## Quantitative determination of creatinine IVD

Store at 2-8 °C

### PRINCIPLE OF THE METHOD

The assay is based on the reaction of creatinine with sodium picrate as described by Jaffé.

Creatinine reacts with alkaline picrate forming a red complex. The time interval chosen for measurements avoids interferences from other serum constituents.

The intensity of the color formed is proportional to the creatinine concentration in the sample<sup>1</sup>.

### CLINICAL SIGNIFICANCE

Creatinine is the result of the degradation of the creatine, component of muscles, it can be transformed into ATP, that a source of high energy for the cells.

The creatinina production depends on the modification of the muscular mass, and it varies little and the levels usually are very stable.

Is excreted by the kidneys. With progressive renal insufficiency there is retention in blood of urea, creatinine and uric acid.

Elevate creatinina level may be indicative of renal insufficiency<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R 1</b> Picric reagent	Picric acid	8.5 mmol/L
<b>R 2</b> Alkaline reagent	Borate buffer	0.29 mol/L
<b>CREATININE CAL</b>	Creatinine aqueous primary standard 2 mg/dL	

### PRECAUTIONS

R1: Corrosive (C); R35: Causes severe burns.

### PREPARATION

All the reagents are ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity. RI (Picric acid) may appear turbid or/and precipitate due to its own composition. However, this will not interfere on the results, as disappears by centrifugation (step 4 from the procedure)
- Blank absorbance (A) at 500 nm  $\geq 1.60$ .

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 500 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

### SAMPLES

- Serum or heparinized plasma<sup>1</sup>.
  - Urine<sup>1</sup>: Dilute sample 1/50 with distilled water. Mix. Multiply results by 50 (dilution factor);
- Stability of the sample: Creatinine is stable for 7 days at 2-8°C.

### PROCEDURE

- Assay conditions:  
Wavelength: ..... 500 nm  
Cuvette: ..... 1 cm. light path  
Temperature: ..... 15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

	Blank	Standard	Sample
R 1 (mL)	--	4.5	4.5
Standard <sup>(Note 1,2)</sup> (mL)	--	0.5	--
Sample (mL)	--	--	0.5

- Shake each tube during 15 sec. And centrifuge at 2500 r.p.m. during 10 min.

### 5. Pipette :

	Blank	Standard	Sample
R 1 (mL)	2.5	--	--
Supernatant (mL)	--	2.5	2.5
R 2 (mL)	1.0	1.0	1.0

- Mix and incubate 20 min. at room temperature.
- Read the absorbance (A) of the samples and calibrator, against the Blank.

### CALCULATIONS

$$\frac{(A)_{\text{Sample}}}{(A)_{\text{Calibrator}}} \times 2 \text{ (Calibrator conc.)} = \text{mg/dL of creatinine in the sample}$$

**Conversion factor:** mg/dL x 88.4 =  $\mu\text{mol/L}$ .

### QUALITY CONTROL

Serum controls are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES<sup>1</sup>

Serum or plasma:

Male 0,7 - 1,4 mg/dL  $\cong$  61.8 – 123.7  $\mu\text{mol/L}$

Female 0,6 - 1,1 mg/dL  $\cong$  53.0 – 97.2  $\mu\text{mol/L}$

Urine:

Male 10 - 20 mg/Kg/24 h  $\cong$  88– 177  $\mu\text{mol/Kg/24 h}$

Female 8 – 18 mg/Kg/24 h  $\cong$  71– 177  $\mu\text{mol/Kg/24 h}$

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 0.12 mg/dL to linearity limit of 10 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

### Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	0.99	2.97	1.06	2.94
SD	0.04	0.08	0.03	0.05
CV (%)	4.09	2.68	2.89	1.70

**Sensitivity:** 1 mg/dL = 0.07 (A)

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples was the following:

Correlation coefficient (r): 0,99.

Regression equation  $y = 0,9274x + 0,0937$

The results of the performance characteristics depend on the used analyzer.

### INTERFERENCES

No interferences were observed to hemoglobin up to 1 g/L and bilirubin up to 10 mg/dL<sup>2</sup>.

A list of drugs and other interfering substances with creatinine determination has been reported by Young et. al<sup>2,3</sup>.

### NOTES

- CREATININE CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- Calibration with the aqueous Standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### BIBLIOGRAPHY

- Murray R.L. Creatinine. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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