
Quantitative determination of urea IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

 Urea in the sample is hydrolyzed enzymatically into ammonia (NH₄⁺) and carbon dioxide (CO₂).

 Ammonia ions formed reacts with α-ketoglutarate in a reaction catalysed by glutamate dehydrogenase (GLDH) with simultaneous oxidation of NADH to NAD⁺:

 The decrease in concentration of NADH, is proportional to urea concentration in the sample¹.

CLINICAL SIGNIFICANCE

Urea is the final result of the metabolism of proteins; It is formed in the liver from their destruction.

 It can appear the urea elevated in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	TRIS pH 7.8	80 mmol/L
Buffer	α-Ketoglutarate	6 mmol/L
	Urease	75000 U/L
R 2	GLDH	60000 U/L
Enzymes	NADH	0.32 mmol/L
	UREA CAL	Urea aqueous primary standard 50 mg/dL

PREPARATION

 Working reagent (WR): Mix **4 vol. R1 Buffer + 1 vol. R2 Substrate**.

The (WR) is stably for 1 month at 2-8°C.

UREA CAL: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1.00.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm..
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment^(Note 1).

SAMPLES

 - Serum or heparinized plasma¹: Do not use ammonium salts or fluoride as anticoagulants.

 - Urine¹: Dilute sample 1/50 in distilled water. Mix. Multiply the results by 50 (dilution factor). Preserve urine samples at pH < 4.

Urea is stable at 2-8°C for 5 days.

PROCEDURE

- Assay conditions:
 Wavelength: 340 nm
 Cuvette: 1 cm light path
 Temperature 37°C / 15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

	Blank	Standard	Sample
WR (mL)	1.0	1.0	1.0
Standard ^(Note 2-3) (μL)	--	10	--
Sample (μL)	--	--	10

- Mix and read the absorbance after 30 s (A₁) and 90 s (A₂).
- Calculate: ΔA= A₁ - A₂ .

CALCULATIONS

$$\frac{(\Delta A)_{\text{Sample}}}{(\Delta A)_{\text{Calibrator}}} \times 50 (\text{Calibrator conc}) = \text{mg/dL urea in the sample}$$

 10 mg/L urea BUN divided by 0.466 = 21 mg/L urea = 0.36 mmol/L urea¹.

Conversion factor: mg/dL x 0.1665 = mmol/L.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES^{4,5}

Serum or plasma:

15-45 mg/dL ≅ 2.5-7.5 mmol/L

Urine:

26 - 43 g/24 h ≅ 428-714 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS
Measuring range: From detection limit 1 mg/dL to linearity limit 350 mg/dL.

If the concentration is greater than linearity limit dilute 1:2 the sample with ClNa 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	40.6	141	42.5	141
SD	1.22	1.03	2.12	1.15
CV (%)	2.99	0.73	4.99	0.81

Sensitivity: 1 mg/dL = 0,00087 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples was the following:

Correlation coefficient (r): 0.99.

Regression equation y= 0.9993x + 0.0394.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

 It is recommended to use heparin as anticoagulant. Do not use ammonium salts or fluoride¹.

 A list of drugs and other interfering substances with urea determination has been reported by Young et. al^{2,3}.

NOTES

- UREA CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- Glassware and distilled water must be free of ammonia and ammonium salts¹.
- Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
- Use clean disposable pipette tips for its dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- Kaplan A. Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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