**Quantitative determination of creatinine IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

In the first reaction, creatinase and sarcosine oxidase were used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide, that is eliminated by catalase. In the second reaction, the catalase is inhibited by sodium azide, and creatinase and 4-aminoantipyrine (4-AA) were added, and only the creatine generated from creatinine by creatinase was hydrolyzed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly-formed hydrogen peroxide was measured in a coupled reaction catalyzed by peroxidase, with N-ethyl-n-sulphopropyl-mtoluidine (TOPS)/4-AA as a chromogen.

**CLINICAL SIGNIFICANCE**

Creatinine measurements are used in the diagnosis and treatment of renal diseases, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes. Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

**REAGENTS**

<b>R 1</b>	MOPS 25 mmol/L, TOPS 0.5 mmol/L, Creatinase 10 KU/L, Sarcosine Oxidase 5 KU/L, Catalase 3 KU/L, EDTA 1mmol/L, pH 7,5.
<b>R 2</b>	MOPS 90 mmol/L, Creatinase 30 KU/L, peroxidase KU/L, pH 7,5. Azida sódica 0,5 g/L.
<b>CREATININE CAL</b>	Creatinine aqueous primary standard

**PREPARATION**

R1 and R2 are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. R1 and R2 are stable 8 weeks after opening bottle.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer measuring at 545±20 nm
- Cell holder thermostatable at 37°C
- General laboratory equipment.

**SAMPLES**

- Serum or plasma<sup>1</sup>.
- Urine<sup>1</sup>: Dilute fresh urine 1/50 with distilled water. Creatinine is stable 1 day at 2-8°C.

**PROCEDURE**

- Assay conditions:  
Wavelength: ..... 545 nm (525-565)  
Temperature: ..... 37°C (±0.1°C)
- Bring the reagents and the instrument to 37°C (±0.1°C).
- Pipette into a thermostated cuvette:

	Blank	Calibrator	Sample
R1 (µL)	450	450	450
H2O/Calibrator/Sample (µL)	10	10	10

- Mix and incubate 5 minutes.
- Read the absorbance (A<sub>1</sub>) of the calibrator and the samples, at 545nm against the blank.
- Add:

	Blank	Calibrator	Sample
R2 (µL)	150	150	150

- Mix and incubate 5 minutes.
- Read the absorbance (A<sub>2</sub>) of the calibrator and the samples, at 545nm against the blank.

**CALCULATIONS**

$$\text{Creatinine} = \frac{\Delta A_{\text{Sample}} \times k - \Delta A_{\text{Blank}} \times k}{\Delta A_{\text{Calibrator}} \times k - \Delta A_{\text{Blank}} \times k} \times C$$

$$K = 0.754 = 276 \mu\text{L} / 366 \mu\text{L}$$

C= Concentration of the calibrator

$$\Delta A = A_2 - A_1$$

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 y 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:

Men 0,9 - 1,3 mg/dL

Women 0,6 - 1,1 mg/dL

Urine:

Men 14- 26 mg/Kg/24 h

Women 11-20 mg/Kg/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 0.03 mg/dL to linearity limit of 150 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	CV (%)	
Mean (mg/dL)	0.982	3.776	0.969	3.753
SD	0.030	0.105	0.052	0.150
CV (%)	3.092	2.768	5.048	4.008

**Sensitivity:** 1 mg/dL = ΔA 0,0015 A/min . mg/dL.

**Accuracy:** Results obtained using SPINREACT these reagents did not show systematic differences when compared with other commercial reagents or with HPLC method.

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0.977.

Regression equation: y= -0.093x + 0.945.

The results of the performance characteristics depend on the analyzer used.

**INTERFERENCES**

No interferences were observed with hemoglobin until 5 g/dL, bilirubin 40 mg/dL.

Other drugs and substances may interfere<sup>2,3</sup>.

**NOTES**

- Calibration with an aqueous standard may cause matrix related bias, it is recommended to calibrate using a serum based calibrator.
- Urine: Multiply the result by 50 (sample dilution factor).

**BIBLIOGRAPHY**

- Fossati et al. Clin Chem 1983;29:1494-1496.
- Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co.,1999.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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