Quantitative determination of total Homocystein (tHcy) IVD

Store at 2-8°C

INTENDED USE

The Homocystein Reagent is intended for in vitro quantitative determination of total homocystein in human serum and plasma. The device can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia and homocystinuria.

PRINCIPLE OF THE METHOD

Bound or dimerised homocystein (oxidised form) is reduced to free homocystein, which then reacts with serine catalysed by cystathionin beta-synthase (CBS) to form cystathionin. Cystathionin in turn is broken down by cystathionin beta-lyase (CBL) to form homocystein, pyruvate and ammonia. Pyruvate is then converted by lactate dehydrogenase (LDH) to lactate with nicotinamide adenine dinucleotide (NADH) as coenzyme. The rate of NADH conversion to NAD⁺ is directly proportional to the concentration of homocystein ($\Delta A_{340\text{ nm}}$).

CLINICAL SIGNIFICANCE

Homocystein is a thiol-containing amino acid produced by the intracellular demethylation of methionin. Total homocystein (tHcy) represents the sum of all forms of Hcy including forms of oxidized, protein bound and free.

Elevated level of tHcy has emerged as an important risk factor in the assessment of cardiovascular disease. Excess Hcy in the blood stream may cause injuries to arterial vessels due to its irritant nature, and result in inflammation and plaque formation, which may eventually cause blockage of blood flow to the heart. Elevated tHcy levels are caused by four major factors, including: a) genetic deficiencies in enzymes involved in Hcy metabolisms such as cystathionin beta-synthase (CBS), methionin synthase (MS) and methyltetrahydrofolate reductase (MTHFR); b) nutritional deficiency in B vitamins such as B6, B12 and folate; c) renal failure for effective for amino acid clearance, and d) drug interactions such as nitric oxide, methotrexate and phenytoin that interfere with Hcy metabolisms.

Elevated levels of tHcy are also linked with Alzheimer's disease and Osteoporosis

REAGENTS

R1	NADH (0.47 mM), LDH (38 KU/L), Serine (0.76 mM), Trizma Base 1-10%, Trizma Hydrochloride 1-10%, Sodium Azide < 1%, Reductant (TCEP:2.9 mM)
R2	Cycling Enzymes CBS (0.748 KU/L) and CBL (16.4 KU/L) Sodium Azide < 1%.
Optional	Ref. 1002276 Calibrator 2 x 3.0 mL Ref. 1002281 Control 3 x 1.5 mL

PRECAUTIONS

R1 and R2: Xn. R22: Harmful if swallowed. R32: Contact with acids liberates very toxic gas. S36/37/39: Wear suitable protective clothing, gloves and eye/face protection. S29/35: Do not empty into drains; dispose of this material and its container in a safe way. S46: If swallowed, seek medical advice immediately and show this container or label.

CALIBRATION

Spinreact recommends using tHcy Calibrator (Ref: 1002276). The calibrators are prepared gravimetrically and are traceable to NIST SRM 1955. It is not recommended the use of other commercially available tHcy calibrators.

STORAGE AND STABILITY

- All the kit components are stable until the expiration date on the label, when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.
- Do not mix the reagents of different lots.
- **Do not freeze reagents.**
- Do not expose Reagent 1 and Reagent 2 to light during on-board use.
- Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.
- The reagents should be clear of particulate material. They should be discarded if they become turbid.

ADDITIONAL EQUIPMENT


- Spectrophotometer or colorimeter measuring at 340 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

- Fresh serum or heparin plasma. It is important to centrifuge blood samples immediately after collection to separate the plasma from the blood cells. If immediate centrifugation is not possible, collected blood specimens should be kept on ice and centrifuged within an hour.
- Hemolysed, turbid or lipemic specimens are not recommended.
- If the assay will be performed within 2 weeks after collection, the specimen should be stored at 2-8°C. If the testing will be delayed more than 2 weeks, the specimen should be stored frozen at -20°C or colder. Specimens have been shown to be stable at -20°C for 8 months.^{1,2}
- Mix specimens thoroughly after thawing by low speed vortexing or by gentle inversion to ensure consistency in results. Avoid repeated freezing and thawing.

PROCEDURE

- All components are ready to use.

CALCULATIONS

Plot (A2- A1) obtained against the tHcy concentration of both calibrators. tHcy concentration in the sample is calculated by interpolation of its (A2- A1) in the linear calibration curve.

QUALITY CONTROL

It is recommended to validate the performance of Hcy reagents with an available set of Hcy Control (Ref. 1002281).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own acceptable tolerances for the Control.

REFERENCE VALUES

Scientific literature reports reference values for adult male and females between 5 and 15 µmol/L, men having higher values than women, and post menopausal woman having higher homocystein values than pre-menopausal women.^{1,3,4}

These values are orientative. The reference range should be determined by each laboratory to confirm the characteristics of the population being tested.

PERFORMANCE CHARACTERISTICS

Limit of detection: 0.33 µmol / L.

Measurement range: The linear range is 1-46 µmol/L. Specimens > 46 µmol/L should be diluted 1 part specimen to 2 parts Cal 0 µmol/L or 1 part specimen to 9 parts Cal 0 µmol/L as appropriate.

Precision:

The study was conducted according to the NCCLS EP-5 protocol. Three HCY controls were assayed, in replicates of two, at two separate times per day for 20 days on one instrument (n=80).

Level of Hcy control	Within Precision			Total Precision		
	Low	Medium	High	Low	Medium	High
Mean (µM)	6.3	12.3	25.5	6.3	12.3	25.5
CV (%)	2.6	1.5	1.5	4.4	3.0	2.5

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 94 samples were the following:

Correlation coefficient r= 0.999.

Regression equation: y= 0.165+ 0.991x.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

The following substances produced less than 10% deviation when tested with the following stated concentrations: bilirubin 20 mg/dL, haemoglobin 500 mg/dL, red blood cell 0.4%, triglyceride 500 mg/dL, glutathione 1000 µM, methionine 800 µM, cysteine 200 µM, pyruvate 1250 µM. Samples with raised protein levels show > 10% difference compared to results obtained from normal samples and should be avoided.

LIMITATIONS

- Cystathionin is measured with homocystein, but in the general population the cystathionin level (0.065 to 0.3 µmol/L) has a negligible effect. In very rare cases, end stage renal disease and patients with severe metabolic disturbances, cystathionin levels may rise dramatically and in severe cases cause greater than 20% interference.
- Hydroxylamine, present in several iron reagents may carryover (reagent probe or reaction cuvette) and cause falsely low results.
- Carbamazepine, methotrexate, phenytoin, nitrous oxide, or 6-azauridine triacetate may affect the homocystein concentration.¹
- Specimens from patients who are on drug therapy involving S-adenosyl-methionin may show falsely elevated levels of homocystein.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

1. Ueland PM, Refsum H, Stabler SP, et al. Total Homocystein in Plasma or Serum: Methods and Clinical Applications. *Clin Chem* 1993;39:1764-1779
2. Fiskerstrand T, Refsum H, Kvalheim G, et al. Homocystein and Other Thiols in Plasma and Urine: Automated Determination and Sample Stability. *Clin Chem* 1993;39:263-271
3. Nehler MR, Taylor LM Jr, Porter JM. Homocysteinemia as a Risk Factor for Atherosclerosis: A Review. *Cardiovascular Pathol* 1997;6:1-9
4. Lussier-Cacan S, Xhignesse M, Piolot A, et al. Plasma Total Homocystein in Healthy Subjects: Sex-Specific Relation with Biological Traits. *Am J Clin Nutr* 1996;64:587-593

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1. Refsum H, Stabler SP, et al. Total Homocysteine in Plasma or Serum: Methods and Clinical Applications. *Clin Chem* 1993;39:1764-1779
2. Fiskerstrand T, Refsum H, Kvalheim G, et al. Homocysteine and Other Thiols in Plasma and Urine: Automated Determination and Sample Stability. *Clin Chem* 1993;39:263-271
3. Nehler MR, Taylor LM Jr, Porter JM. Homocysteinemia as a Risk Factor for Atherosclerosis: A Review. *Cardiovascular Pathol* 1997;6:1-9
4. Lussier-Cacan S, Xhignesse M, Pilot A, et al. Plasma Total Homocysteine in Healthy Subjects: Sex-Specific Relation with Biological Traits. *Am J Clin Nutr* 1996;64:587-593

PRESENTACION

Ref:41057

Cont.

R1: 1x 30 mL

R2: 1x 5 mL