## Diluent for Spincell 5 Series IVD

### INTENDED USE

This product is designed to be used on a Spinreact 5-Part-Diff haematology analyser along with Lyse and Sheath reagents produced by SPINREACT, for testing the cell content in human blood.

### FUNCTIONS

Diluent is a buffer solution for diluting cells and preventing accumulation. Its osmotic concentration, ionic strength and conductivity allows the preservation of cells integrity and volume.

### INGREDIENTS

Anhydrous Sodium Sulfate	5~15 g/L
Preservative Solution	<20 g/L
Muriate	2~6 g/L
Edetate disodium	<0.5g/L

### STORAGE & STABILITY

This product will be stable for 24 months when stored unopened at 2°C~40°C and away from sunlight. Once opened, this product should be consumed within 60 days and stored below 35°C

### BLOOD SAMPLES

1. The blood samples should be venous and capillary specimens collected in EDTA-K<sub>2</sub> anticoagulant. Coagulated, hemolysis and lipidaemia samples cannot be measured.
2. Blood samples should be stored at 2°C~8°C.

### OPERATIONS

1. Place the Diluent container beneath the analyzer to avoid direct sunlight. Keep away from cooling or heating air outlet.
2. Use it at 15°C~30°C.
3. Connect Diluent container with the analyzer according to Operation Manual.

### TEST RESULT


If the test result is out of the range of reference, flagged results will be marked with letter "L" or "H". "L" means the result is below the lower limit; "H" means the result is higher than the upper limit.

### RESTRICTIONS, PRECAUTIONS & NOTICE

1. Use the reagent under the professional technician's guidance.
2. Make sure that the reagent matches with your analyzer. Any mismatch will lead to misleading results.
3. Avoid contacting with skin and eyes. In case of contact, rinse immediately with plenty of water and seek medical advice.
4. Avoid inhaling the evaporated reagent.
5. Tight up the container's cap to prevent from volatilization and contamination.
6. Once opened, dispose the remains over 60 days.

7. If product has been frozen, let product warm to room temperature, and then gently invert product to mix it completely. Verify background results before use.

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