Lyse reagent for Spincell 5 Series IVD

INTENDED USE

This product is designed to be used on a Spinreact 5-Part-Diff haematology analyser along with Lyse and Sheath produced by SPINREACT, for testing the cell content in human blood.

FUNCTIONS

Lyse does not contain trinitride or cyanide. It has the functions of lysing the red cells and releasing the hemoglobin to protect the RBC count from being affected by RBC fragments. At the same time it converts hemoglobin to steady colored hemo-compound for detecting the concentration of hemoglobin. The Lyse preserves the integrity of WBC nucleus allowing the counting and classification of WBCs

INGREDIENTS

Dodecyl trimethyl ammonium chloride	25~30 g/L
Myristyltrimethylammonium bromide	1,0~2,0g/L
Hydroxylamine hydrochloride	1,0~2,0 g/L

PRECAUTIONS

H290-May be corrosive to metals. EUH208-Contains Hydroxylamine hydrochloride. May produce an allergic reaction.

Follow the precautionary statements given in MSDS and label of the product.

STORAGE & STABILITY

This product will be stable for 24 months when stored unopened at 2°C~40°C and away from sunlight. Once opened, this product should be consumed within 60 days and stored below 35°C

BLOOD SAMPLES

1. The blood samples should be venous and capillary specimens collected in EDTA-K₂ anticoagulant. Coagulated, hemolysis and lipidaemia samples cannot be measured.
2. Blood samples should be stored at 2°C~8°C.

OPERATIONS

1. Place the Lyse container beneath the analyzer to avoid direct sunlight. Keep away from cooling or heating air outlet.
2. Use it at 15°C~30°C.
3. Connect Lyse container with the analyzer according to Operation Manual.

TEST RESULT


If the test result is out of the range of reference value, flagged results will be marked with letter "L" or "H". "L" means the result is below the lower limit; "H" means the result is higher than the upper limit.

RESTRICTIONS, PRECAUTIONS & NOTICE

1. Use the reagent under the professional technician's guidance.
2. Make sure that the reagent matches with your analyzer. Any mismatch will lead to misleading results.
3. Avoid contacting with skin and eyes. In case of contact, rinse immediately with plenty of water and seek medical advice.

4. Avoid inhaling the evaporated reagent.
5. Tight up the container's cap to prevent from volatilization and contamination.
6. Once opened, dispose the remains over 60 days.
7. If product has been frozen, let product warm to room temperature, and then gently invert product to mix it completely. Verify background results before use.

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