## Reticulocyte reagent for Spincell 5 Series IVD

### INTENDED USE

Stain is used for measuring the content of reticulocyte (RETIC) in human blood together with SPINREACT Sheath on adaptable hematology analyzers.

### RESTRICTIONS, PRECAUTIONS & NOTICES

1. Use the detergent under the professional technician's guidance.
2. This product should not be consumed, if so, please see a doctor immediately.
3. Avoid contact with skin and eyes. In case of contact, rinse immediately with plenty of water or seek medical advice.
4. Deal with the waste according to the local regulations.

### INGREDIENTS

Sodium chloride, potassium oxalate, buffer solution, preservative, surfactant and new methylene blue.

### STORAGE & STABILITY

1. Storage Condition: keep away from direct sunlight and store between 2°C and 40°C. Store it in the temperature between 15°C and 30°C once opening.
2. Shelf Life: 24 months.

### APPLICABILITY

It's appropriate to series of Spinreact 5-Part-Diff Auto Hematology Analyzer.

### BLOOD SAMPLES

1. The blood samples should be venous and capillary specimens collected in EDTA-K<sub>2</sub> anticoagulant. Coagulated, hemolysis and lipidaemia samples cannot be measured.
2. Blood samples should be stored at 2°C~8°C.

### OPERATION PROCEDURE

1. Add 20µL whole blood into one tube of RETIC Stain. Mix it upside down, and place the mixture in the temperature between 15°C to 30°C for 15 minutes.
2. Operate it according to the Instruction of 5-Part-Diff Auto Hematology Analyzer.

### PERFORMANCE

1. pH value: 7.20 ± 0.20 at (25 ± 1) °C.
2. Osmotic concentration (560 ± 10) mOsm / kg.
3. Peak wavelength of absorption:  $\lambda_{max} = (675 \pm 10)nm$ .
4. Absorbance value: 3.20 ± 0.30.
5. Blank value  $\leq 0.01 \times 10^{12}/L$ .
6. Repeatability  $\leq 20\%$ .
7. Scatter plot: test the fresh blood by applicable analyzer, the scatter plot which shall show the categories of scatter group of reticulocytes and RBC.
8. Batch differences:  $\Delta PH$  value  $\leq 0.40$ , osmotic concentration  $\leq 10mOsm/kg$ ,  $\Delta\lambda_{max} \leq 10nm$ .

### TEST RESULTS

1. Test results please take reference of the normal range values.
2. If the test result is out of the range of reference value, flagged results will be marked with letter "L" or "H". "L" means the result is below the lower limit; "H" means the result is higher than the upper limit.

### CHANGES ON PROCESSES AND PERFORMANCE

1. Test the processed sample within 2 hours.
2. Please note that the reagent shall match with analyzer model. Please do not use the reagent produced by other company, otherwise it may cause measurement errors.

### QUALITY CONTROL

Detect three pieces of products in each production lot, two of them should be the first box and the last box in manufacturing process. The production lot shall be determined to be unqualified if there one does not meet inspection standards.


### QUALITY CONTROL AND TRACEABILITY OF CALIBRATORS

Not applicable.

### REFERENCE VALUE

ITEM	Adult Male	Adult Female	Child
RETIC%	0.5%~2.0%		
RETIC ( $\times 10^{12}/L$ )	34.8~148	21.1~123	49.5~94.6
IRF	0.1~0.38		0.02~0.50

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