Quantitative determination of fecal occult blood IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

FOB-Turbilatex is a quantitative turbidimetric test for the measurement of human hemoglobin in human feces.

Latex particles coated with specific anti-human hemoglobin are agglutinated when mixed with samples containing hemoglobin or hemoglobin/haptoglobin complexes. The agglutination causes an absorbance change, dependent upon the hemoglobin contents of the patient sample that can be quantified by comparison from a calibrator of known hemoglobin concentration.

CLINICAL SIGNIFICANCE

The presence of hemoglobin in stool is associated with lower intestinal tract like colon polyps, Crohn's disease, ulcerous colitis and colorectal cancer. An early diagnostic is the key for a successful therapy. FOB-Turbilatex allows detecting amounts of hemoglobin that are invisible to the naked eye. The test can be performed with an automatic analyzer using directly a stool sample collection tube. These features make this the ideal test for screening programs.

REAGENTS

Diluent (R1)	Buffered solution. Preservative
Latex (R2)	Latex particles coated with polyclonal antibody anti-human hemoglobin. Preservative.
FOB-CAL	Calibrator. Hemoglobin concentration is stated on the vial label.

PRECAUTIONS

Contains products of animal origin. Handle cautiously as potentially infectious.

CALIBRATION

Use FOB Calibrator Ref. 1107182. The value of the calibrator is traceable to Reference Material WHO International Standard Haemoglobinocyanide NIBSC (98/708). Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

The reagents are ready to use. FOB Calibrator and controls: reconstitute (→) lyophilized products with 2 mL of purified water.

Calibration curve: Prepare the following FOB calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the FOB calibrator by the corresponding factor stated in table below to obtain the Hemoglobin concentration of each dilution:

Calibrator dilution	0	1/16	1/8	1/4	1/2	1
FOB Calibrator (µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	-
Factor	0	0,062	0,125	0,25	0,5	1,0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date. Do not freeze. Reagent deterioration: Presence of particles and turbidity.

Calibrator/controls stability: after reconstitution, 24h at 2-8°C or 1 month at -20°C.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E autoanalyzer.
- Laboratory equipment.

SAMPLES

Use only Spinreact's sample collection tube, following the instructions included in the package. Store stool samples in collection tube during maximum 7 days at 2-8°C.

REFERENCE VALUES

Each laboratory should determine its own normal reference value. The sensitivity and specificity of the results with the cutoff of 100ng/mL was 76% and 95%, respectively.

QUALITY CONTROL

Control solutions are recommended to monitor the performance of assay procedures (Control 1 and Control 2). If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

The clinical diagnostic should not be based on a single test result; the complete medical history of the patient should be taken into consideration.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS

Test	FOB	R1	180
No.	**	R2	60
Full Name	FOB	Sample volume	8.5
Standard No.	6	R1 Blank	
Reac. Type	Fixed Time	Mixed Rgt Blank	
Pri. Wavelength	510 nm	Linearity Range	0 ng/mL 1100 ng/mL
Sec. Wavelength		Linearity Limit	*
Direction	Increase	Substrate Limit	*
Reac. Time	1 17	Factor	
Incuba. Time		Prozone check	*
Unit	ng/mL	q1	q2
Precision	0.01	q3	q4
		PC	Abs

CALIBRATION (Cal + Rgt BIK)

Rule	Spline
Sensitivity	1
Replicates	2
Interval (days)	0
Difference Limit	
SD	
Blank Response	
Error Limit	
Correlation Coefficient	

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **30 days**. After this period the blank parameter must be performed again in order to validate the calibration.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to the maximum value of the calibrator (aprox. 1000 ng/mL), under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in sample diluent (obtained from sample collection tube). The linearity limit depends on the sample/reagent ratio and the analyzer used.
- Detection limit:** Values less than 10 ng/mL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 10000 ng/mL.
- Precision:** The reagent has been tested during 19 days on a Spintech240, using three different hemoglobin concentrations in a EP5 (CLSI)-based study.

	CV (%)	
	75 ng/mL	232 ng/mL
Total	5,2	5,0
Within Run	2,4	0,9
Between Run	3,4	2,2
Between Days	3,0	4,5

- Accuracy:** 108 fecal samples were tested, including 14 positives. All results were reported correctly, at the cutoff of 100 ng/mL. Results were compared against a validated immunocromatographic method.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Weizman AV. et al., Minerva Gastroenterol Dietol., 2010.
2. Burt RW et al., J. Natl. Compr. Canc. Netw., 2013
3. Kuipers EJ et al., Nat.Rev.Clin.Oncol., 2014.
4. Vilkin A et al., Am. J. Gastroent., 2005.

PACKAGING

Ref: MI1107180	R1: 2 x 15 mL
	R2: 1 x 10 mL
Cont.	FOB-CAL :1 x 2 mL
Ref: MI1107181	R1: 4 x 30 mL
	R2: 2 x 20 mL
	FOB-CAL :1 x 2 mL

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