

**Quantitative determination of triglycerides IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Sample triglycerides incubated with lipoproteinlipase (LPL), liberate glycerol and free fatty acids. Glycerol is converted to glycerol-3-phosphate (G3P) and adenosine-5-diphosphate (ADP) by glycerol kinase and ATP. Glycerol-3-phosphate (G3P) is then converted by glycerol phosphate dehydrogenase (GPO) to dihydroxyacetone phosphate (DAP) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>).

In the last reaction, hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) reacts with 4-aminophenazone (4-AP) and p-chlorophenol in presence of peroxidase (POD) to give a red colored dye:

The intensity of the color formed is proportional to the triglycerides concentration in the sample<sup>1,2,3</sup>.

**CLINICAL SIGNIFICANCE**

Triglycerides are fats that provide energy for the cell. Like cholesterol, they are delivered to the body's cells by lipoproteins in the blood. A diet with a lot of saturated fats or carbohydrates will raise the triglyceride levels. The increases in serum triglycerides are relatively non-specific. For example liver dysfunction resulting from hepatitis, extra hepatic biliary obstruction or cirrhosis, diabetes mellitus is associated with the increase<sup>3,6,7</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R</b>	GOOD pH 6.3	50 mmol/L
	p-Chlorophenol	2 mmol/L
	Lipoprotein lipase (LPL)	150000 U/L
	Glycerol kinase (GK)	500 U/L
	Glycerol-3-oxidasa (GPO)	3500 U/L
	Peroxidase (POD)	440 U/L
	4 - Aminophenazone (4-AP)	0,1 mmol/L
ATP	0,1 mmol/L	
<b>Optional</b>	SPINTROL H CAL	

**PREPARATION**

Ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm ≥ 0,40.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or colorimeter measuring at 505 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

**SAMPLES**

 Serum or plasma<sup>1</sup>.

Stability of the sample: Triglycerides are stable for 5 days at 2-8°C.

**APPLICATION SPINLAB 180**

Name	TRIGLYCERIDES	Ref. male low	40 mg/dL
Abbr. Name	TRIG	Ref. male high	160 mg/dL
Mode	Endpoint	Ref. female low	35 mg/dL
Wavelength	505 nm	Ref. female high	135 mg/dL
Units	mg/dL	Ref. Ped. Low	_ mg/dL
Decimals	0	Ref. Ped. High	_ mg/dL
Low Conc.	0	Panic value low	-
High Conc.	1000	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
Normal volume	300 µL		
Rerun volume	300 µL		
Sample			
Normal volume	3,0 µL		
Rerun volume	2,0 µL		
R2 bottle (mL)	0 mL		
Normal volume	0 µL		
Rerun volume	0 µL		
Predilution	No		
Slope blank	No		
Incubation	4,5 min.		
<b>Factor</b>			
Reagent blank	Yes		
Low Absorbance	-0,100 Abs		
High Absorbance	3,000 Abs		
R. Abs. L. Limit	-0,100 Abs		
R. Abs. H. Limit	3,000 Abs		
R. Abs. Deviation	3,000 Abs		

**QUALITY CONTROL**

Control Sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES**

Men	40 – 160 mg/dL
Women	35 – 135 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

**NOTES**

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

**BIBLIOGRAPHY**

1. Buccolo G et al. Quantitative determination of serum triglycerides by use of enzymes. Clin Chem 1973; 19 (5): 476-482.
2. Fossati P et al. Clin. Chem 1982; 28(10): 2077-2080.
3. Kaplan A et al. Tryglicerides. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 437 and Lipids 1194-1206.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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1. Buccolo G et al. Quantitative determination of serum triglycerides by use of enzymes. Clin Chem 1973; 19 (5): 476-482.
2. Fossati P et al. Clin. Chem 1982; 28(10): 2077-2080.
3. Kaplan A et al. Tryglycerides. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 437 and Lipids 1194-1206.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
6. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP41031

Cont.

R: 10 x 25 mL