Quantitative determination of Ferritin IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Ferritin-turbilatex is a quantitative turbidimetric test for the measurement of ferritin in human serum or plasma.

Latex particles coated with specific anti-human ferritin are agglutinated when mixed with samples containing ferritin. The agglutination causes an absorbance change, dependent upon the ferritin contents of the sample that can be quantified by comparison from a calibrator of known ferritin concentration.

CLINICAL SIGNIFICANCE

Serum ferritin concentration usually reflects body iron stores and is considered one of the most reliable indicators of iron status of patients. Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, pH 8.2. Sodium azide 0.95 g/L.
Latex (R2)	Latex particles coated with anti human ferritin, pH 8.2. Sodium azide 0.95 g/L.
FERR-CAL	Calibrator. Ferritin concentration is stated on the vial label.
Optional	Ref.:1107044 Ferritin control.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Use Ferritin Calibrator Reference 1107042.

The sensitivity of the assay and the target value of the calibrator have been standardized against the 2nd International Standard of Ferritin (80/578, 1992 WHO). Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ferritin Calibrator: Reconstitute (→) with 3.0 mL of distilled water. Mix gently and incubate at room temperature for about 10 minutes before testing.

Calibration curve: Prepare the following FERR calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the FERR calibrator by the corresponding factor stated in table below to obtain the concentration of each dilution.

Calibrator dilution	1	2	3	4
FERR Calibrator (µL)	--	33.3	66.6	100
NaCl 9 g/L (µL)	100	66.6	33.3	--
Factor	0	1/3	2/3	1

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: Presence of particles and turbidity.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

SPINLAB 180 APPLICATION

Name	Ferritin	Ref. male low	*
Abbr. Name	FERR	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	µg/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	0 µg/L	Panic value low	*
High Conc.	400 µg/L	Panic value high	*
Calibrator name	CAL FERR(4 pts)	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	180 µL		
rerun volume	180 µL		
Sample			
normal volume	20.0 µL		
rerun volume	20.0 µL		
R2 bottle (mL)	5 mL		
normal volume	45 µL		
rerun volume	45 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Ferritin Control (Ref.: 1107044).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Men: 30 – 300 µg/L.

Women: 30 – 300 µg/L /Pregnant: 12-160 µg/L.

Children: 15-120 µg/L.

Each laboratory should establish its own reference range.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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