**Quantitative determination of Rheumatoid Factors (RF)
 IVD**

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The RF-Turbilatex is a quantitative turbidimetric test for the measurement of RF in human serum or plasma.

Latex particles coated with human gammaglobulin are agglutinated when mixed with samples containing RF. The agglutination causes an absorbance change, dependent upon the RF contents of sample that can be quantified by comparison from a calibrator of known RF concentration.

CLINICAL SIGNIFICANCE

Rheumatoid factors are a group of antibodies directed to determinants in the Fc portion of the immunoglobulin G molecule. Although rheumatoid factors are found in a number of rheumatoid disorders, such as systemic lupus erythematosus (SLE) and Sjögren's syndrome, as well as in nonrheumatic conditions, its central role in clinic lies its utility as an aid in the diagnosis of rheumatoid arthritis (RA).

A study of the "American College of Rheumatology" shows that the 80.4% of RA patients were RF positive.

REAGENTS

| | |
|---------------------|---|
| Diluent (R1) | Tris buffer 20 mmol/L, pH 8.2. Sodium azide 0.95 g/L. |
| Latex (R2) | Latex particles coated with human gammaglobulin, pH 7.4. Sodium azide 0.95 g/L. |
| RF-CAL | Calibrator. Human serum. The RF concentration is stated on the vial label. |
| Optional | Ref.:1102114 Control serum ASO/CRP/RF Level L Ref.:1102115 Control serum ASO/CRP/RF Level H. |

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Use RF Calibrator Reference 1107007.

The sensitivity of the assay and the target value of the calibrator have been standardized against the International Reference NIBSC 64/2 (Rheumatoid Arthritis Serum) WHO.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use

RF Calibrator: Reconstitute (→) with 2.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.

Calibration Curve: Prepare the following RF calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the RF calibrator by the corresponding factor stated in table below to obtain the RF concentration of each dilution.

| Calibrator dilution | 1 | 2 | 3 | 4 | 5 | 6 |
|---------------------|-----|--------|-------|------|-----|-----|
| Calibrator RF (µL) | -- | 25 | 50 | 100 | 200 | 400 |
| NaCl 9 g/L (µL) | 400 | 375 | 350 | 300 | 200 | - |
| Factor | 0 | 0,0625 | 0,125 | 0,25 | 0,5 | 1,0 |

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: Presence of particles and turbidity.

Reconstituted calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

Do not freeze; frozen latex and diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolized or lipemic samples.

SPINLAB 180 APPLICATION

| | | | |
|-------------------|-----------|------------------|-------|
| Name | RF | Ref. male low | * |
| Abbr. Name | RF | Ref. male high | * |
| Mode | Twopoints | Ref. female low | * |
| Wavelength | 620 nm | Ref. female high | * |
| Units | IU/mL | Ref. Ped. Low | * |
| Decimals | 2 | Ref. Ped. High | * |
| Low Conc. | 6 IU/mL | Panic value low | * |
| High Conc. | 160 IU/mL | Panic value high | * |
| Calibrator name** | CAL RF** | Control 1 | * |
| Prozone check | No | Control 2 | * |
| | | Control 3 | * |
| | | Correlat. factor | 1.000 |
| | | Correlat. offset | 0.000 |

DUAL MODE

| | |
|------------------|--------------|
| Sample blank | No |
| R1 bottle (mL) | 25 mL |
| normal volume | 240 µL |
| rerun volume | 240 µL |
| Sample | |
| normal volume | 3.0 µL |
| rerun volume | 3.0 µL |
| R2 bottle (mL) | 5 mL |
| normal volume | 60 µL |
| rerun volume | 60 µL |
| Predilución | No |
| Slope blank | No |
| Point one,two | -3, 130 sec. |
| Reagent blank | Yes |
| Low Absorbance | -0.100 Abs |
| High Absorbance | 3.000 Abs |
| R. Abs. L. Limit | -0.100 Abs |
| R. Abs. H. Limit | 3.000 Abs |
| Substr.Depletion | 3.000 Abs |

** Adjust maximum inaccuracy calibration % until the curve cover all the Abs found.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used SPINREACT Control ASO/CRP/RF Level L (Ref.: 1102114) and Level H (Ref.: 1102115). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Normal values up to 20 IU/mL. Each laboratory should establish its own reference range.

INTERFERENCES

Hemoglobin (10 g/L), bilirubin (20 mg/dL) and lipemia (10 g/L), do not interfere. Other substances may interfere⁶.

NOTES

1. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951- 960.
2. Robert W Dorner et al. Clinica Chimica Acta 1987; 167: 1-21.
3. Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528 - 534.
4. Vladimir Muié et al. Scand J Rheumatology 1972; 1: 181 - 187.
5. Paul R et al. Clin Chem 1979; 25/11: 1909 - 1914.
6. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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