Stool collection tube

Stool collection tube for FOB-turbidimetric assay IVD

Store at 2-30°C

PRINCIPLE OF THE METHOD

Stool collection tube for fecal occult blood testing to be used with Spinreact's FOB-turbidimetric assay. The buffered solution contained in the tube stabilizes the hemoglobin present in the sample during up to 7 days; time enough from sampling to testing. The tube design standardizes the amount of stool collected for an easier interpretation of the results.

REAGENTS

Sample Buffer	1,7 mL of buffered solution, stabilizers, preservatives.	
Related products	FOB-turbi reagent	Ref.: MI1107180
		Ref.: MI1107181
	Control 1 and Control 2	Ref.: 1107183
	Calibrator	Ref.: 1107182

STORAGE AND STABILITY

The buffered solution contained in the tube is stable until the expiry date printed on the tube, if stored at 2-30°C. The hemoglobin obtained with the sample is stable 7 days if the tube is stored at 2-8°C.

PROCEDURE

1.

Place toilet paper in the toilet bowl.

2.


Defecate directly onto the toilet paper.

3.

Unscrew the blue cap of the tube.

4.

Dip the stick into 4 different zones before placing it back into the tube.

5.

Screw the blue dip-stick back into the tube.

6.

Write the patient ID on the tube.

7.

Keep the tube in the refrigerator. Do not freeze. Return to the laboratory as soon as possible.

NOTES FOR THE LABORATORY

- Before testing, shake vigorously the tube and wait at least 1 hour to help the feces sample dissolve in the buffered solution.
- Place the tube in the autoanalyzer after releasing the white cap.
- Verify that the tube does not collide with the sampling probe; if necessary transfer the sample to the instrument sample cup.
- If the tube seal is totally broken discard the sample.

LIMITATIONS

- Avoid sample collection during menstrual period, and in case of hemorrhoids or blood in urine.
- Blood, serum or other bodily fluids containing traces of blood processed with the same autoanalyzer may interfere with the assay.

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