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

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref:1001050

Cont.


R1: 20 x 2,5 mL, R2: 20 → 2,5 mL

Quantitative determination of creatine kinase (CK) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Creatine kinase (CK) catalyses the reversible transfer of a phosphate group from phosphocreatine to ADP. This reaction is coupled to those catalysed by hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6P-DH):

The rate of NADPH formation, measured photometrically, is proportional to the catalytic concentration of CK present in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Creatine kinase is a cellular enzyme with wide tissue distribution in the body. Its physiological role is associated with adenosine triphosphate (ATP) generation for contractile or transport systems.

Elevated CK values are observed in diseases of skeletal muscle and after myocardial infarction^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

| | | |
|-------------------------|--|------------|
| R 1 Buffer | Imidazol pH 7,0 | 100 mmol/L |
| | Glucose | 20 mmol/L |
| | Magnesium acetate | 10 mmol/L |
| | EDTA | 2 mmol/L |
| R 2 Substrate | ADP | 2 mmol/L |
| | AMP | 5 mmol/L |
| | di-Adenosine-5- pentaphosphate | 10 mmol/L |
| | NADP ⁺ | 2 mmol/L |
| | Hexokinase (HK) | 2500 U/L |
| | Glucose-6-phosphate dehydrogenase (G6P-DH) | 1500 U/L |
| | N-acetyl cysteine | 20 mmol/L |
| | Creatine phosphate | 30 mmol/L |

Optional

| | | |
|-------------------------------|-------------------------|--------------|
| CK-Nac / CK-MB CONTROL | Lyophilized human serum | Ref: 1002260 |
|-------------------------------|-------------------------|--------------|

PRECAUTIONS

R1: H360- May damage fertility or the unborn child.
Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Working reagent (WR):
Ref: 1001050. Dissolve (→) 1 tablet of R 2 Substrate with one vial of R 1.
Ref: 1001051. Dissolve (→) 1 tablet of R 2 Substrate with 15 mL of R 1.
Cap vial and mix gently to dissolve contents.
Stability: 5 days at 2-8°C or 24 hours at room temperature (15-25°C).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use the tablets if appears broken.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1,00.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C ó 37° C (± 0,1°C).
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C, protected from light.
The creatin kinase activity decreases 10% after 1 day at 2-5°C or after 1 hour at 15-25°C.

PROCEDURE

- Assay conditions:
Wavelength: 340 nm
Cuvette: 1 cm light path
Constant temperature: 25°C / 30°C / 37°C
- Adjust the instrument to zero with distilled water or air.
- Pipette into a cuvette:

| | | |
|-------------|-----------|------|
| | 25 - 30°C | 37°C |
| WR (mL) | 1,0 | 1,0 |
| Sample (µL) | 40 | 20 |

- Mix, incubate for 2 minutes.
- Read initial absorbance (A) of the sample, start the stopwatch and read absorbances at 1 minute intervals thereafter for 3 minutes.
- Calculate the difference between absorbances and the average absorbance differences per minute (ΔA/min).

CALCULATIONS

$$25^{\circ} - 30^{\circ}\text{C} \quad \Delta A / \text{min} \times 4127 = \text{U/L CK}$$

$$37^{\circ}\text{C} \quad \Delta A / \text{min} \times 8095 = \text{U/L CK}$$

Units: One international unit (IU) is the amount of enzyme that transforms 1 µmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

Temperature conversion factors

To correct results to other temperatures multiply by:

| Assay temperature | Conversion factor to | | |
|-------------------|----------------------|------|------|
| | 25°C | 30°C | 37°C |
| 25°C | 1,00 | 1,56 | 2,44 |
| 30°C | 0,64 | 1,00 | 1,56 |
| 37°C | 0,41 | 0,63 | 1,00 |

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

| | | | |
|--------------|--------|---------|---------|
| | 25°C | 30°C | 37°C |
| Men, up to | 80 U/L | 130 U/L | 195 U/L |
| Women, up to | 70 U/L | 110 U/L | 170 U/L |

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,878 U/L to linearity limit of 1300 U/L.
If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

| Mean (U/L) | Intra-assay (n=20) | | Inter-assay (n=20) | |
|------------|--------------------|------|--------------------|------|
| | 144 | 478 | 146 | 494 |
| SD | 3,88 | 6,98 | 4,55 | 9,57 |
| CV (%) | 2,69 | 1,49 | 3,11 | 1,94 |

Sensitivity: 1 U/L = 0,0001 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,997.

Regression equation: y = 1,031x - 0,5355.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

No interferences were observed with bilirubin up to < 20 mg/dL and hemoglobin up to 10 g/L^{1,2}. A list of drugs and other interfering substances with CK determination has been reported^{3,4}.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAPHIE

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
- Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: 1001050

Cont.

R1: 20 x 2,5 mL, R2: 20 → 2,5 mL