**CK-MB Control serum**
**IVD**

Store at 2 - 8°C

**PRODUCT CHARACTERISTICS**

CK-NAC / CK-MB CONTROL is a human lyophilised control serum.

Is intended for control of accuracy for use with manual and automated analytical procedures.

**REAGENTS**

Human serum. Biological additives. Bacteriostatics agents.

The concentration of the components are lot-specific.

See chart below for the values of ingredients.

**PRECAUTIONS**

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and HCV found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

**PREPARATION**

-Reconstitute (→) with 2.0 mL of distilled water.

-Mix thoroughly, avoiding foam forming.

-Bring to room temperature for about 30 min. before use.

Improper handling and/or storage can affect results.

Inaccurate reconstitution and errors in assay technique can cause erroneous results.

**STORAGE AND STABILITY**

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:

At 15°C to 25°C 8 hours

At 2°C to 8°C 5 days

At -25°C to -15°C 1 month

**PACKAGING**

REF	1002260	Cont.	4 x 2 mL
-----	---------	-------	----------

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REF	1002260	Cont.	4 x 2 mL
-----	---------	-------	----------

LOT		⌚	
-----	--	---	--

COMPONENTE COMPONENT	METODO METHOD	T <sup>a</sup>	VALOR VALUE	INTERVALO RANGE	UNIDADES UNITS
<b>CK-NAC</b>	DGKC	37°C			U / L
<b>CK-MB</b>	Método inmunológico Immunological method	37°C			U / L

\*This Control data sheet for lot is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.

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