CK-MB Control serum
IVD

Store at 2 - 8°C

PRODUCT CHARACTERISTICS

CK-Nac / CK-MB CONTROL is a human lyophilised control serum.

Is intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Biological additives. Bacteriostatics agents.

The concentration of the components are lot-specific.

See chart below for the values of ingredients.

PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and HCV found negative / non reactive by an FDA approved test. However, handle cautiously as potentially infectious.

PREPARATION

-Reconstitute (→) with 2.0 mL of distilled water.

-Mix thoroughly, avoiding foam forming.

-Bring to room temperature for about 30 min. before use.

Improper handling and/or storage can affect results.

Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth.

Store tightly capped when not use.

- After reconstitution is stable for:

At 15°C to 25°C 8 hours

At 2°C to 8°C 5 days

At -25°C to -15°C 1 month

PACKAGING

REF

1002262

Cont.

Level 1: 2 x 2 mL

Level 2: 2 x 2 mL

LOT

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REF

1002262

Cont.

Nivel 1: 2 x 2 mL

Nivel 2: 2 x 2 mL

	COMPONENTE COMPONENT	METODO METHOD	T ^a	VALOR VALUE	INTERVALO RANGE	UNIDADES UNITS
NIVEL 1 LEVEL1	CK-Nac	DGKC	37°C			U / L
	CK-MB	Método inmunológico Immunological method	37°C			U / L
NIVEL 2 LEVEL2	CK-Nac	DGKC	37°C			U / L
	CK-MB	Método inmunológico Immunological method	37°C			U / L

*This Control data sheet for lot is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.

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