ADA CONTROL

Adenosine Deaminase Control/
Control Adenosina Deaminasa

**Adenosine Deaminase (ADA) Control serum
IVD**

Store at 2-8°C

PRODUCT CHARACTERISTICS

The Adenosine Deaminase Control is used for quality control procedures in examining the accuracy of quantitative adenosine deaminase assays. The adenosine deaminase control is prepared in a bovine serum base, provided in lyophilized powder.

PRECAUTIONS

Product contains highly purified bovine source material from non-BSE-countries. The manufacturing facility does not receive, store or process ruminant materials from restricted countries.

Possible infectious agents in the materials have been inactivated. Because no method can offer complete assurance as to the absence of infectious agents, this material should be handled as thought capable of transmitting infectious disease and disposed as biohazard waste or medical waste according to applicable local and national laws.

PREPARATION

Open one vial carefully to avoid any loss of material and reconstitute with exactly 1 mL of distilled water. Close the vial and let stand for 30 min. at room temperature, dissolving contents completely by gently swirling or rotating.

STORAGE AND STABILITY

The lyophilized powder (not reconstituted) is stable until the expiration date indicated on the label when stored at 2-8°C. Once reconstituted, the components of the serum are stable for 1 week at 2-8°C.

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