Cholinesterase - Dibucaine


Ancillary reagent. Dibucaine inhibition

Quantitative determination of cholinesterase with dibucaine inhibition. IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Butyrylthiocholine kinetic method using dibucaine as inhibitor. Cholinesterase hydrolysed butyrylthiocholine to butyrate and thiocholine. Thiocholine reacts with 5,5'-dithiobis-2-nitrobenzoic acid (DTNB) to form 5-mercapto-2-nitrobenzoic acid (5-MNBA), according the following reactions:

The rate of 5-MNBA formation, measured photometrically, is proportional to the enzymatic activity of cholinesterase in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Cholinesterase is an enzyme present in plasma and synthesized by the liver. Its true physiological function is unknown, so its function may be to hydrolyze choline in plasma. Cholinesterase activity is usually measured for liver function, is a sensitive test of exposure to pesticides organophosphorus and identification of patients with the atypical form of enzyme whose presents high sensitivity to succinylcholine.

The dibucaine percent inhibition is determined in addition to total activity cholinesterase to detect the presence of an abnormal genetic variant of cholinesterase^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

CHOLINESTERASE Ref: 1001100		
R 1 Buffer	Phosphate pH 7.7	50 mmol/L
R 2 Substrate	5,5-dithiobis-2-nitrobenzoic acid (5,5 DTNB) Butyrylthiocholine	0.25 mmol/L 7 mmol/L

CHE-DIBUCAINE		
R 3 Inhibitor	Dibucain Chlorhidrate	2.6 mmol/L

PREPARATION

Working reagent (WR):
Dissolve (→) one tablet of R 2 Substrate in one vial of R 1.
Cap vial and mix gently to dissolve contents.
Inhibition (WR-I):
Add 105 µL of R 3 Inhibitor to working reagent (WR).
Stability: 2 hours at 2-8°C.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.
Do not use the tablets if appears broken.
Do not use reagents over the expiration date.
Signs of reagent deterioration:
- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm ≥ 0.70.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 405 nm.
- Thermostatic bath at 25°C, 30°C o 37°C (± 0.1°C)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment (Note 1).

SAMPLES

Serum or heparinized plasma¹: Stability 7 days at 2-8°C.

PROCEDURE

- Assay conditions:
Wavelength: 405 nm
Cuvette: 1 cm light path
Constant temperature 25°C /30°C / 37°C
- Adjust the instrument to zero with distilled water or air.

- Pipette into a cuvette:

	25 - 30°C	37°C
Without inhibitor WR (mL)	1.5	1.5
With inhibitor WR-I (mL)	1.5	1.5
Sample (µL)	10	--
Sample diluted 1/2 with NaCl 9 g/L (µL)	--	10

- Mix. Read initial absorbance (A) of the sample, start the stopwatch and read absorbances at 30 seconds intervals thereafter for 1.5 minutes.
- Calculate the difference between absorbances and the average absorbance differences per 30 seconds (ΔA/30 s).

CALCULATIONS

$$25^{\circ} - 30^{\circ}\text{C} \quad \Delta A / 30 \text{ s} \times 22537^* = \text{U/L}$$

$$37^{\circ}\text{C} \quad \Delta A / 30 \text{ s} \times 45074^* = \text{U/L}$$

* $\frac{\text{Tv} \times 1000}{\epsilon \times \text{LP} \times \text{Sv}}$	Tv= Total volumen in mL ε DTNB = 13.5 at 405 nm LP= Light path Sv= Sample volum in mL
--	--

$$\text{Percent inhibition (\%): } 100 - \frac{\text{Activity with inhibitor (U/L)}}{\text{Activity without inhibitor (U/L)}} \times 100$$

Units: One international unit (IU) is the amount of enzyme that transforms 1 µmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES^{1,2}

Percent inhibition of cholinesterase activity by dibucaine:

Normal 70-90% Heterozygotes 30-70% Homozygotes 0-20%

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: Up to linearity limit of 0.250 ΔA/30 s.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

Moderate hemolytic will not interfere in the results^{1,2}. A list of drugs and other interfering substances with cholinesterase determination has been reported by Young et. al^{3,4}.

NOTES

- Use glass containers only. Inhibitors may be extracted from some plastics.

BIBLIOGRAPHY

- King M. Cholinesterase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1108-1111.
- Whittaker M. et al. Comparison of a Commercially Available Assay System with Two Reference Methods for the Determination of Plasma Cholinesterase Variants. Clin. Chem 1983;(29/10); 1746-1760.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

- King M. Cholinesterase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1108-1111.
- Whittaker M. et al. Comparasion of a Commercially Available Assay System with Two Reference Methods for the Determination of Plasma Cholinesterase Variants. Clin. Chem 1983;(29/10); 1746-1760.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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