

Quantitative determination of Lipase IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The pancreatic lipase in presence of colipase, desoxycholate and calcium ions, hydrolyses the substrate 1-2-O-dilauryl-rac-glycero-3-glutaric acid-(6' -methylresorufin)-ester. The sequence of reactions involved in the enzymatic direct lipase determination is the following:

The rate of methylresorufin formation, measured photometrically, is proportional to the catalytic concentration of lipase present in the sample.

CLINICAL SIGNIFICANCE

Lipase (LPS) is a pancreatic enzyme necessary for the absorption and digestion of nutrients that catalyzes the hydrolysis of glycerol esters of fatty acids. Determination of LPS is used for diagnosis of diseases of pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct^{1,7,8}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 8.3	40 mmol/L
	Colipase	≥ 1 mg/L
	Desoxycholate	1,8 mmol/L
	Taurodesoxycholate	≥ 20 mmol/L
R 2 Substrate (micro-emulsion)	Tartrate pH 4,0	15 mmol/L
	Lipase Substrate	≥ 0,7 mmol/L
	Calcium chloride (CaCl ₂)	≥ 0,05 mmol/L

PREPARATION

The reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidimetry (note 1)
- Blank absorbance (A) at 570 nm ≥ 1,00.
- R2 is a turbid orange-colored micro-emulsion, discard if turning to red.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 570 nm.
- Thermostatic bath at 37° C (± 0,1°C)
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment (note 2)

SAMPLES

Fresh serum samples are recommended.

Spinreact does not recommend the use of EDTA and citrate plasma samples. These sample types could only be used if different reference ranges are established by each laboratory using different sample types other than fresh serum.

Avoid repeated freezing and thawing. Stability: 14 days at 2-8°C and at -20°C

PROCEDURE

- Assay conditions:
Wavelength: 570 nm
Cuvette: 1 cm light path
Constant temperature 37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette (Note 2):

	Blank	Standard / Sample
R 1 (mL)	1,0	1,0
R 2 (µL)	200	200
Distilled water (µL)	10	--
Standard / Sample (µL)	--	10

- Mix, incubate at 37°C for 1 minute.
- Read initial absorbance (A) of the sample, start the stopwatch and read absorbances at 1 minute intervals thereafter for 2 minutes.
- Calculate the difference between absorbances and the average absorbance differences per minute (ΔA/min).

CALCULATIONS

$$(\Delta A/\text{min})_{\text{Sample}} - (\Delta A/\text{min})_{\text{Blank}} = (\Delta A/\text{min})_{\text{of sample}}$$

$$(\Delta A/\text{min})_{\text{Calibrator}} - (\Delta A/\text{min})_{\text{Blank}} = (\Delta A/\text{min})_{\text{of Calibrator}}$$

$$\frac{\Delta A/\text{min Sample}}{\Delta A/\text{min Standard}} \times \text{Calibrator activity} = \text{Lipase activity of the sample (U/L)}$$

Units: One international unit (IU) is the amount of enzyme that transforms 1 µmol of substrate per minute, in standard conditions. The concentration is expressed in units per litre of sample (U/L).

Conversion factor: LPS [U/L] x 0,01667= LPS [µkal/L]

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

the normal range for pancreatic lipase when measured by the DGGR method is 13-60 U/L⁸

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

1. Trueness: Trueness has been performed using Certified Reference Material. Results obtained with Spinreact reagents show an acceptable proportional systematic deviation at low levels compared to the reference standard.

Theoretical conc. (U/L)	Measured conc. (U/L)	Bias (%)
10,04	10,03	-0,03
20,07	20,77	3,49
40,15	42,83	6,68
80,3	83,79	4,34
160,6	163,47	1,77

2. Precision: repeatability and reproducibility were evaluated based on the CLSI protocol EP05-A3:

Repeatability (n= 84)		Reproducibility (n= 90)	
Mean (U/L)	CV (%)	Mean (U/L)	CV (%)
46.8	3.3	49.25	5.0
93.46	2.4	103.2	4.3

3. Analytical sensitivity: the observed analytical sensitivity at 570nm was 0,378 mAbs/U/L.

4. Limit of detection (LoD): The Limit of Detection (LoD) for LIPASE-LQ. DGGR is 1,591 U/L, determined based on the CLSI EP17-A2 protocol.

5. Measuring range: From *quantification limit* (LoQ) of 4.32 U/L to *linearity limit* of 700 U/L. LoQ determined based on the CLSI EP17-A2 guidelines. And Linearity determined base on CLSI EP06-A 2nd Ed.

6. Method comparison: Results obtained using SPINREACT reagents (y) compared with other commercial reagents (x) did show slight proportional bias (less than 9%).

The results obtained using 50 samples were the following:

Correlation coefficient r= 0,994.

Ordinary least-squares fit: y= 1,313 + 1,076x

The results of the correlation may vary between analyzers and on the samples used.

INTERFERENCES

Bilirubin up to 75 mg/dL, hemoglobin up to 110 mg/dL and ascorbic acid up to 120 mg/dL do not interfere at the tested Lipase concentrations (40 U/L). 600 mg/dL Intralipid (1300 mg/dL triglycerides) do not interfere.

A list of drugs and other interfering substances with lipase determination has been reported^{2, 3}.

NOTES

- In some storage conditions (i.e. storage at a temperature lower that the one indicate) a precipitate may appear in the vial, it is not recommended to resuspend the precipitate as it could alter the functionality.
- In order to avoid contamination, it is recommended to use disposable material.

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