# Total and Direct Bilirubin

Jendrassik – Grof. Colorimetric

## Quantitative determination of bilirubin IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically.

Of the two fractions presents in serum, bilirubin-glucuronide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (bilirubin direct), while free bilirubin requires solubilization with caffeine to react (bilirubin indirect). In the determination of indirect bilirubin the direct is also determined, the results correspond to total bilirubin. The intensity of the color formed is proportional to the bilirubin concentration in the sample<sup>1,2,3</sup>.

### CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of hemoglobin.

It is transported from the spleen to the liver and excreted into bile.

Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma. Causes of hyperbilirubinemia:

Total bilirubin: Increased hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs.

Direct bilirubin: Hepatic cholestasis, genetic errors, hepatocellular damage<sup>1,6,7</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R 1</b>	Sulfanilic acid Hydrochloric acid (HCl)	30 mmol/L 400 mmol/L
<b>R 2</b>	Sodium nitrite	50 mmol/L
<b>R 3</b>	Caffeine	100 mmol/L
<b>Optional</b>	<b>BILIRUBIN CAL</b>	Ref: 1002250

### PRECAUTIONS

R1: H314-Causes severe burns and eye damage. EUH208-Contains sulfanilic acid. May produce an allergic reaction.

R3: H302.Harmful if swallowed. H371.May cause damage to organs.

Follow the precautionary statements given in MSDS and label of the product.

### PREPARATION

All the reagents are ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Color development in R 2.

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 540 nm.
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

### SAMPLES

Serum or plasma, free of hemolysis<sup>1</sup>.

Protect samples from direct light.

Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at -20°C.

### PROCEDURE

- Assay conditions:  
Wavelength: ..... 540 nm  
Cuvette: ..... 1 cm light path  
Temperature ..... 15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

	Total B.	Direct B.	Blank
R 1 (µL)	200	200	200
R 2 (drop)	1	1	--
NaCl 9 g/L (mL)	--	2,0	2,0
R 3 (mL)	2,0	--	--
Sample / Calibrator (µL) (Note 1)	200	200	200

- Mix and incubate for exactly **5 minutes** at 15-25°C.
- Measure the absorbance (A) of calibrator and sample.

### CALCULATIONS

- With Calibrator:

$$\frac{(A) \text{ Sample} - (A) \text{ Sample Blank}}{(A) \text{ Calibrator} - (A) \text{ Calibrator Blank}} \times \text{Conc. Calibrator} = \text{mg/dL bilirubin}$$

- **Theoretical calculation:** In the absence of a suitable calibrator, theoretical factors can be used, only as guidelines. The calculation with calibrator is the most appropriate.

Calculation: (A) Sample - (A) Blank Sample x Factor = mg / dL bilirubin

Total Bilirubin Factor: 13,8

Direct Bilirubin Factor: 10,4

**Conversion factor:** mg/dL x 17,1 = µmol/L.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES<sup>1</sup>

Total Bilirubin in adults Up to 1,10 mg/dL ≅ 18,81 µmol/L

Total Bilirubin in newborns <12 mg/dL ≅ < 205,2 µmol/L

Direct Bilirubin Up to 0,25 mg/dL ≅ 4,27 µmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of (T) 0,0926 mg/L (D) 0,0453 to *linearity limit* of 20 mg/dL. (T&D).

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

T Bilirubin	Intraserie (n= 20)	
Mean (mg/dL)	1,16	4,21
SD	0,02	0,04
CV (%)	2,03	1,06

Interserie (n= 20)	
1,15	4,27
0,02	0,13
1,91	3,10

D Bilirubin	Intra-assay (n=20)	
Mean (mg/dL)	0,78	2,28
SD	0,01	0,01
CV (%)	1,28	0,65

Inter-assay (n=20)	
0,80	2,18
0,01	0,03
1,63	1,53

**Sensitivity:** (T) 1 mg/dL = 0,07880 (A). (D) 1 mg/dL = 0,12772 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

#### TOTAL BILIRUBIN

Correlation coefficient (r)<sup>2</sup>: 0,9894

Regression equation: y= 0,9832x + 0,0224.

#### DIRECT BILIRUBIN

Correlation coefficient (r)<sup>2</sup>: 0,9867

Regression equation: y= 0,9923x + 0,0048.

The results of the performance characteristics depend on the analyzer used.

### INTERFERENCES

Hemolysis causes decreased bilirubin values<sup>1,2,3</sup>.

A list of drugs and other interfering substances with bilirubin has been reported by Young et. al<sup>4,5</sup>.

### NOTES

- For bilirubin determination in newborns, pipette 50 µL of sample. Multiply the result by 4.

### BIBLIOGRAPHY

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241, 436 and 650.
- Malloy H T et al. The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937; 112 (2): 481-491.
- Jendrassik L et al. Biochemische Zeitschrift Band 1938; 297:80-89.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241, 436 and 650.
- Malloy H T et al. The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937; 112 (2): 481-491.
- Jendrassik L et al. Biochemische Zeitschrift Band 1938; 297:80-89.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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