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### BIBLIOGRAFIA

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinez R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

### PRESENTACION

Ref: 1001043 Cont. R 1: 2 x 150 mL  
R 2: 1 x 10 mL

**Quantitative determination of bilirubin IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Bilirubin is converted to colored azobilirubin by diazotized sulfanilic acid and measured photometrically. Of the two fractions presents in serum, bilirubin-glucuronide and free bilirubin loosely bound to albumin, only the former reacts directly in aqueous solution (bilirubin direct), while free bilirubin requires solubilization with dimethylsulphoxide (DMSO) to react (bilirubin indirect). In the determination of indirect bilirubin the direct is also determined, the results correspond to total bilirubin. The intensity of the color formed is proportional to the bilirubin concentration in the sample<sup>1,2,3</sup>.

**CLINICAL SIGNIFICANCE**

Bilirubin is a breakdown product of hemoglobin. It is transported from the spleen to the liver and excreted into bile. Hyperbilirubinemia results from the increase of bilirubin concentrations in plasma. Causes of hyperbilirubinemia: Total bilirubin: Increase hemolysis, genetic errors, neonatal jaundice, ineffective erythropoiesis, and drugs. Direct bilirubin: Hepatic cholestasis, genetic errors, hepatocellular damage<sup>1,6,7</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	Sulfanilic acid (C <sub>6</sub> H <sub>7</sub> NO <sub>3</sub> S) Hydrochloric acid (HCl)	30 mmol/L 150 mmol/L
<b>R 2</b>	Sodium nitrite	29 mmol/L

**PRECAUTIONS**

R1: H314-Causes severe burns and eye damage. EUH208-Contains sulphanic acid (C<sub>6</sub>H<sub>7</sub>NO<sub>3</sub>S). May produce an allergic reaction. Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Color development in R 2.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or colorimeter measuring at 555 nm (530-580).
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

**SAMPLES**

Serum or plasma, free of hemolysis<sup>1</sup>. Protect samples from direct light. Stability: Bilirubin is stable at 2-8°C for 4 days and 2 months at -20°C.

**PROCEDURE**

- Assay conditions:  
Wavelength: ..... 555 nm (530-580)  
Cuvette: ..... 1 cm light path  
Temperature: ..... 15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette<sup>(Note 2)</sup>:

	Blank	B. Direct
R 1 (mL)	1,5	1,5
R 2 (µL)	--	50
Sample / Calibrator (µL) <sup>(Note 1)</sup>	100	100

- Mix and incubate for exactly **5 minutes** at 15-25°C.
- Read the absorbance (A).

**CALCULATIONS**

**- With Calibrator:**

$$\frac{(A) \text{ Sample} - (A) \text{ Sample Blank}}{(A) \text{ Calibrator} - (A) \text{ Calibrator Blank}} \times \text{Conc. Calibrator} = \text{mg/dL bilirubin}$$

**- With Factor:**

$$((A) \text{ Sample} - (A) \text{ Sample Blank}) \times \text{Factor}^* = \text{mg/dL bilirubin in the sample}$$

**\*Factor:**  $\frac{\text{Concentration of Calibrator}}{(A) \text{ Calibrator} - (A) \text{ Calibrator Blank}}$

**Conversion factor:** mg/dL x 17,1 = µmol/L.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**REFERENCE VALUES<sup>1</sup>**

Bilirubin Direct: Up to 0,25 mg/dL ≅ 4,27 µmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From *detection limit* of 0,07 mg/dL to *linearity limit* of 20 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	0,96	2,48	0,96	2,50
SD	0,024	0,051	0,043	0,035
CV (%)	2,52	2,06	4,49	1,41

**Sensitivity:** 1 mg/dL = 0,006856 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,96.

Regression equation: y=0,71177x -0,05267.

The results of the performance characteristics depend on the analyzer used.

**INTERFERENCES**

Hemolysis causes decreased bilirubin values<sup>1,2,3</sup>.

A list of drugs and other interfering substances with bilirubin has been reported<sup>4,5</sup>.

**NOTES**

- For bilirubin determination in newborns, pipette 50 µL of sample. Multiply the result by 2.
- Use clean disposable pipette for the dispensation.
- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

**BIBLIOGRAPHY**

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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### BIBLIOGRAPHIE

- Kaplan A et al. Bilirubin. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1238-1241. 436 and 650.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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