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- Franco R S. Hemoglobin. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1294-1296 and 418.
- Van Kampen EJ et al. Standardization of hemoglobinometry Clin. Chim 1961;6: 438-544.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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 Ref: 1001230 Cont. R: 4 x 5 mL
 Ref: 1001230S Cont. R: 4 x 5 mL, CAL: 1 x 1 mL

Quantitative determination of hemoglobin IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Hemoglobin is oxidized by potassium ferricyanide into metahemoglobin, which is converted into cyanometahemoglobin, by potassium cyanide.

 The intensity of the color formed is proportional to the hemoglobin concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

The hemoglobin is a protein that contains iron and that the red color to the blood. The hemoglobin is in red globules and it is the one in charge of oxygen transport by the blood from the lungs to weaves.

When the level of hemoglobin appears underneath the normal levels is describing an anemia that can be of different origins: primary anemia, cancer, pregnancy, renal diseases, and hemorrhages.

 If the hemoglobin levels appear high it can be due to: cardiopathies, dehydration and stays in places of much altitude^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

HEMOGLOBIN 50x	Potassium ferricyanide	0,60 mmol/L
	Potassium cyanide	77 mmol/L
	Dihydrogen potassium phosphate	2 mmol/L

Optional

HEMOGLOBIN CAL Ref.1001232	Hemoglobin Standard Animal origin	15 g/dL
----------------------------	-----------------------------------	---------

PRECAUTIONS

R: H302-Harmful if swallowed. H311+H331-Toxic in contact with skin or if inhaled. H412-Harmful to aquatic life with long lasting effects.

 CAL: H302-Harmful if swallowed. H315-Causes skin irritation. H319-Causes serious eye irritation. H331-Toxic if inhaled. H336-May cause drowsiness or dizziness. H351-Suspected of causing cancer. H361-Suspected of damaging fertility or the unborn child. H372-Causes damage to organs through prolonged or repeated exposure. Contains chloroform (CHCl₃). Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Working reagent (WR):

- For 5 mL 4,9 mL of distilled water + 2 drops of Reagent
- For 250 mL 245 mL of distilled water + 1 vial (5 mL) of Reagent

Mix well.

Stability: 2 months at 2-8°C, protected from the sunlight.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 540 nm \geq 0,012.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 540 nm.
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES

 Venous or capillary blood¹.

Use anticoagulants like EDTA, heparin or oxalate.

Stability of the sample: 1 week at 2-8°C.

PROCEDURE

- Assay conditions:
 Wavelength: 540 nm
 Cuvette: 1 cm. light path
 Temperature 15-25°C

2. Adjust the instrument to zero with distilled water.

3. Pipette into a cuvette:

A) MACRO METHOD:

	Blank	Standard	Sample
WR (mL)	5,0	5,0	5,0
Calibrator (μL)	--	20	--
Sample (μL)	--	--	20

B) MICRO METHOD:

	Blank	Standard	Sample
WR (mL)	2,5	2,5	2,5
Calibrator (μL)	--	10	--
Sample (μL)	--	--	10

4. Mix and incubate for 3 min. at room temperature (15-25°C).

5. Read the absorbance (A) of the samples and calibrator, against the Blank.

CALCULATIONS

 - With factor²:

$$(A) \text{ Sample} \times 36,77 = \text{g/dL hemoglobin in the sample}$$

- With calibrator:

$$\frac{(A) \text{ Sample} - (A) \text{ Blank}}{(A) \text{ Standard} - (A) \text{ Blank}} \times 15 (\text{Standard conc.}) = \text{g/dL hemoglobin in the sample}$$

QUALITY CONTROL

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES¹

 Men 14 - 18 g/dL \cong 8,7 - 11,2 mmol/L

 Women 12 - 16 g/dL \cong 7,5 - 9,9 mmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS
Measuring range: From detection limit of 0,108 g/dL to linearity limit of 20 g/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (g/dL)	SD	CV (%)	CV (%)
Mean (g/dL)	8,00	15,2	7,81	15,1
SD	0,29	0,33	0,19	0,26
CV (%)	3,59	2,19	2,51	1,74

Sensitivity: 1 g/dL = 0,027 (A).

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents. The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

 A list of drugs and other interfering substances with hemoglobin determination has been reported by Young et. al^{3,4}.

BIBLIOGRAPHY

- Franco R S. Hemoglobin. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1294-1296 and 418.
- Van Kampen EJ et al. Standardization of hemoglobinometry Clin. Chim 1961;6: 438-544.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
- Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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BIBLIOGRAFIA

1. Franco R S. Hemoglobin. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1294-1296 and 418.
2. Van Kampen EJ et al. Standardization of hemoglobinometry Clin. Chim 1961;6: 438-544.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
5. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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Ref: 1001230

Ref: 1001230S

Cont.

R: 4 x 5 mL

R: 4 x 5 mL, CAL: 1 x 1 mL