## Quantitative determination of potassium ion IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Potassium ions in a protein-free alkaline medium react with sodium tetraphenylboron to produce a finely dispersed turbid suspension of potassium tetraphenylboron.<sup>(1,2)</sup> The turbidity produced is proportional to the potassium concentration and read photometrically.

### CLINICAL SIGNIFICANCE

Potassium (K<sup>+</sup>) is the major positive ion within cells and is particularly important for maintaining the electric charge on the cell membrane. This charge allows nerves and muscles to communicate and is necessary for transporting nutrients into cells and waste products out of the cell. The concentration of potassium inside cells is about 30 times that in the blood and other fluids outside of cells. Potassium levels are mainly controlled by the steroid hormone aldosterone. Aldosterone is secreted from the adrenal gland when levels of potassium increase. Aldosterone, in turn, causes the body to rid itself of the excess potassium. Metabolic acidosis (for example, caused by uncontrolled diabetes) or alkalosis (for example, caused by excess vomiting) can affect blood potassium.

In normal people, taking potassium supplements or potassium-containing drugs is of no consequences, because the kidneys efficiently dispose of excess potassium.

### REAGENTS

<b>R1</b> TPB-Na	Sodium tetraphenylboron (TPB-Na)	0,2 mol/L
<b>R2</b> NaOH	Sodium hydroxide	2,0 mol/L
<b>R3</b> PREC	Trichloroacetic acid (TCA)	0,3 mol/L
<b>K-p CAL</b>	Potassium aqueous primary standard 5,0 mmol/L	

### PRECAUTIONS

R2: H314-Causes severe skin burns and eye damage.  
R3 /CAL: H314-Causes severe skin burns and eye damage. H335- May cause respiratory irritation. H411-Toxic to aquatic life with long lasting effects.

Follow the precautionary statements given in MSDS and label of the product.

### PREPARATION (Note 5)

Working reagent (WR):  
Mix equal volumes of R1 TPB-Na and R2 NaOH (Shake before to use).  
Don't use before 30 min. after its mixing. The working reagent must be shaken before each use.  
The working reagent is stable for 7 days at 15-25°C and 30 days at 2-8°C.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

**Do not use reagents over the expiration date.**

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 578 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment <sup>(Note 1, 2, 3)</sup>.

### SAMPLES

- Serum and lithium heparin plasma <sup>(Note 2)</sup>

### PROCEDURE

- Assay conditions:  
Wavelength: ..... 578 nm  
Cuvette: ..... 1 cm. light path  
Temperature ..... 37°C /15-25°C
- Adjust the instrument to zero with distilled water.
- Pipette into a tube <sup>(Note 3)</sup>:

Sample (µL)	50
R3-PREC (µL)	500

- Mix carefully. Centrifuge at high speed for 5-10 min.
- Separate the clear supernatant and pipette on another cuvette:

	Standard	Sample
Working reagent (mL) <sup>(Note 5)</sup>	1,0	1,0
Standard (µL) <sup>(Note 1,4)</sup>	100	--
Supernatant (µL)	--	100

To produce an homogeneous turbidity, the standard or the clear supernatant must be added to the center of the surface of the working reagent in the cuvette. Mix each cuvette carefully before proceeding to the next sample.

- Read the absorbance (A) of standard and samples against working reagent blank after 5 min. Color is stable up to 30 minutes.

### CALCULATIONS

$$\frac{(A) \text{ Sample} - (A) \text{ Blank}}{(A) \text{ Standard} - (A) \text{ Blank}} \times 5 (\text{Standard conc.}) = \text{mmol/L potassium in the sample}$$

**Conversion factor:** mmol/L = mEq/L.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).  
If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.  
Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

Serum:	3.60 – 5.50 mmol/L
Plasma:	4.00 – 4.80 mmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 2 mmol/L to linearity limit of 10 mmol/L.  
If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mmol/L)	SD	CV (%)	
Mean (mmol/L)	4,64	7,60	4,61	7,63
SD	0,095	0,10	0,113	0,148
CV (%)	2,05	1,32	2,45	1,94

**Sensitivity:** 1 mmol/L = 0.537A.

**Accuracy:** Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents.

Correlation coefficient: (r)<sup>2</sup> : 0,997

Regression equation: y = 0,988x + 0,489

The results of the performance characteristics depend on the analyzer used.

### INTERFERENCES

No significant interfering effect has been up to following concentrations: Bilirubin 40 mg/dL, Hemoglobin 450 mg/dL, Triglycerides 2500 mg/dL and Ascorbate 20 mg/dL.

### NOTES

- K-p CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
- As red blood cells contain about 25 times the amount of potassium, they have to be separated from the serum within one hour after blood collection. Otherwise, falsely elevated potassium concentrations will be found.
- Traces of detergents produce turbidity which leads to falsely elevated potassium concentrations. They therefore have to be avoided.
- The R2 (NaOH) and the working reagent must be shaken before their use.

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