### Quantitative determination of direct bilirubin IVD

Store at 2-8°C

#### PRINCIPLE OF THE METHOD

Direct bilirubin (conjugated) couples with the diazo reagent in the presence of sulfamic acid to form azobilirubin. The intensity of color formed is proportional to the bilirubin concentration in the sample tested. The increase of absorbance at 546 nm is directly proportional to the direct bilirubin concentration.

#### CLINICAL SIGNIFICANCE

Bilirubin is caused by the degradation of hemoglobin and exists in two forms. Unconjugated bilirubin is transported to the liver bound by albumin where it becomes conjugated (direct) with glucuronic acid and excreted. Hyperbilirubinemia is the result of an increase of bilirubin in plasma. Possible causes:

**Total bilirubin:** Increase hemolysis, genetic, neonatal jaundice, ineffective erythropoiesis and presence of drugs.

**Direct bilirubin:** Hepatic cholestasis, genetic, hepatocellular damage. Clinical diagnosis should not be made based on a single test result; it should integrate clinical and other laboratory data.

#### REAGENTS

R 1	Sulfamic acid	100 mM
R 2	2,4-DPD Hydrochloric acid (HCl)	0,5 mM 0,3 M

#### PRECAUTIONS

R1/R2: H314 Causes severe skin burns and eye damage. Follow the safety advice given in MSDS and product label.

#### PREPARATION

The reagents are provided in a ready to use format.

#### STORAGE AND STABILITY

The reagents are stable until the expiry date stated on the label when stored at 2-8°C, protected from light and contaminations are prevented during their use. Do not use reagents over the expiration date.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.

#### ADDITIONAL EQUIPMENT

- Spectrophotometer or analyzer capable of measuring absorbance at 546 nm.
- Cuvettes with 1,0 cm light path.
- General laboratory equipment.

#### SAMPLES

Serum or plasma, free of hemolysis. Protect samples from light. Stability of the sample: 4 days at 2-8°C or 2 month at -20°C.

#### PROCEDURE

- Assay conditions:  
Wavelength: .....546nm (530-580)  
Cuvette: .....1cm light path  
Temperature: .....37°C
- Adjust the instrument to zero with distilled water.
- Pipette into a cuvette:

	Calibrator blank	Sample blank
R 1 (μL)	800	800
Calibrator (μL)	50	-
Sample (μL)	-	50

- Mix and incubate for **5 minutes** at 37°C.
- Read the absorbance (A1) of the sample and calibrator.
- Add:

	Calibrator	Sample
R 2 (μL)	200	200

- Mix and incubate for **5 minutes** at 37°C.
- Read the absorbance (A2) of the sample and calibrator against the blank.
- Calculate the increase of the absorbance:  $\Delta A = A2 - A1$ .

#### CALCULATIONS:

- **With calibrator:**

$$\frac{(\Delta A)_{\text{Sample}}}{(\Delta A)_{\text{Calibrator}}} \times \text{Calibrator conc.} = \text{mg/dL of bilirubin in the sample}$$

- **With Factor:**  $(\Delta A)_{\text{Sample}} \times \text{Factor}^* = \text{mg/dL bilirubin in the sample}$

$$\text{*Factor: } \frac{\text{Calibrator concentration}}{(\Delta A)_{\text{Calibrator}}}$$

**Conversion factor:** mg/dL x 17,1 = μmol/L.

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210). If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

#### REFERENCE VALUES

Direct bilirubin 0 – 0,2 mg/dL (0 – 3,42 μmol/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

#### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 0,03 mg/dL to *linearity limit* of 9 mg/dL. If the results obtained are greater than the linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

	Inter assay (n= 40)		Intra assay (n= 80)	
	Mean (mg/dL)	SD	CV (%)	
Mean (mg/dL)	0,7458	2,444	0,7458	2,444
SD	0,05868	0,0550	0,0276	0,024
CV (%)	7,9	2,2	3,7	1,0

**Sensitivity:** 1 mg/dL = 0,040 Abs. units

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x) on a Spintech 240 analyzer. The results obtained using 53 samples ranging from 0,06 a 9 mg/dL (1,02 to 153,9 μmol/L) were:

Correlation coefficient (r): 0,9986

Regression equation:  $y = 1,0056 x - 0,1046$

The results of the performance characteristics depend on the analyzer used.

#### INTERFERENCES

Interferences from hemolysis, lipemia and ascorbic acid were evaluated for this direct bilirubin method on a Spintech 240 analyzer. Two concentrations of direct bilirubin were evaluated. No interferences were observed for lipemia (Intralipid) up to 350 mg/dL and ascorbic acid up to 40 mg/L. Hemolysis causes decreased direct bilirubin values, therefore hemolytic samples should be discarded.

A list of drugs and other interfering substances with bilirubin has been reported by Young et. al <sup>4,5</sup>.

#### NOTES

- SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

#### BIBLIOGRAPHY

- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2: 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966: Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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### BIBLIOGRAPHIE

- David G Levitt and Michael D Levitt. Quantitative assessment of the multiple processes responsible for bilirubin homeostasis in health and disease. Clin Exp Gastroenterol. 2014; 7: 307-328.
- Malloy H T. et al. The determination of bilirubin with the photoelectric colorimeter. J. Biol Chem 1937; 112, 2; 481-491.
- Martinek R. Improved micro-method for determination of serum bilirubin. Clin Chim 1966; Acta 13: 61-170.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Text book of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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