Quantitative determination of Fibrinogen IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Fibrinogen in presence of an excess of thrombin concentration, changes into Fibrin.

The time for clot formation in dilute plasma is inversely proportional to the fibrinogen concentration in the sample.

The thrombin clotting time fibrinogen assay is based on the method originally described by Clauss.¹ In the presence of high concentrations of thrombin, the time required for clot formation in dilute plasma is inversely proportional to the fibrinogen concentration.

CLINICAL SIGNIFICANCE

Fibrinogen (Factor I), protein synthesized by the liver, is the substance used in the blood to form a clot. Its determination is used to evaluate abnormal blood clotting.

Elevated Fibrinogen levels are observed in acute inflammations and in pregnancy; low values are observed in trombotic therapy, in hepatic disease, in the congenital non fibrinogen, in DIC (Disseminated Intravascular Coagulation) and in pancreatitis (low values)¹.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Bovine thrombin	≈ 100 NIH u /mL
R 2	Imidazole Buffer	
R3	Caolin Solution	
Optional	COAGULATION CAL	REF: 1709101
	CONTROL NORMAL	REF: 1709104
	CONTROL PATHOLOGIC	REF: 1709106

PREPARATION

R1: Dissolve (→) the vial content with 1.5 mL of distilled water and 0.5 mL of caolin solution, previously mixed. Cap and mix thoroughly avoiding foam forming to dissolve the total content, 15 minutes.. Stability: 7 days at 2-8°C or 1 month at -20°C, if stored in the original container and immediately frozen once reconstituted. Do not re-freeze once thawed.

Introduce the magnetic stirrer in the vial before to use the reagent in the analyzer.

R2: Ready for use. Mix thoroughly before use.

Both reagent positions are defined in the reagent's paragraph of Biobas 1000.

TRACEABILITY

The Coagulation Calibrator is traceable for Fibrinogen to the WHO standard, 2nd International Standard for Fibrinogen, plasma (98/612).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Quality control values outside established ranges.
- Product colour variations.

ADDITIONAL EQUIPMENT

- BIOBAS 1000.
- BIOBAS 1000 consumables.
- General laboratory equipment^(Note 1).

SAMPLES

Plasma from venous puncture diluted 1/10 in trisodium citrate solution 3.8% (105 mmol/L).

Mixing immediately the blood with anticoagulant. Avoid foaming the specimen. Centrifuge the sample at 3000 x g for 10 min and transfer the plasma to siliconized glass or plastic containers.

Turbid, icteric, lipemic or hemolyzed samples may generate erroneous results. The sample is stable for 4 hours at room temperature (15-25°C) or 28 days if immediately frozen at below 20°C.

NOTES

1. All labware must be clean and free of trace amounts of detergents.
2. Always follow instrument manufacturer's instructions; the results must be validated by the test laboratory.

PROCEDURE

The procedure to use is defined in the analyzer. Follow preparation and reagents assignment instructions.

CALIBRATION

Use the Spinreact coagulation calibrator, Ref. 1709101.

Apply to calibrator the samples dilution factor, multiplying the value of calibrator by 10. Use the value of calibrator calculated and assign it in the analyzer's calibration paragraph.

Introduce the values of the calibration curve indicated in the table:

Concentration (mg/dL)	600	400	300	200	100	0
-----------------------	-----	-----	-----	-----	-----	---

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures:

CONTROL NORMAL REF: 1709104

CONTROL PATHOLOGIC REF: 1709106

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

200 - 400 mg/dL¹ (2.00 – 4.00 g/L)

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Precision:

	Inter-assay (n=30)		
Mean (U/L)	144	294	488
CV (%)	5.9	3.4	2.9

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents.

INTERFERENCES AND LIMITATIONS

Has been observed interferences in samples with fibrinogen degradation. Acute inflammatory reactions can elevate circulating fibrinogen. Hemolysis can cause clotting factor activation and end point detection interference. High paraprotein levels, and drugs that activate the fibrinolytic system can interfere with fibrinogen assays. A list of drugs and other interfering substances with the determination has been reported^{2,3}.

BIBLIOGRAPHY

1. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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