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COIS05 13/06/14

SPINREACT, S.A./S.A.U. Ctra.Santa Coloma, 7 E-17176 SANT ESTEVE DE BAS (GI) SPAIN
Tel. +34 972 69 08 00 Fax +34 972 69 00 99 e-mail: spinreact@spinreact.com**COAGULATION CONTROL N****(CONTROL COAGULACION N)
CONTROL NORMAL / NORMAL CONTROL****Quantitative determination of coagulation factors****IVD**
Store at 2 - 8°C.**PRODUCT CHARACTERISTICS**

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REAGENTSHuman plasma collected with <0.4% sodium citrate anticoagulant, with a NORMAL concentration of coagulation factors.
Coagulation factors concentration is indicated below.**PRECAUTIONS**

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Actual results depend on many factors, including lot number, type of reagent, and instrument. Ranges must be determined in each laboratory with each change of lot number, reagent, or instrument.

This Control data sheet is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.

REF 1709104

4 x 1 mL

LOT

COIS05-F 03/03/16

SPINREACT,S.A./S.A.U. Ctra.Santa Coloma, 7 E-17176 SANT ESTEVE DE BAS (GI) SPAIN
Tél. +34 972 69 08 00 Fax +34 972 69 00 99 e-mail : spinreact@spinreact.com

COAGULATION CONTROL N

CONTROL NORMAL / NORMAL CONTROL
Coagulation / Coagulación**Quantitative determination of coagulation factors IVD**

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

The Control is a lyophilised human plasma, used to evaluate the precision and accuracy of PT, APTT and Fibrinogen determinations in human plasma.

REAGENTS

Human plasma collected with <0.4% sodium citrate anticoagulant, with a NORMAL concentration of coagulation factors.
Coagulation factors concentration is indicated below.

PRECAUTIONS

Each unit of source material used in the preparation of this product has been tested by an FDA licensed method and found non-reactive for HBsAg and negative for antibodies to HIV and HCV. However, as no known test method can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product must be handled as potentially infectious biological material.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Swirl gently (do not invert vial or mix vigorously) and let stand undisturbed for 15 minutes at room temperature before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date. After reconstitution, it's stable for 8 hours tightly closed at 2-8°C. Gently mix contents prior to each use.

Erratic values, product color variations, or lack of vacuum in the vials could indicate product deterioration. However, poor control performance could also be due to other factors within the test.

PROCEDURE

The Control should be run as a sample, at the initiation of testing each day and at least once each shift, or with each group of assays. Controls should also be tested with each reagent change or major instrument adjustment. Compare test results obtained to the expected results for the test method and control plasma.

COMPONENT / COMPOSANT	METHOD / MÉTHODE	RANGE / PLAGE
PT (Prothrombin Time / Temps de Prothrombine)	According to Instructions Sheet of corresponding Spinreact reagent/ Selon fiche d'instructions du réactif Spinreact correspondante	14.0 (11.9 – 16.1) s
APTT Activated Partial Thromboplastin Test Temps de Thromboplastine Partielle Activée		28.6 (24.3– 32.9) s
Fibrinogen / Fibrinogène		202 (172- 232) mg/dL

Actual results depend on many factors, including lot number, type of reagent, and instrument. Ranges must be determined in each laboratory with each change of lot number, reagent, or instrument.

This Control data sheet is applicable to sublots. Sequential alphabetical letter (e.g. A, B, C etc.) following the lot number.

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COAGULATION CONTROL N

(CONTROL COAGULAÇÃO P)
CONTROL NORMAL / NORMAL CONTROL**Quantitative determination of coagulation factors****IVD**

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

The Control is a lyophilised human plasma, used to evaluate the precision and accuracy of PT, APTT and Fibrinogen determinations in human plasma.

REAGENTS

Human plasma collected with <0.4% sodium citrate anticoagulant, with a NORMAL concentration of coagulation factors.
Coagulation factors concentration is indicated below.

PRECAUTIONS

Each unit of source material used in the preparation of this product has been tested by an FDA licensed method and found non-reactive for HBsAg and negative for antibodies to HIV and HCV. However, as no known test method can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product must be handled as potentially infectious biological material.

PREPARATION

Reconstitute (→) with 1.0 mL of distilled water. Swirl gently (do not invert vial or mix vigorously) and let stand undisturbed for 15 minutes at room temperature before use.

STORAGE AND STABILITY

The control is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date. After reconstitution, it's stable for 8 hours tightly closed at 2-8°C. Gently mix contents prior to each use.

Erratic values, product color variations, or lack of vacuum in the vials could indicate product deterioration. However, poor control performance could also be due to other factors within the test.

PROCEDURE

The Control should be run as a sample, at the initiation of testing each day and at least once each shift, or with each group of assays. Controls should also be tested with each reagent change or major instrument adjustment. Compare test results obtained to the expected results for the test method and control plasma.

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