Ecstasy Assay for Urine IVD

Store at 2 - 8°C.

INTENDED USE

The Ecstasy Enzyme Immunoassay (EIA) is a homogeneous enzyme immunoassay system intended for use in the qualitative and semi-quantitative analysis of MDMA and related compounds in human urine.

The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method (1, 2). Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.

PRINCIPLE OF THE METHOD

The ecstasy immunoassay is a homogeneous enzyme immunoassay (4) with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labelled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, ecstasy-labelled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to free drug; the unbound ecstasy-labelled G6PDH then exhibits its maximal enzyme activity.

Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

SIGNIFICANCE OF THE TEST

Ecstasy drugs are a group of amphetamine derivatives, including MDMA (3,4-methylenedioxyamphetamine), MDA (3,4-methylenedioxymphetamine), and MDEA (3,4-methylenedioxyethylamphetamine). They are central nervous system (CNS) stimulants. At light dose, ecstasy drugs produce euphoria, increase self-awareness. However, they are popularly abused for their psychotropic effects at high doses and become hallucinogenic and loss control of behaviour.

Toxic overdose causes depression, uncontrolled body fluid excretion, cardiac arrhythmias, and sleep disorder. Since there is no known medical application of ecstasy drugs with high abuse potential, US DEA list both MDMA and MDA as schedule I drugs. After ingestion of the drug, MDMA is known to metabolize to MDA by de-methylation.

Within the human body, most of the drug is eliminated through urinary excretion. Most of the urinary excretion is unchanged MDMA with a small fraction of MDA. Other urinary excretions include mono- and dihydroxy derivatives which appear as glucuronide conjugates (3). Detection of MDMA or its metabolites in urine indicates use of ecstasy.

REAGENTS

Antibody/Substrate Reagent (R₁): Contains monoclonal antibodies to ecstasy, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide as preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains ecstasy-labelled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide as preservative.

Avoid prolonged exposures of the reagent at temperatures higher than 25°C.

CALIBRATORS AND CONTROLS

Negative Human Urine (Level 0): Contains negative human urine with sodium azide as preservative (Ref.: 933010 o Ref.:933015).

Ecstasy Calibrator for Urine (Ref.: 933030):

Ecstasy Calibrator Level 0: Contains negative human urine with sodium azide as preservative.

Ecstasy Calibrator Level 1: Contains 100 ng/mL ecstasy in human urine with sodium azide as preservative.

Ecstasy Calibrator Level 2: Contains 500 ng/mL ecstasy in human urine with sodium azide as preservative.

Ecstasy Calibrator Level 3: Contains 750 ng/mL ecstasy in human urine with sodium azide as preservative.

Ecstasy Calibrator Level 4: Contains 1000 ng/mL ecstasy in human urine with sodium azide as preservative.

Ecstasy Calibrator (Cut-off 500 ng/mL): Contains 500 ng/mL ecstasy in human urine with sodium azide as preservative (Ref.: 932935).

Ecstasy Control: Contains 375 and 625 ng/mL ecstasy in human urine with sodium azide as preservative (Ref.: 935935).

PREPARATION

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

WARNING AND PRECAUTIONS

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

SPECIMEN COLLECTION

Urine samples may be collected in plastic or glass containers. Some plastics may adsorb drugs. Use a fresh urine specimen for the test. If the sample can not be analyzed immediately, it may be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

INSTRUMENTATION REQUIRED

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

PROCEDURE

Analyzers with above indicated specifications are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameter used for each analyzer before performing the assay. Typical assay parameters used for the analyzers include a sample to antibody reagent (R₁) to enzyme conjugate reagent (R₂) ratio of 1:10:3,75 respectively; a 37°C incubation temperature, 2-4 min reading frames, and 340nm primary wavelength.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

CALIBRATION

For qualitative determinations, the reagent must be calibrated with the selected cut-off calibrator. For semi-quantitative determinations, the reagent can be calibrated with a 5 point calibration curve.

Calibrator levels (ng/mL)		Control levels (ng/mL)
QUALITATIVE	SEMI-QUANTITATIVE	
Level 0 (0) Level 2 (500)	Level 0 (0)	Control -25% (375) Control +25% (625)
	Level 1 (100)	
	Level 2 (500)	
	Level 3 (750)	
	Level 4 (1000)	

The reagent should be recalibrated every month, when the controls are out of specifications (see Quality Control information), and when changing the reagent lot or the instrument settings.

INTERPRETATION

For qualitative determinations, the cut-off calibrator (500ng/mL) of ecstasy is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance (Δ mA/min) equal to, or greater than, that obtained with the cut-off calibrator is considered positive. A sample with a change in absorbance value lower than that obtained with the cut-off calibrator is considered negative.

For semi-quantitative determinations, a calibration curve with multiple calibrators is required. The concentration of the ecstasy in the sample may then be estimated from the calibration curve.

QUALITY CONTROL

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

The calibration curve can be validated with the controls level 375 and 625 ng/mL or with commercial controls.

LIMITATIONS

1. A positive result from the assay indicates only the presence of ecstasy.
2. Refer to the intended use statement for details of recommended confirmation testing methods.
3. The test is designed for use with human urine only.

PRECISION

According to the EP5 standards (CLSI), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80).

Qualitative analysis:

	CV (%)		
	Cut-off 500 ng/mL		
	375 ng/mL	500 ng/mL	625 ng/mL
Mean (mAU/min)	366.1	394.5	422.1
Total	1,3%	0,9%	0,9%
Within Run	0,3%	0,2%	0,2%
Between Run	0,9%	0,7%	0,5%
Between Day	0,8%	0,6%	0,7%

	Cut-off 500 ng/mL	
	0 ng/mL	1000 ng/mL
Mean (mAU/min)	200.2	452.3
Total	1%	0,8%
Within Run	0,3%	0,2%
Between Run	0,7%	0,4%
Between Day	0,7%	0,7%

Semi-quantitative analysis:

	CV (%)		
	Cut-off 500 ng/mL		
	375 ng/mL	500 ng/mL	625 ng/mL
Mean (ng/mL)	375.7	493.4	682.5
Total	2.8%	2.6%	4%
Within Run	0.8%	0.9%	1.1%
Between Run	2.2%	2.1%	2.8%
Between Day	1.4%	1.3%	2.7%

ACCURACY

52 urine samples from an external quality assurance scheme (UKNEQAS) were tested using GC-MS. The results from all participating laboratories were averaged to give the definitive concentration of drug in the sample. The same 52 samples were analysed using the Liquid EIA and the results compared.

3 samples screened positive by the Liquid EIA, and the same 3 samples were confirmed positive by GC-MS.

49 samples screened negative by the Liquid EIA, and the same 49 samples were confirmed negative by GC-MS.

There were no mismatching data sets.

SPECIFICITY

Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cut-off calibrator (500ng/mL).

It is possible that other substances and/or factors not listed below may interfere with the test and cause false positive results.

The following compounds were found to be non-cross reacting at 100,000ng/mL

11-hydroxy-delta9-THC	Diazepam
11-nor9-carboxy-delta9-THC	Dihydrocodeine
6-Acetyl-Morphine	Ecgonine Methyl Ester
Amitriptyline	EDDP
Amobarbital	EMDP
Amphetamine	Ephedrine
Aspirin	Heroin
Benzoyllecgonine	LAAM

β-phenylethylamine	Methadone
Cannabidiol	Methamphetamine
Chlorpheniramine	Morphine
Cocaeethylene	Oxycodone
Cocaine	Paracetamol
Codeine	Pseudoephedrine
Cotinine	Temazepam
Delta9-THC	

The following compounds will produce a positive response relative to the 500ng/mL cut-off at the following concentration:

Compound	Concentration (ng/mL)
MBDB	10,000
MDA	1,000
MDEA	1,000

BIBLIOGRAPHY

1. "Urine Testing for Drug of Abuse", National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. "Mandatory Guidelines for Federal Workplace Drug Testing Program", National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988)
3. Baselt RC and Cravey RH. "Disposition of Toxic Drugs and Chemicals in Man". Chemical Toxicology Institute, 4th ed. Foster City, CA., 1995
4. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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