Multidrug / *Multidroga* Low Calibrator Urine / *Calibrador Orina Bajo* DoA Liquid Test

Multidrug Calibrator Low Cut-off for Urine IVD

Store at 2 - 8°C.

INTENDED USE

The Multidrug calibrator is used for the calibration of the following Spinreact DoA-EIA qualitative liquid tests at low cutoff levels: Amphetamine, Barbiturates, Cocaine metabolite, Methadone, Opiates.

REAGENTS

Ref. 932910	Multidrug Calibrator Low Cut-off Human Urine. Sodium azide 0.5 g/L Methamphetamine, Secobarbital, Benzoyllecgonine, Methadone, Morphine.
-------------	--

The concentration of each drug in this calibrator is:

	ng/mL
Methamphetamine	300
Secobarbital	200
Benzoyllecgonine	150
Methadone	150
Morphine	300

PREPARATION AND STABILITY

The calibrator is ready to use. No calibrator preparation is required. The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

WARNING AND PRECAUTIONS

- The calibrator is for in vitro diagnostic use only. Harmful if swallowed.
- The calibrator contains sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such calibrator or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the calibrator beyond their expiration date.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix calibrators from different manufacturers.
- Do not freeze the calibrator.

Handle all urine specimens as if they were potentially infectious.

PROCEDURE

The Multidrug calibrator should be used in the following Urine reagents: Amphetamine, Barbiturates, Cocaine metabolite, Methadone, Opiates.

REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Spinreact DoA- EIA package inserts.

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REFERENCIAS

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Hoja de instrucciones de los kits Spinreact DoA-EIA.

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Ref. 932910 CONT. 1x10 mL