## Negative Calibrator / *Calibrador Negativo* DoA Liquid Test

### Negative Calibrator for Urine EIA IVD

Store at 2 - 8°C.

#### INTENDED USE

The EIA Negative Calibrator is for use as a calibrator in the EIA semiquantitative determination of drugs of abuse in human urine on automatic clinical chemistry analyzers.

#### REAGENTS

Refs. 933010, 933015	The calibrators contain Human Urine and preservative
-------------------------	--

#### PREPARATION AND STABILITY

The calibrator is ready to use. No calibrator preparation is required. The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

#### WARNING AND PRECAUTIONS

- This calibrator is for in vitro diagnostic use only. Harmful if swallowed.
- The calibrator contains sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such calibrators or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the calibrator beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix calibrators from different manufacturers.
- Do not freeze the calibrator.

Handle all urine specimens as if they were potentially infectious.

#### PROCEDURE

The calibrator is sold separately and may be used with any reagent lot.

#### REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.

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