BZO CAL 200 Urine

# Benzodiazepines / Benzodiazepinas

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DoA Liquid Test

### Benzodiazepine Calibrator Low Cut-off (200 ng/mL) for Urine IVD

Store at 2 - 8°C.

#### INTENDED USE

The Benzodiazepine Calibrator is used for the calibration of Spinreact DoA-EIA qualitative liquid test for Benzodiazepines at the low cutoff level (200ng/mL).

#### REAGENTS

Ref. 933081	Benzodiazepine Calibrator Low Cut-off (200 ng/mL) <b>Human Urine. Sodium azide 0.5 g/L</b> <b>Oxazepam</b>
-------------	--

The concentration of drug in this calibrator is:

	ng/mL
Oxazepam	200

#### PREPARATION AND STABILITY

The calibrator is ready to use. No calibrator preparation is required. The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

#### WARNING AND PRECAUTIONS

- The calibrator is for in vitro diagnostic use only. Harmful if swallowed.
- The calibrator contains sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such calibrator or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the calibrator beyond their expiration date.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix calibrators from different manufacturers.
- Do not freeze the calibrator.

Handle all urine specimens as if they were potentially infectious.

#### PROCEDURE

The calibrator should be used in the Benzodiazepine reagent for urine.

#### REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Spinreact Benzodiazepines EIA (Refs #930018/930020) package inserts.

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#### REFERENCIAS

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Hoja de instrucciones del kit Spinreact Benzodiazepinas-EIA (Refs # 930018/930020).

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Ref. 933081 **CONT.** 1X10 mL