Cannabinoids / *Cannabinoides* Multilevel Calibrator Urine / *Calibrador Multinivel Orina* DoA Liquid Test

Cannabinoids Calibrator for Urine IVD

Store at 2 - 8°C.

INTENDED USE

The Cannabinoids Multilevel Calibrator is used for the calibration of Spinreact DoA-EIA semi-quantitative liquid test for Cannabinoids.

REAGENTS

5 Levels Ref.: 933070	Cannabinoids Calibrator for Urine
Level 0	Negative Human Urine. Sodium azide 0.5 g/l.
Level 1	Cannabinoids Calibrator Level 1. Human Urine. Sodium azide 0.5 g/l. THC
Level 2	Cannabinoids Calibrator Level 2. Human Urine. Sodium azide 0.5 g/l. THC
Level 3	Cannabinoids Calibrator Level 3. Human Urine. Sodium azide 0.5 g/l. THC
Level 4	Cannabinoids Calibrator Level 4. Human Urine. Sodium azide 0.5 g/l. THC

Each level has the following concentrations:

ng/ml	Level 0	Level 1	Level 2	Level 3	Level 4
11-nor-Δ^9-THC-9-COOH	0	20	50	100	200

PREPARATION AND STABILITY

The calibrators are ready to use. No calibrator preparation is required. The calibrators are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

WARNING AND PRECAUTIONS

- The calibrators are for in vitro diagnostic use only. Harmful if swallowed.
- The calibrators contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such calibrators or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the calibrators beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix calibrators from different manufacturers.
- Do not freeze calibrators.

Handle all urine specimens as if they were potentially infectious.

PROCEDURE

The calibrators should be used in the Cannabinoids reagent for urine.

REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Spinreact THC EIA (Refs #930028/930030) package inserts.

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REFERENCIAS

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Hoja de instrucciones del kit Spinreact THC-EIA (Refs # 930028/930030).

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Ref. 933070 CONT. 5x4 mL