# Benzodiazepines / Benzodiazepinas

## Multilevel Calibrator Urine / *Calibrador Multinivel Orina*

DoA Liquid Test

### Benzodiazepine Calibrator for Urine IVD

Store at 2 - 8°C.

#### INTENDED USE

The Benzodiazepines Multilevel Calibrator is used for the calibration of Spinreact DoA-EIA semi-quantitative liquid test for Benzodiazepines.

#### REAGENTS

5 Levels Ref.: 933080	Benzodiazepine Calibrator for Urine
Level 0	<b>Negative Human Urine. Sodium azide 0.5 g/l.</b>
Level 1	<b>Benzodiazepine Calibrator Level 1. Human Urine. Sodium azide 0.5 g/L. Oxazepam</b>
Level 2	<b>Benzodiazepine Calibrator Level 2. Human Urine. Sodium azide 0.5 g/L. Oxazepam</b>
Level 3	<b>Benzodiazepine Calibrator Level 3. Human Urine. Sodium azide 0.5 g/L. Oxazepam</b>
Level 4	<b>Benzodiazepine Calibrator Level 4. Human Urine. Sodium azide 0.5 g/L. Oxazepam</b>

Each level has the following concentrations:

ng/ml	Level 0	Level 1	Level 2	Level 3	Level 4
Oxazepam	0	100	200	300	1000

#### PREPARATION AND STABILITY

The calibrators are ready to use. No calibrator preparation is required. The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

#### WARNING AND PRECAUTIONS

- The calibrators are for in vitro diagnostic use only. Harmful if swallowed.
- The calibrators contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such calibrators or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the calibrators beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix calibrators from different manufacturers.
- Do not freeze calibrators.

Handle all urine specimens as if they are potentially infectious.

#### PROCEDURE

The calibrators should be used in the Benzodiazepine reagent for urine.

#### REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Spinreact Benzodiazepines EIA (Refs #930018/930020) package inserts.

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#### REFERENCIAS

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Hoja de instrucciones del kit Spinreact Benzodiazepinas-EIA (Refs # 930018/930020).

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Ref. 933080

<b>CONT.</b>
--------------

5x10 mL