## Cannabinoids / *Cannabinoides* High Cut-off Urine Control 50 / *Control Cut-off Alto Orina 50* DoA Liquid Test

### Cannabinoids Control High Cut-off (50 ng/mL) for Urine IVD

Store at 2 - 8°C.

#### INTENDED USE

The Cannabinoids Control High Cutoff is consisting of two levels  $\pm 25\%$  from the high cutoff (50ng/mL) for the control of Spinreact DoA-EIA liquid test for Cannabinoids.

#### REAGENTS

|             |  |
|-------------|--|
| Level - 25% | Cannabinoids Control High Cut-off (50 ng/mL)<br><b>Human Urine. Sodium azide 0.5 g/L THC</b> |
| Level + 25% | Cannabinoids Control High Cut-off (50 ng/mL)<br><b>Human Urine. Sodium azide 0.5 g/L THC</b> |

Each level has the following concentrations:

| ng/mL                          | - 25 % | + 25 % |
|--------------------------------|--------|--------|
| 11-nor- $\Delta^9$ -THC-9-COOH | 37,5   | 62,5   |

#### PREPARATION AND STABILITY

The controls are ready to use. No control preparation is required. The controls are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

#### WARNING AND PRECAUTIONS

- The controls are for in vitro diagnostic use only. Harmful if swallowed.
- The controls contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such controls or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the controls beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix controls from different manufacturers.
- Do not freeze controls.

Handle all urine specimens as if they were potentially infectious.

#### PROCEDURE

The controls should be used in the Cannabinoids reagent for urine.

#### REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Spinreact THC EIA (Refs #930028/930030) package inserts.

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#### REFERENCIAS

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTS, 2002.
3. Hoja de instrucciones del kit Spinreact THC-EIA (Refs # 930028/930030).

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|             |              |              |        |             |        |
|-------------|--------------|--------------|--------|-------------|--------|
| Ref. 936072 | <b>CONT.</b> | Nivel - 25%: | 2x4 mL | Nivel +25%: | 2x4 mL |
|-------------|--------------|--------------|--------|-------------|--------|