## Barbiturates Assay for Urine IVD

Store at 2 - 8°C

### INTENDED USE

The Barbiturates Enzyme Immunoassay (EIA) is a homogenous enzyme immunoassay system intended for use in the qualitative and semi-quantitative analysis of barbiturates in human urine.

**The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method (1, 2). Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.**

### PRINCIPLE OF THE METHOD

The barbiturates immunoassay is a homogeneous enzyme immunoassay (3) with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labelled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, barbiturates-labelled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to the free drug; the unbound barbiturates-labelled G6PDH then exhibits its maximal enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

### SIGNIFICANCE OF THE TEST

Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of Barbiturates leads to tolerance and physical dependence. Short acting Barbiturates taken at 400 mg/day for 2-3 months produces a clinically significant degree of physical dependence.

Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death. Only a small amount (less than 5%) of most Barbiturates are excreted unaltered in the urine. The detection period for the Barbiturates in the urine is 4-7 days.

### REAGENTS

**Antibody/Substrate Reagent (R<sub>1</sub>):** Contains monoclonal antibodies to secobarbital, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide as preservative.

**Enzyme-drug Conjugate Reagent (R<sub>2</sub>):** Contains secobarbital-labelled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide as preservative.

**Avoid prolonged exposures of the reagent at temperatures higher than 25°C.**

### CALBRATORS AND CONTROLS

**Negative Human Urine (Level 0):** Contains negative human urine with sodium azide as preservative (Ref.: 933010 and 933015).

**Multidrug Calibrator Urine (Ref.: 933020):**

**Multidrug Calibrator Level 0:** Contains negative human urine with sodium azide as preservative.

**Multidrug Calibrator Level 1:** Contains 100ng/mL secobarbital in human urine with sodium azide as preservative.

**Multidrug Calibrator Level 2:** Contains 200ng/mL secobarbital in human urine with sodium azide as preservative.

**Multidrug Calibrator Level 3:** Contains 1000ng/mL secobarbital in human urine with sodium azide as preservative.

**Multidrug Calibrator Low Cut-Off:** Contains 200 ng/mL Secobarbital in human urine with sodium azide as preservative (Ref.: 932910).

**Multidrug Calibrator High Cut-Off:** Contains 200 ng/mL Secobarbital in human urine with sodium azide as preservative (Ref.: 932915).

**Multidrug Control Low Cut-Off:** Contains 150 and 250 ng/mL Secobarbital in human urine with sodium azide as preservative (Ref.: 935910).

**Multidrug Control High Cut-Off:** Contains 150 and 250 ng/mL Secobarbital in human urine with sodium azide as preservative (Ref.: 935915).

### PREPARATION

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

### WARNING AND PRECAUTIONS

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

### SPECIMEN COLLECTION

Urine samples may be collected in plastic or glass containers. Some plastics may adsorb drugs. Use a fresh urine specimen for the test. If the sample can not be analyzed immediately, it may be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

### INSTRUMENTATION REQUIRED

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

### PROCEDURE

Analyzers with above indicated specifications are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameter used for each analyzer before performing the assay. Typical assay parameters used for the analyzers include a sample to antibody reagent (R<sub>1</sub>) to enzyme conjugate reagent (R<sub>2</sub>) ratio of 1:10:3,75 respectively; a 37°C incubation temperature, 2-4 min reading frames, and 340nm primary wavelength.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

### CALIBRATION

For qualitative determinations, the reagent must be calibrated with the selected cut-off calibrator. For semi-quantitative determinations, the reagent can be calibrated with a 4 point calibration curve.

Calibrator levels (ng/mL)		Control levels (ng/mL)
QUALITATIVE	SEMI-QUANTITATIVE	
Level 0 (0)	Level 0 (0)	Control - 25% (150) Control + 25% (250)
Level 2 (200)	Level 1 (100)	
	Level 2 (200)	
	Level 3 (1000)	

Recalibrate each month, when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

### INTERPRETATION

For qualitative determinations, the cut-off calibrator (200ng/mL) of secobarbital is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance ( $\Delta$ mA/min) equal to, or greater than, that obtained with the cut-off calibrator is considered positive. A sample with a change in absorbance value lower than that obtained with the cut-off calibrator is considered negative.

For semi-quantitative determinations, a calibration curve with multiple calibrators is required. The concentration of secobarbital in the sample may then be estimated from the calibration curve.

### QUALITY CONTROL

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

The calibration curve can be validated with the Control levels 150 and 250 ng/mL (ref. 935910 and ref. 935915), or with commercial controls.

**LIMITATIONS**

1. A positive result from the assay indicates only the presence of secobarbital.
2. Refer to the intended use statement for details of recommended confirmation testing methods.
3. The test is designed for use with human urine only.

**PRECISION**

According to the EP5 standards (NCCLS), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80).

**Qualitative analysis:**

	CV (%)		
	Cut-off 200 ng/mL		
	150 ng/mL	200 ng/mL	250 ng/mL
<b>Mean (mAU/min)</b>	375.9	401.1	403.9
<b>Total</b>	1.1%	1.1%	1%
<b>Within Run</b>	0.2%	0.2%	0.2%
<b>Between Run</b>	0.6%	0.6%	0.6%
<b>Between Day</b>	0.9%	0.9%	0.8%

	Cut-off 200 ng/mL	
	0 ng/mL	1000 ng/mL
<b>Mean (mAU/min)</b>	332.8	463.1
<b>Total</b>	1.3%	1%
<b>Within Run</b>	0.2%	0.2%
<b>Between Run</b>	0.8%	0.7%
<b>Between Day</b>	1%	0.8%

**Semi-quantitative analysis:**

	CV (%)		
	Cut-off 200 ng/mL		
	150 ng/mL	200 ng/mL	250 ng/mL
<b>Mean (ng/mL)</b>	127.8	203.1	216
<b>Total</b>	2.4%	3.5%	3%
<b>Within Run</b>	1.5%	2.1%	1.8%
<b>Between Run</b>	1.8%	2.5%	2.3%
<b>Between Day</b>	0.6%	1.2%	0.8%

**SPECIFICITY**

Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cut-off calibrator (200ng/mL).

It is possible that other substances and/or factors not listed below may interfere with the test and cause false positive results.

The following compounds were found to be non-cross reacting at 100,000ng/mL:

11-hydroxy-delta9-THC	EDDP
11-nor9-carboxy-delta9-THC	EMDP
6 Acetyl Morphine	Ephedrine
Amphetamine	Heroin
Aspirin	LAAM
β-phenylethylamine	MBDB
Benzoyllecgonine	MDA
Cannabidiol	MDEA
Chlorpheniramine	MDMA
Cocaethylene	Methadone
Cocaine	Methamphetamine
Codeine	Morphine
Cotinine	Oxycodone
Delta9-THC	Paracetamol
Diazepam	Pseudoephedrine
Dihydrocodeine	Temazepam
Ecgonine methyl ester	

The following barbiturates will produce a positive response relative to the 200ng/mL cut-off at the following concentrations:

Compound	Concentration(ng/mL)
Amobarbital	1200
Butalbital	500
Hexobarbital	25.000
Phenobarbital	500
Pentobarbital	900

**BIBLIOGRAPHY**

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988).
3. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, Biochem Biophys Res Commun, 47, 846 (1972)

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	CV (%)		
	Cut-off 200 ng/mL		
	150 ng/mL	200 ng/mL	250 ng/mL
<b>Media (mAU/min)</b>	375.9	401.1	403.9
<b>Total</b>	1.1%	1.1%	1%
<b>Within Run</b>	0.2%	0.2%	0.2%
<b>Between Run</b>	0.6%	0.6%	0.6%
<b>Between Day</b>	0.9%	0.9%	0.8%

	Cut-off 200 ng/mL	
	0 ng/mL	1000 ng/mL
<b>Media (mAU/min)</b>	332.8	463.1
<b>Total</b>	1.3%	1%
<b>Within Run</b>	0.2%	0.2%
<b>Between Run</b>	0.8%	0.7%
<b>Between Day</b>	1%	0.8%

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	CV (%)		
	Cut-off 200 ng/mL		
	150 ng/mL	200 ng/mL	250 ng/mL
<b>Media (ng/mL)</b>	127.8	203.1	216
<b>Total</b>	2.4%	3.5%	3%
<b>Within Run</b>	1.5%	2.1%	1.8%
<b>Between Run</b>	1.8%	2.5%	2.3%
<b>Between Day</b>	0.6%	1.2%	0.8%

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1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988)
3. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, Biochem Biophys Res Commun, 47, 846 (1972)

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Ref.: 930098 200 tests	Cont.	: 2 x 21 mL R1 : 2 x 8 mL R2
Ref.: 930100 500 tests	Cont.	: 1 x 105 mL R1 : 1 x 38 mL R2
Ref.: 930102 2500 tests	Cont.	: 2 x 250 mL R1 : 1 x 190 mL R2