Benzodiazepines Urine

DoA Liquid Test

Benzodiazepines Assay for Urine IVD

Store at 2 - 8°C.

INTENDED USE

The Benzodiazepine Enzyme Immunoassay (EIA) is a homogeneous enzyme immunoassay system intended for use in the qualitative and semi-quantitative analysis of benzodiazepines in human urine.

The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method (1, 2). Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.

PRINCIPLE OF THE METHOD

The benzodiazepine assay is a homogeneous enzyme immunoassay (3) with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labelled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, benzodiazepine-labelled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to free drug; the unbound benzodiazepine-labelled G6PDH then exhibits its maximal enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

SIGNIFICANCE OF THE TEST

The benzodiazepine family has a diverse chemical structure which gives rise to a range of physio-chemical and pharmacological properties. They form a large group and are used as sedatives and in the treatment of anxiety. It is estimated that over 50 benzodiazepine drugs are in common use world-wide with approximately a dozen marketed within a give country. Benzodiazepines are strongly associated with opioid users (4).

They are usually taken orally (although intravenous and intramuscular routes are also used) and as would be expected from such a diverse group, the pharmacological effects are broad and are also related with the absorption rate. In general, they are extensively metabolised by dealkylation, hydroxylation, oxidation and glucuronidation in the liver.

Some metabolites have biological activity; for example nordiazepam from diazepam and oxazepam from temazepam. Detection of benzodiazepines and their metabolites in urine can be used as an indication for use of benzodiazepines.

REAGENTS

Antibody/Substrate Reagent (R₁): Contains monoclonal antibodies benzodiazepines, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide as preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains benzodiazepine-labelled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide as preservative.

Avoid prolonged exposures of the reagent at temperatures higher than 25°C.

CALIBRATORS AND CONTROLS

Negative Human Urine (Level 0): Contains negative human urine with sodium azide as preservative (Ref.: 933010 and 933015).

Benzodiazepine Calibrator (Ref.: 933080)

Benzodiazepine Calibrator Level 0: Contains negative human urine with sodium azide as preservative.

Benzodiazepine Calibrator Level 1: Contains 100 ng/mL benzodiazepine in human urine with sodium azide as preservative.

Benzodiazepine Calibrator Level 2: Contains 200 ng/mL benzodiazepine in human urine with sodium azide as preservative.

Benzodiazepine Calibrator Level 3: Contains 300 ng/mL benzodiazepine in human urine with sodium azide as preservative.

Benzodiazepine Calibrator Level 4: Contains 1000 ng/mL benzodiazepine in human urine with sodium azide as preservative.

Benzodiazepine Calibrator Low Cut-Off 200: Contains 200 ng/mL Benzodiazepine in human urine with sodium azide as preservative. (Ref.: 933081)

Benzodiazepine Calibrator High Cut-Off 300: Contains 300 ng/mL Benzodiazepine in human urine with sodium azide as preservative. (Ref.: 933082)

Benzodiazepine Control Low Cut-Off 200: Contains 150 and 250 ng/mL Benzodiazepine in human urine with sodium azide as preservative. (Ref.: 936081).

Benzodiazepine Control High Cut-Off 300: Contains 225 and 375 ng/mL Benzodiazepine in human urine with sodium azide as preservative (Ref.: 936082).

PREPARATION

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

WARNING AND PRECAUTIONS

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

SPECIMEN COLLECTION

Urine samples may be collected in plastic or glass containers. Benzodiazepines are very thermolabile. If the sample cannot be analyzed immediately, they should be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

INSTRUMENTATION REQUIRED

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

PROCEDURE

Analyzers with above indicated specifications are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameter used for each analyzer before performing the assay. Typical assay parameters used for the analyzers include a sample to antibody reagent (R₁) to enzyme conjugate reagent (R₂) ratio of 1:10:3,75 respectively; a 37°C incubation temperature, 2-4 min reading frames, and 340nm primary wavelength.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

CALIBRATION

For qualitative determinations, the reagent must be calibrated with the selected cut-off calibrator. For semi-quantitative determinations, the reagent can be calibrated with a 5 point calibration curve.

CUT-OFF	Calibrator levels (ng/mL)		Control levels (ng/mL)
	QUALITATIVE	SEMI-QUANTITATIVE	
LOW 200 ng/mL	Level 0 (0) Level 2 (200)	Level 0 (0) Level 1 (100) Level 2 (200) Level 3 (300) Level 4 (1000)	Control - 25% (150) Control + 25% (250)
HIGH 300 ng/mL	Level 0 (0) Level 3 (300)		Control - 25% (225) Control + 25% (375)

Recalibrate each month, when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

Benzodiazepines Urine

DoA Liquid Test

INTERPRETATION

For qualitative determinations, the cut-off calibrator (200 or 300ng/mL) of benzodiazepine is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance ($\Delta A/\text{min}$) equal to, or greater than, that obtained with the cut-off calibrator is considered positive. A sample with a change in absorbance value lower than that obtained with the cut-off calibrator is considered negative.

For semi-quantitative determinations, a calibration curve with multiple calibrators is required. The concentration of benzodiazepines in the sample may then be estimated from the calibration curve.

QUALITY CONTROL

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

The calibration curve can be validated with the Control levels 150 and 250 ng/mL (ref. 936081) or levels 225 and 375 ng/mL (ref. 936082), or with commercial controls.

LIMITATIONS

1. A positive result from the assay indicates only the presence of benzodiazepines.
2. Positive results should be confirmed by other affirmative, analytical methods (e.g. chromatography), and preferably GC/MS.
3. The test is designed for use with human urine only.

PRECISION AND SENSITIVITY

Details of the precision and sensitivity studies are available on request.

ACCURACY

52 urine samples from an external quality assurance scheme (UKNEQAS) were tested using GC-MS. The results from all participating laboratories were averaged to give the definitive concentration of drug in the sample. The same 52 samples were analysed using the Liquid EIA and the results compared.

13 samples screened positive by the Liquid EIA, and the same 13 samples were confirmed positive by GC-MS.

38 samples screened negative by the Liquid EIA, and the same 38 samples were confirmed negative by GC-MS.

One sample was positive in the screen, but negative when confirmed. The donor in question was taking clonazepam which was not looked for by the confirmation method. This would explain the discrepancy in results.

SPECIFICITY

Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cut-off calibrator (300ng/mL).

It is possible that other substance and/or factors not listed below may interfere with the test and cause false positive results.

The following compounds were found to be non-cross reacting at 100,000ng/mL:

11-hydroxy-delta9-THC	EDDP
11-nor9-carboxy-delta9-THC	EMDP
6-Acetyl-Morphine	Ephedrine
Amitriptyline	Heroin
Amobarbital	LAAM
Amphetamine	MBDB
Aspirin	MDA
Benzoylcegonine	MDEA
β -phenylethylamine	MDMA
Cannabidiol	Metadone
Cocaine	Methamphetamine
Codeine	Morphine
Cotinine	Benzodiazepine
Delta9-THC	Paracetamol
Dihydrocodeine	Pseudoephedrine
Ecgonine Methyl Ester	

The following benzodiazepines will produce a positive relative to the 300ng/mL cut-off at the following concentrations:

Compound	Concentration (ng/mL)
Alprazolam	300
Bromazepam	2,000
Chlordiazepoxide	300
Clobazam	2,000
Clonazepam	300
Diazepam	300
Flunitrazepam	300
Flurazepam	300
Lormetazepam	300
Lorazepam	300
Medazepam	300
Midazolam	400
Nitrazepam	300
Prezepam	300
Temazepam	300
Triazolam	300

BIBLIOGRAPHY

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988).
3. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, Biochem Biophys Res Commun, 47, 846 (1972).
4. The Forensic Pharmacology of Drugs of Abuse, Olaf Drummer. Arnold publication ISBN 034076257

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