Cannabinoid Assay for Urine IVD

Store at 2 - 8°C

INTENDED USE

The Cannabinoid Enzyme Immunoassay (EIA) is a homogenous enzyme immunoassay system intended for use in the qualitative and semi-quantitative analysis of cannabinoids in human urine.

The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method (1, 2). Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.

PRINCIPLE OF THE METHOD

The cannabinoids immunoassay is a homogeneous enzyme immunoassay (8) with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labelled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, Δ^9 -THC-labelled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to the free drug; the unbound Δ^9 -THC-labelled G6PDH then exhibits its maximal enzyme activity.

Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

SIGNIFICANCE OF THE TEST

The principal, active constituent in marijuana (or hashish), obtained from *Cannabis sativa* plant, is the Δ^1 -3,4-*trans* tetrahydrocannabinol, frequently referred as Δ^9 -tetrahydrocannabinol or Δ^9 -THC. Cannabis has been used for its euphoric effects for over 4000 years (3). It is one of the most commonly used illicit drugs in the United States.

Marijuana is frequently self-administered for its mood-altering properties. Chronic use has been shown to cause reversible psychological impairment, an abstinence syndrome, and can cause users to develop tolerance (4). At low doses, it produces mixed depressant and stimulant effects; while at higher dose it acts as a CNS depressant (5,6,7).

Δ^9 -THC is easily absorbed by inhalation (smoking) or ingestion through the gastrointestinal tract. Due to its highly fat-soluble nature, Δ^9 -THC is readily deposited in fatty tissues, where it may remain for days or even weeks (5). It is primarily metabolized in the liver to a variety of compounds, the primary one being the 11-nor- Δ^9 -tetrahydrocannabinol-9-carboxylic acid (11-nor- Δ^9 -THC-9-COOH) (6, 7). Approximately 70% of THC is excreted in faeces and urine within 72 hours of administration.

REAGENTS

Antibody/Substrate Reagent (R₁): Contains monoclonal antibodies to Δ^9 -THC, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide as preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains Δ^9 -THC-labelled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide as preservative.

Avoid prolonged exposures of the reagent at temperatures higher than 25°C.

CALIBRATORS AND CONTROLS

Negative Human Urine (Level 0): Contains negative human urine with sodium azide as preservative (Ref.: 933010 or Ref.:933015).

THC Calibrator (Ref. 933070):

THC Calibrator Level 0: Contains human urine with sodium azide as preservative.

THC Calibrator Level 1: Contains 20 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative.

THC Calibrator Level 2: Contains 50 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative.

THC Calibrator Level 3: Contains 100 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative.

C THC Calibrator Level 4: Contains 200 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative.

THC Calibrator Low Cut-Off 20: Contains 20 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative (Ref.: 933071).

THC Calibrator High Cut-Off 50: Contains 50 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative (Ref.: 933072).

THC Control Low Cut-Off 20: Contains 15 and 25 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative (Ref.: 936071).

THC Control High Cut-Off 50: Contains 37,5 and 62,5 ng/mL 11-nor- Δ^9 -THC-9-COOH in human urine with sodium azide as preservative (Ref.: 936072).

PREPARATION

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

WARNING AND PRECAUTIONS

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

SPECIMEN COLLECTION

Urine samples may be collected in plastic or glass containers. Cannabinoids bind easily to plastic/glass surfaces; minimize sample handling. If the sample cannot be analyzed immediately, it may be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

INSTRUMENTATION REQUIRED

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

PROCEDURE

Analyzers with above indicated specifications are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameter used for each analyzer before performing the assay. Typical assay parameters used for the analyzers include a sample to antibody reagent (R₁) to enzyme conjugate reagent (R₂) ratio of 1:10:3,75, a 37°C incubation temperature, 2-4 min reading frames, and 340nm primary wavelength.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

CALIBRATION

For qualitative determinations, the reagent must be calibrated with the selected cut-off calibrator. For semi-quantitative determinations, the reagent can be calibrated with a 4 or 5 point calibration curve:

CUT-OFF	Calibrator levels (ng/mL)		Control Levels (ng/mL)
	QUALITATIVE	SEMI-QUANTITATIVE	
LOW 20 ng/mL	Level 0 (0) Level 1 (20)	Level 0 (0) Level 1 (20) Level 2 (50) Level 3 (100)	Control -25% (15) Control +25% (25)
HIGH 50 ng/mL	Level 0 (0) Level 2 (50)	Level 0 (0) Level 1 (20) Level 2 (50) Level 3 (100) Level 4 (200)	Control -25% (37,5) Control +25%(62,5)

Recalibrate each month, when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

INTERPRETATION

For qualitative determinations, the cut-off calibrator (20 or 50 ng/mL) of THC is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance ($\Delta A/\text{min}$) equal to, or greater than, that obtained with the cut-off calibrator is considered positive. A sample with a change in absorbance value lower than that obtained with the cut-off calibrator is considered negative.

For semi-quantitative determinations, a calibration curve with multiple calibrators is required. The concentration of THC in the sample may then be estimated from the calibration curve.

QUALITY CONTROL

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

The calibration curve can be validated with the controls level 15 and 25 ng/mL, or 37,5 and 62,5 ng/mL (Ref. 936071 and 936072, respectively) or with commercial controls.

LIMITATIONS

1. A positive result from the assay indicates only the presence of 11-nor⁹-carboxy- Δ^9 -THC
2. Refer to the intended use statement for details of recommended confirmation testing methods.
3. The test is designed for use with human urine only.

PRECISION AND SENSITIVITY

According to the EP5 standards (NCCLS), the reagent has been tested for 20 days, measuring several levels per duplicate twice a day (n=80).

Mean (ng/mL)	CV (%)		
	14.5	18.7	22.2
Total	7.2%	7.1%	5.3%
Within Run	6.0%	3.7%	2.9%
Between Run	1.7%	3.7%	4.1%
Between Day	3.6%	4.8%	1.7%
Mean (mAU/min)	358.6	364	368.9
Total	1.1%	1.1%	1.2%
Within Run	0.3%	0.3%	0.3%
Between Run	0.7%	0.8%	0.8%
Between Day	0.8%	0.7%	0.7%

Mean (ng/mL)	CV (%)		
	28.4	44.5	60.8
Total	7.7%	5.2%	3.3%
Within Run	4.2%	2.7%	1.5%
Between Run	4.8%	3.3%	1.4%
Between Day	4.3%	2.9%	2.6%
Mean (mAU/min)	387.6	407.3	435
Total	1.4%	1.5%	1.4%
Within Run	0.4%	0.5%	0.4%
Between Run	1.0%	1.1%	0.8%
Between Day	0.8%	0.8%	1.1%

Mean (ng/mL)	CV (%)		
	67.8	88.2	109.8
Total	4.4%	8.6%	9.4%
Within Run	2.6%	4%	4%
Between Run	2.7%	5.6%	5.6%
Between Day	2.4%	5.1%	5.1%
Mean (mAU/min)	447	466.2	476.3
Total	1.5%	1.4%	1.2%
Within Run	0.5%	0.4%	0.4%
Between Run	1.1%	1.1%	0.9%
Between Day	0.9%	0.8%	0.8%

ACCURACY

51 urine samples from an external quality assurance scheme (UKNEQAS) were tested using GC-MS. The results from all participating laboratories were averaged to give the definitive concentration of drug in the sample. The same 51 samples were analysed using the Liquid EIA and the results compared.

10 samples screened positive by the Liquid EIA, and the same 10 samples were confirmed positive by GC-MS.

41 samples screened negative by the Liquid EIA, and the same 41 samples were confirmed negative by GC-MS.

There were no mismatching data sets.

SPECIFICITY

Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cut-off calibrator (50ng/mL).

It is possible that other substances and/or factors not listed below may interfere with the test and cause false positive results.

The following compounds were found to be non-cross reacting at 100,000ng/mL:

6-Acetyl-Morphine	EDDP
Amitriptyline	EMDP
Amobarbital	Ephedrine
Amphetamine	Heroin
Aspirin	LAAM
Benzoylcegonine	MBDB
β -phenylethylamine	MDA
Cannabidiol	MDEA
Chlorpheniramine	MDMA
Cocaehtylene	Methadone
Cocaine	Methamphetamine
Codeine	Morphine
Cotinine	Oxycodone
Diazepam	Paracetamol
Dihydrocodeine	Pseudoephedrine
Ecgonine Methyl Ester	Temazepam

The following cannabinoids will produce a positive response relative to the 50ng/mL cut-off at the following concentrations:

Compound	Concentration(ng/mL)
11-hydroxy- Δ^9 -THC	80
Δ^9 -THC	200

BIBLIOGRAPHY

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988)
3. Marilyn A Huestics, "Marijuana", in Contemporary Practice in Clinical Toxicology, 2nd edition, Leslie M. Shaw, editor-in-chief. AACC, 2000.
4. Nahas, G.G., Cannabis: Toxicological Properties and Epidemiological Aspects, *Med J. Aust.*, 145:82 (1986)
5. Baselt, R.C. and R.H. Cravey, Disposition of Toxic Drugs and Chemicals in Man, 3rd Edition, pp. 780-783, Chicago, IL. Year Book Medical Publishers Inc. 1990.
6. Wall, M.E., D.R. Brine and M. Peres-Reyes, Metabolism of Cannabinoids in Man., in The Pharmacology of Marijuana, Brande, M.C. and S. Szara, editors, Raven Press, p. 93, 1976.
7. Chiang, C.N., and G. Barnett, Marijuana Pharmacokinetics and Pharmacodynamics, in Cocaine, Marijuana, Designer Drugs: Chemistry, Pharmacology, and Behavior, Redda. K.K., C.A. Walker, and G. Barnett, editors, CRC Press, Boca Raton, FL. 1989
8. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988)
3. Marilyn A Huestics, "Marijuana", in Contemporary Practice in Clinical Toxicology, 2nd edition, Leslie M. Shaw, editor-in-chief. AACCC, 2000.
4. Nahas, G.G., Cannabis: Toxicological Properties and Epidemiological Aspects, *Med J. Aust.*, 145:82 (1986)
5. Baselt, R.C. and R.H. Cravey, Disposition of Toxic Drugs and Chemicals in Man, 3rd Edition, pp. 780-783, Chicago, IL. Year Book Medical Publishers Inc. 1990.
6. Wall, M.E., D.R. Brine and M. Peres-Reyes, Metabolism of Cannabinoids in Man., in The Pharmacology of Marijuana, Brande, M.C. and S. Szara, editors, Raven Press, p. 93, 1976.
7. Chiang, C.N., and G. Barnett, Marijuana Pharmacokinetics and Pharmacodynamics, in Cocaine, Marijuana, Designer Drugs: Chemistry, Pharmacology, and Behavior, Redda. K.K., C.A. Walker, and G. Barnett, editors, CRC Press, Boca Raton, FL. 1989
8. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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Ref.: 930028 200 tests	Cont.	R1: 2 x 21 mL R2: 2 x 8 mL
Ref.: 930030 500 tests	Cont.	R1: 1 x 105 mL R2: 1 x 38 mL
Ref.: 930032 2500 tests	Cont.	R1: 2 x 250 mL R2: 1 x 190 mL