Cocaine Metabolite Urine

DoA Liquid Test

Cocaine Metabolite Assay for Urine IVD

Store at 2 - 8°C.

INTENDED USE

The Cocaine Metabolite Enzyme Immunoassay (EIA) is a homogeneous enzyme immunoassay system intended for use in the qualitative and semi-quantitative analysis of Benzoylcegonine (cocaine metabolite) in human urine.

The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method (1, 2). Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.

PRINCIPLE OF THE METHOD

The benzoylcegonine immunoassay is a homogeneous enzyme immunoassay (8) with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labelled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, benzoylcegonine-labelled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to free drug; the unbound benzoylcegonine-labelled G6PDH then exhibits its maximal enzyme activity.

Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

SIGNIFICANCE OF THE TEST

Cocaine (methylbenzoylcegonine) is an alkaloid found in the plant *Erythroxylon coca*, which is principally grown in South America. It is chemically, but not pharmacologically, related to atropine.

Cocaine is a CNS stimulant; however, it also exhibits numerous undesirable side effects including cardiac toxicity and behaviour responses such as paranoia and hallucinations. The most important action of cocaine clinically is its ability to block nerve conduction upon local application (3).

Cocaine sold on the street includes the hydrochloride salt and crack cocaine. The salt is frequently abused by inhalation or intravenous injection. Crack is a free base form of cocaine that produces a characteristic cracking sound when smoked.

Cocaine is rapidly metabolized, with less than 5% excreted unchanged in urine. The major metabolite is benzoylcegonine. Other notable metabolites are methylecgonine and ecgonine. Cocaine metabolites are detectable in urine for 1-3 days after moderate use (4, 5). However, for long term, heavy use, the metabolites may be found in urine for up to 3 weeks (6, 7). Cocaine readily passes through the placenta into the foetus. Cocaine abuse during pregnancy can adversely affect the foetal development and cause serious problems in the neonate (5).

REAGENTS

Antibody/Substrate Reagent (R₁): Contains monoclonal antibodies to benzoylcegonine, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide as preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains benzoylcegonine-labelled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide as preservative.

Avoid prolonged exposures of the reagent at temperatures higher than 25°C.

CALIBRATORS AND CONTROLS

Negative Human Urine (Level 0): Contains negative human urine with sodium azide as preservative (Ref.: 933010 and 933015).

Multidrug Calibrator Urine (Ref.: 933020):

Multidrug Calibrator Level 0: Contains negative human urine with sodium azide as preservative.

Multidrug Calibrator Level 1: Contains 150 ng/mL benzoylcegonine in human urine with sodium azide as preservative.

Multidrug Calibrator Level 2: Contains 300 ng/mL benzoylcegonine in human urine with sodium azide as preservative.

Multidrug Calibrator Level 3: Contains 500 ng/mL benzoylcegonine in human urine with sodium azide as preservative.

Multidrug Calibrator Level 4: Contains 1000 ng/mL benzoylcegonine in human urine with sodium azide as preservative.

Multidrug Calibrator Low Cut-Off: Contains 150 ng/mL benzoylcegonine in human urine with sodium azide as preservative (Ref.: 932910).

Multidrug Calibrator High Cut-Off: Contains 300 ng/mL benzoylcegonine in human urine with sodium azide as preservative (Ref.: 932915).

Multidrug Control Low Cut-Off: Contains 112,5 and 187,5 ng/mL benzoylcegonine in human urine with sodium azide as preservative (Ref.: 935910).

Multidrug Control High Cut-Off: Contains 225 and 375 ng/mL benzoylcegonine in human urine with sodium azide as preservative (Ref.: 935915).

PREPARATION

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

WARNING AND PRECAUTIONS

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

SPECIMEN COLLECTION

Urine samples may be collected in plastic or glass containers. Some plastics may adsorb drugs. Use a fresh urine specimen for the test. If the sample can not be analyzed immediately, it may be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

INSTRUMENTATION REQUIRED

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

PROCEDURE

Analyzers with above indicated specifications are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameter used for each analyzer before performing the assay. Typical assay parameters used for the analyzers include a sample to antibody reagent (R₁) to enzyme conjugate reagent (R₂) ratio of 1:10:3,75 respectively; a 37°C incubation temperature, 2-4 min reading frames, and 340nm primary wavelength.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

CALIBRATION

For qualitative determinations, the reagent must be calibrated with the selected cut-off calibrator. For semi-quantitative determinations, the reagent can be calibrated according to the following table:

CUT-OFF	Calibrator levels (ng/mL)		Control levels (ng/mL)
	QUALITATIVE	SEMI-QUANTITATIVE	
LOW 150 ng/mL	Level 0 (0) Level 1 (150)	Level 0 (0) Level 1 (150) Level 2 (300) Level 3 (500) Level 4 (1000)	Control - 25% (112,5) Control + 25% (187,5)
HIGH 300 ng/mL	Level 0 (0) Level 2 (300)		Control - 25% (225) Control + 25% (375)

The reagent should be recalibrated every month, when the controls are out of specifications (see Quality Control information), and when changing the reagent lot or the instrument settings.

Cocaine Metabolite Urine

DoA Liquid Test

INTERPRETATION

For qualitative determinations, the cut-off calibrator (150 or 300ng/mL) of benzoylecgonine is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance (Δ mA/min) equal to, or greater than, that obtained with the cut-off calibrator is considered positive. A sample with a change in absorbance value lower than that obtained with the cut-off calibrator is considered negative.

For semi-quantitative determinations, a calibration curve with multiple calibrators is required. The concentration of the Cocaine metabolite in the sample may then be estimated from the calibration curve.

QUALITY CONTROL

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

The calibration curve can be validated with the Control levels 112,5 and 187,5 ng/mL (ref. 935910) or levels 225 and 375 ng/mL (ref. 935915), or with commercial controls.

LIMITATIONS

1. A positive result from the assay indicates only the presence of cocaine metabolite.
2. Refer to the intended use statement for details of recommended confirmation testing methods
3. The test is designed for use with human urine only.

PRECISION

According to the EP5 standards (NCCLS), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80).

Qualitative analysis:

	CV (%)		
	Cut-off 150 ng/mL		
	112,5 ng/mL	150 ng/mL	187,5 ng/mL
Mean (mAU/min)	354,8	363,1	372,4
Total	1,2%	1,3%	1,2%
Within Run	0,4%	0,4%	0,3%
Between Run	0,6%	0,6%	0,6%
Between Day	1,0%	1,1%	1,0%

	Cut-off 300 ng/mL		
	225 ng/mL	300 ng/mL	375 ng/mL
Mean (mAU/min)	377,4	387	404,3
Total	1,3%	1,2%	1,2%
Within Run	0,4%	0,4%	0,3%
Between Run	0,6%	0,7%	0,6%
Between Day	1,1%	1,0%	1,0%

	Cut-off 150 and 300 ng/mL	
	0 ng/mL	1000 ng/mL
Mean (mAU/min)	299	441
Total	1,2%	1,1%
Within Run	0,5%	0,4%
Between Run	0,5%	0,5%
Between Day	1,0%	0,9%

Semi-quantitative analysis:

	CV (%)		
	Cut-off 150 ng/mL		
	112,5 ng/mL	150 ng/mL	187,5 ng/mL
Mean (ng/ml)	121.2	149.9	188.7
Total	7.6%	6.6%	6.8%
Within Run	4.3%	3.9%	2.5%
Between Run	5%	4.1%	5.5%
Between Day	3.9%	3.5%	3.2%

	Cut-off 300 ng/mL		
	225 ng/mL	300 ng/mL	375 ng/mL
Mean (ng/ml)	212.9	266	390.5
Total	6.7%	5.7%	6.9%
Within Run	3.3%	2.7%	3.2%
Between Run	4.5%	3.7%	4.7%
Between Day	3.7%	3.4%	3.9%

ACCURACY

52 urine samples from an external quality assurance scheme (UKNEQAS) were tested using GC-MS. The results from all participating laboratories were averaged to give the definitive concentration of drug in the sample. The same 52 samples were analysed using the Liquid EIA and the results compared.

13 samples screened positive by the Liquid EIA, and the same 13 samples were confirmed positive by GC-MS.

39 samples screened negative by the Liquid EIA, and the same 39 samples were confirmed negative by GC-MS.

There were no mismatching data sets.

SPECIFICITY

Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cut-off calibrator (300ng/mL).

It is possible that other substances and/or factors not listed below may interfere with the test and cause false positive results.

The following compounds were found to be non-cross reacting at 100,000ng/mL:

11-hydroxy- Δ 9-THC	EDDP
11-nor9-carboxy- Δ 9-THC	EMDP
6-Acetyl-Morphine	Ephedrine
Amitriptyline	Heroin
Amobarbital	LAAM
Amphetamine	MBDB
Aspirin	MDA
β -phenylethylamine	MDEA
Cannabidiol	MDMA
Chlorpheniramine	Methadone
Codeine	Methamphetamine
Cotinine	Morphine
Delta9-THC	Oxycodone
Diazepam	Paracetamol
Dihydrocodeine	Pseudoephedrine
Ecgonine Methyl Ester	Temazepam

The following cocaine will produce a positive response relative to the 300ng/mL cut-off at the following concentrations:

Compound	Concentration (ng/mL)
Cocaeethylene	100,000
Cocaine	100,000

BIBLIOGRAPHY

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988)
3. Goodman, L.S., and A. Gilman. The Pharmacological Basis of Therapeutics, 4th edition, 380, The MacMillan Co., 1970
4. Bouknight, L.G. and R.R. Bouknight. Cocaine- A Particularly Addictive Drug, *Postgrad. Med.* 83, 115, 1988
5. Benowitz, N.L. Clinical Pharmacology and Toxicology of Cocaine, *Pharmacol. & Toxicol.* 72, 3, 1993
6. Weiss, R.D. and F.H. Gawin. Protracted Elimination of Cocaine Metabolites in Long-term, High-dose Cocaine Abusers. *Am. J. Med.*, 85, 879, 1988
7. Burke, W.M. et al, Prolong Presence of Metabolite in Urine after Compulsive Cocaine Use, *J. Clin. Psychiatry*, 51, 145, 1990
8. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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- Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
- Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988)
- Goodman, L.S., and A. Gilman. The Pharmacological Basis of Therapeutics, 4th edition, 380, The MacMillan Co., 1970
- Bouknight, L.G. and R.R. Bouknight. Cocaine- A Particularly Addictive Drug, *Postgrad. Med.* 83, 115, 1988
- Benowitz, N.L. Clinical Pharmacology and Toxicology of Cocaine, *Pharmacol. & Toxicol.* 72, 3, 1993
- Weiss, R.D. and F.H. Gawin. Protracted Elimination of Cocaine Metabolites in Long-term, High-dose Cocaine Abusers. *Am. J. Med.*, 85, 879, 1988
- Burke, W.M. et al, Prolong Presence of Metabolite in Urine after Compulsive Cocaine Use, *J. Clin. Psychiatry*, 51, 145, 1990
- Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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Ref.: 930038
200 tests

Cont. : 2 x 21 mL R1
: 2 x 8 mL R2

Ref.: 930040
500 tests

Cont. : 1 x 105 mL R1
: 1 x 38 mL R2

Ref.: 930042
2500 tests

Cont. : 2 x 250 mL R1
: 1 x 190 mL R2