**Ethyl Alcohol Assay for Urine****IVD**

Store 2 - 8°C.

**INTENDED USE**

The Ethyl Alcohol is a homogenous enzymatic assay system intended for use in the quantitative analysis of ethyl alcohol in human urine samples, as well as human serum and plasma.

**Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.**

**PRINCIPLE OF THE METHOD**

The ethyl alcohol enzymatic assay is a homogeneous ready-to-use liquid reagent. The Alcohol Dehydrogenase (ADH) oxidizes ethyl alcohol to acetaldehyde, and reduces nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm. The concentration of ethyl alcohol is directly proportion to the ADH activity, measured at 340 nm wavelength.

**SIGNIFICANCE OF THE TEST**

Ethyl alcohol is present in some drinks, liquors, medicinal preparations and foods. After ethyl alcohol ingestion, the alcohol quickly spread to whole body and the most of it is metabolized in liver and excreted. Then, ethyl alcohol can be found in human urine and serum. Ethyl alcohol acts as a Central Nervous System depressant, and can lead to loss of alertness, stupor, coma and death and frequently cause public safety issues. Alcohol consumption during pregnancy carries a heavy risk of permanent mental and physical defects in the child, known as fetal alcohol syndrome (1,2,3,4,5). The concentration of ethyl alcohol is used for measuring legal impairment and forensic judgment. Ethyl alcohol measurements are also used in the diagnosis and treatment of alcohol intoxication and poisoning.

**REAGENTS**

**Buffer Reagent (R<sub>1</sub>)**: Contains Tris buffer with sodium azide as preservative.

**Enzyme Reagent (R<sub>2</sub>)**: Contains alcohol dehydrogenase and nicotinamide adenine dinucleotide, stabilizers and sodium azide as preservative.

**Avoid prolonged exposures of the reagent at temperatures higher than 25°C.**

**CALIBRATORS/CONTROLS**

**Negative Human Urine (Level 0)**: Contains negative human urine with sodium azide as preservative (Ref.: 933010).

**Ethyl Alcohol 50 mg/dL Control**: Contains 50 mg/dL ethyl alcohol in phosphate buffer with sodium azide as preservative (Ref.: 932995).

**Ethyl Alcohol 100 mg/dL Calibrator**: Contains 100 mg/dL ethyl alcohol in phosphate buffer with sodium azide as preservative (Ref.: 933000).

**Ethyl Alcohol 300 mg/dL Control**: Contains 300 mg/dL ethyl alcohol in phosphate buffer with sodium azide as preservative (Ref.: 933005).

**Ethyl Alcohol Calibrator/Control Urine**: (Ref.: 933060).

**PREPARATION**

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

**WARNING AND PRECAUTIONS**

- This test is for *in vitro* diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

**SPECIMEN COLLECTION**

Urine, plasma and serum samples may be collected in plastic or glass containers. They should have a tight cap to prevent alcohol evaporation. Anticoagulants citrate, EDTA, fluoride-oxalate, and heparin for plasma specimens may be used for this assay. Use a fresh urine specimen for the test. If the sample can not be analyzed immediately, it may be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing.

Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

**INSTRUMENTATION REQUIRED**

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting sample, mixing reagents, measuring enzyme rates at 340 nm and timing the reaction accurately can be used to perform this homogeneous immunoassay.

**PROCEDURE**

Analyzers with above indicated specifications are suitable for performing this homogeneous enzyme immunoassay. Refer to the specific parameter used for each analyzer before performing the assay. Typical assay parameters used for the analyzers include a sample to buffer reagent (R<sub>1</sub>) to enzyme reagent (R<sub>2</sub>) ratio of 1:10:3,75 respectively; a 37°C incubation temperature, 2-4 min reading frames, and 340nm primary wavelength.

Spinreact has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

**CALIBRATION**

The reagent has to be calibrated with the Negative human urine and the Ethyl Alcohol 100 mg/dl calibrator.

The reagent should be recalibrated when the controls are out of specifications (see Quality Control information), and when changing the reagent lot or the instrument settings.

**INTERPRETATION**

Level (in Blood)	Sporadic Drinkers	Chronic Drinkers
100 mg/dL(0.1%)	Legally intoxicated	Minimal signs
200-250 mg/dL (0.2-0.25%)	Alertness lost, lethargic	Effort needed to maintain control
300-350mg/dL (0.30-0.35%)	Stupor to coma	Drowsy, Slow
>500 mg/dL (>0.50%)	Death possible	Coma

**QUALITY CONTROL**

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

Use controls (50 and 300 mg/dl) to validate the calibration.

**LIMITATIONS**

1. The legal alcohol intoxication levels vary. The test result should be interpreted in light of clinical signs and symptoms.
2. Ethyl Alcohol is volatile. Precaution as suggested in the specimen collection and required to prevent alcohol evaporation from calibrators, controls and samples.
3. The test is designed for use with human urine, serum and plasma only.

**PRECISION AND SENSITIVITY**

According to the EP5 standards (CLSI), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80).

**Qualitative method:**

	CV (%)		
	50 mg/dL	75 mg/dL	100 mg/dL
Mean (mAU/min)	148,7	221,7	280
Total	2,1	1,8	1,5
Within Run	1,0	1,2	0,8
Between Run	1,2	1,2	1,3
Between Day	1,5	0,6	0,3

	CV (%)	
	125 mg/dL	250 mg/dL
Mean (mAU/min)	360,3	705,9
Total	1,5	1,7
Within Run	1,1	1,5
Between Run	0,7	0,8
Between Day	0,8	0

**Semi-Quantitative method:**

	CV (%)		
	50 mg/dL	75 mg/dL	100 mg/dL
Mean (ng/mL)	52,6	78,6	99,4
Total	2,2	1,8	1,5
Within Run	1,0	1,2	0,8
Between Run	1,2	1,2	1,3
Between Day	1,5	0,6	0,3

	CV (%)	
	125 mg/dL	250 mg/dL
Mean (ng/mL)	128	251,2
Total	1,5	1,7
Within Run	1,1	1,5
Between Run	0,7	0,8
Between Day	0,8	0

**ACCURACY**

Urine and serum specimens containing ethyl alcohol were analyzed by the commercial Ethyl Alcohol Assay (A), Ethyl Alcohol Test (B), and by the Spinreact Ethyl Alcohol Enzymatic Assay. Linear regression analyses of the results are summarized in the following table.

	<b>Spinreact vs A</b>		<b>Spinreact vs B</b>	
	Serum	Urine	Serum	Urine
N =	55	110	55	110
Slope	1.059	1.08	0.991	1.003
Intercept	21.03	-7.65	11.03	-0.95
Correlation	0.971	0.998	0.995	0.999

**BIBLIOGRAPHY**

1. Baselt RC, Disposition of Toxic Drugs and Chemicals in Man. 3<sup>rd</sup> edition, Chicago, IL. Year Book Medical Publishers Inc.; 1989:322-324.
2. Beutler HO, Ethanol. In: Bergmeyer HU, ed Methods of Enzyme Analysis, Vol VI. 3<sup>rd</sup>, New York: Academic Press, 1984:598-606.
3. Wyngaerde JB, Smith LH Jr, eds. Cecil Textbook of Medicine. Philadelphia, PA: WB Saunders Co; 1988:48-52.
4. Ellenhorn MJ, Barceloux DG, Medical Toxicology. New York, NY: Elsevier Science publishing Company, Inc. 1988:525-526, 782-796.
5. Tietz NW, ed Textbook of Clinic Chemistry. Philadelphia, PA: WB Saunders Co; 1986:1704-1706, 1692-1694.

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1. Baselt RC, Disposition of Toxic Drugs and Chemicals in Man. 3<sup>rd</sup> edition, Chicago, IL. Year Book Medical Publishers Inc.; 1989:322-324.
2. Beutler HO, Ethanol. In: Bergmeyer HU, ed Methods of Enzyme Analysis, Vol VI. 3<sup>rd</sup>, New York: Academic Press, 1984:598-606.
3. Wyngaarde JB, Smith LH Jr, eds. Cecil Textbook of Medicine. Philadelphia, PA: WB Saunders Co; 1988:48-52.
4. Ellenhorn MJ, Barceloux DG, Medical Toxicology. New York, NY: Elsevier Science publishing Company, Inc. 1988:525-526, 782-796.
5. Tietz NW, ed Textbook of Clinic Chemistry. Philadelphia, PA: WB Saunders Co; 1986:1704-1706, 1692-1694.

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Ref.: 930108  
200 tests

Cont. R1: 2 x 21 mL  
R2: 2 x 8 mL

Ref.: 930110  
500 tests

Cont. R1: 1 x 105 mL  
R2: 1 x 38 mL

Ref.: 930112  
2500 tests

Cont. R1: 2 x 250 mL  
R2: 1 x 190 mL