Creatine kinase

NAC. Kinetic UV. Liquid


Quantitative determination of creatine kinase liquid (CK) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Kinetic determination of the creatine kinase based upon IFCC and DGKC recommendations.

Creatine kinase (CK) catalyses the reversible transfer of a phosphate group from phosphocreatine to ADP. This reaction is coupled to those catalysed by hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6P-DH):

The rate of NADPH formation, measured photometrically, is proportional to the catalytic concentration of CK present in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Creatine kinase is a cellular enzyme with wide tissue distribution in the body. Its physiological role is associated with adenosine triphosphate (ATP) generation for contractile or transport systems. Elevated CK values are observed in diseases of skeletal muscle and after myocardial infarction^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Imidazol, pH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
	N-Acetyl-L-Cysteine	25 mmol/L
	Magnesium acetate	12,5 mmol/L
	NADP	2,52 mmol/L
	EDTA	2,02 mmol/L
	Hexokinase	≥6 800 U/L
R 2	ADP	15,2 mmol/L
	AMP	25 mmol/L
	di-Adenosine-5- penta-P	103 mmol/L
	Glucose-6-phosphate DH	≥8 800 U/L
	Creatine phosphate	250 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use the tablets if appears broken.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1,60.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum free of hemolysis or heparin plasma.

Stability 7 days at 2-8°C, protected from light.

The creatin kinase activity decreases 10% after 1 day at 2-5°C or after 1 hour at 15-25°C.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	240
Test	CKNAC	Vol. R2	60
Full Name	Creatin Kinase	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	5	Decimal	1 Slope 1
Vol. Sample Increas.		Unit	U/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	340	Reagent Blank	0-0
Sec. Wave.		React. Time	57-77

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	CK	No.	**
Full name	CK-NAC	Print name	CK
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	U/L	Decimal	0
Reagent Blank	-	React. Time	62 - 82
Vol. Sample	5 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until 32 days. After this period the Calibration must be performed again in order to obtain good results.

REFERENCE VALUES¹

Men, up to 195 U/L
Women, up to 170 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 2,12 U/L to linearity limit of 2000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay		Inter-assay	
	Mean (U/L)	CV (%)	Mean (U/L)	CV (%)
Mean (U/L)	147	0,84	145	2,01
SD	1,23	0,73	485	8,97
CV (%)	0,84	0,73	2,91	1,85

Sensitivity: 1 U/L = 0,00012 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained were the following:

Correlation coefficient (r)²: 0,9995

Regression equation: y = 1,0846x - 0,3512.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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Ref: MD41250	Cont.	R1: 4 x 40 mL
		R2: 2 x 20 mL

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BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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Ref: MD41250	Cont.	R1: 4 x 40 mL
		R2: 2 x 20 mL

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BIBLIOGRAFIA

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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Ref: MD41250	Cont.	R1: 4 x 40 mL
		R2: 2 x 20 mL