# Alkaline phosphatase


p-Nitrophenylphosphate. kinetic. Liquid. DGKC

## Quantitative determination of alkaline phosphatase (ALP) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Alkaline phosphatase (ALP) catalyses the hydrolysis of p-nitrophenyl phosphate at pH 10.4, liberating p-nitrophenol and phosphate, according to the following reaction:

The rate of p-Nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of alkaline phosphatase present in the sample<sup>1,2</sup>.

### CLINICAL SIGNIFICANCE

Alkaline phosphatase is an enzyme present in almost all weaves of the organism, being particularly high in bone, liver, placenta, intestine and kidney. Both increases and decreases of plasma ALP are of importance clinically.

Causes of increased plasma ALP: Paget's disease of bone, obstructive liver disease, hepatitis, hepatotoxicity caused by drugs or osteomalacia.

Causes of decreased plasma ALP: Cretinism and vitamin C deficiency<sup>1,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R 1</b>	Diethanolamine (DEA) pH 10.4	1 mmol/L
Buffer	Magnesium chloride	0,5 mmol/L
<b>R 2</b>	p-Nitrophenylphosphate (pNPP)	10 mmol/L
Substrate		

### PREPARATION

All the reagents are ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm  $\geq 1,30$ .

### ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

### SAMPLES

Serum or heparinized plasma<sup>1</sup>. Use unhemolyzed serum, separated from the clot as soon as possible. Stability: 3 days at 2-8°C.

### QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

## SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	240
Test	ALP	Vol. R2	60
Full Name	Alkaline Phosp.	Vol. R3	
Standard nº		Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	5	Decimal	1 Slope 1
Vol. Sample Increas.		Unit	U/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	412	Reagent Blank	0-0
Sec. Wave.		React. Time	50-70

## SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	ALP	No.	**
Full name	ALK PHOSPH	Print name	ALP
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	412	Sec. Wave.	
Unit	U/L	Decimal	0
Reagent Blank	0 - 0	React. Time	57 - 77
Vol. Sample	5 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

*The Calibration is stable until 30 days. After this period the Calibration must be performed again in order to obtain good results.*

### REFERENCE VALUES<sup>1</sup>

Children (1-14 years) < 645 U/L

Adults 98 - 279 U/L

Factors affecting ALP activities in a normal population include exercise, periods of repaid growth in children and pregnancy.

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 0,6845 U/L to linearity limit of 1200 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

### Precision:

Mean (U/L)	Intra-assay (n=20)		Inter-assay (n=20)	
	174	443	175	434
SD	0,72	1,56	6,88	11,93
CV (%)	0,41	0,35	3,93	2,75

**Sensitivity:** 1 U/L = 0,0003 ΔA/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following: Correlation coefficient (r)<sup>2</sup>:0,99938.

Regression equation: y= 1,025x - 1,105.

The results of the performance characteristics depend on the analyzer used.

### BIBLIOGRAPHY

1. Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
2. Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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