
Quantitative determination of chloride ion IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of a red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed is proportional to the chloride ion concentration in the sample^{1,2,3,4}.

CLINICAL SIGNIFICANCE

It is important clinically the determination of chloride due regulation of osmotic pressure of extra cellular fluid and to its significant role in acid-base balance. Increases in chloride ion concentration may be found in severe dehydration, excessive intake of chloride, severe renal tubular damage and in patients with cystic fibrosis.

Decrease in chloride ion concentration may be found in metabolic acidosis, loss from prolonged vomiting and chronic pyelonephritis^{2,7,8}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R Thiocyanate-Hg	Mercuric thiocyanate	4 mmol/L
	Ferric nitrate	40 mmol/L
	Mercuric nitrate	2 mmol/L
	Nitric acid	45 mmol/L

PRECAUTIONS

R: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm \geq 0.15.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment (Note 1, 2, 3).

SAMPLES

- Serum or plasma free of hemolysis and separated from cells as rapidly as possible. Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.

- Urine¹: Collect 24-hour urine specimen in chloride free containers. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the sample: Ion chloride is stable 1 week at room temperature (15-25°C), in refrigerator (2-8°C) or frozen (-20°C) temperatures.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	300
Test	CL	Vol. R2	
Full Name	Chloride	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	3	Decimal	0,1 Slope 1
Vol. Sample Increas.		Unit	mmo/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	505	Reagent Blank	10-11
Sec. Wave.		React. Time	26-27

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	CL	No.	**
Full name	CHLORIDE	Print name	CL
Reac. Type	End Point	Direction	Increase
Pri. Wave.	505	Sec. Wave.	
Unit	mmol/L	Decimal	0.01
Reagent Blank	10 - 11	React. Time	46 - 47
Vol. Sample	3 ul	R1	300 ul
Increased		R2	
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

REFERENCE VALUES¹

Serum or plasma:	95 - 115 mmol/L	CSF:	95 - 110 mmol/L
Urine:	110 - 250 mmol/24h	Sweat:	Up to 60 mmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0,454 mmol/L to *linearity limit* of 190 mmol/L.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with distilled water and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mmol/L)	84,2	114	82,5
SD	0,81	0,62	1,07	1,87
CV (%)	0,96	0,55	1,30	1,68

Sensitivity: 1 mmol/L = 0,00471 A.

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents. The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,96731.

Regression equation: y=0,990x + 0,100.

The results of the performance characteristics depend on the analyzer used.

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