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
BIBLIOGRAFIA

- Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
- Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

PRESENTACION

Ref: MD1001247

Cont.


R1: 4 x 30 mL, R2: 4 x 500 mg, R3: 1 x 10 mL

Quantitative determination of iron IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The iron is dissociated from transferrin-iron complex in weakly acid medium. Liberated iron is reduced into the bivalent form by means of ascorbic acid. Ferrous ions give with FerroZine a coloured complex:

 The intensity of the color formed is proportional to the iron concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

The iron is the component of a great number of enzymes. The myoglobin, muscular protein, contains iron, as well as the liver. Iron is necessary for the hemoglobin production, molecule that transports oxygen inside red globules. Their deficit in the last causes the ferropenic anemia. High levels of iron are found in hemochromatosis, cirrhosis, hepatitis and in increased transferrin levels.

 The variation day to day is quite marked in healthy people^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Buffer	Acetate pH 4,9	100 mmol/L
R 2	Reductant	Ascorbic acid	99,7%
R 3	Color	FerroZine	40 mmol/L

PREPARATION

Dissolve (→) the contents of one tube of reductant (R2) in one bottle of R 1 Buffer. Cap and mix gently to dissolve contents. Stability: 3 months at 2-8°C or 1 month at 15-25°C. R3 is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 562 nm ≥ 0,020.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment

SAMPLES

Serum or heparinized plasma.

Free of hemolysis and separated from cells as rapidly as possible.

 Stability of the sample: 2-8°C for 7 days¹.

REFERENCE VALUES⁵

Male 65 - 175 µg/dL ≅ 11,6 - 31,3 µmol/L (Note 4)

Female 40 - 150 µg/dL ≅ 7,16 - 26,85 µmol/L (Note 4)

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	240
Test	IRON	Vol. R2	20
Full Name	Iron	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	40	Decimal	1 Slope
	1		
Vol. Sample Increas.		Unit	µg/dL Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	End Point	Direction	
	Increase		
Pri. Wave.	570	Reagent Blank	39-40
Sec. Wave.		React. Time	75-76

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	FE	No.	**
Full name	IRON	Print name	FE
Reac. Type	End Point	Direction	Increase
Pri. Wave.	570	Sec. Wave.	
Unit	ug/dL	Decimal	0.1
Reagent Blank	46 - 47	React. Time	81 - 82
Vol. Sample	40 ul	R1	240 ul
Increased		R2	20 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until 35 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0,850 µg/dL to *linearity limit* of 1000 µg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (µg/dL)	113	250	111	249
SD	0,89	0,72	3,51	6,29
CV (%)	0,79	0,29	3,17	2,52

Sensitivity: 1 µg/dL = 0,00104 A.

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents.

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,9934.

Regression equation: y = 1,0243x - 3,877.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. It is recommended to use disposable material. If glassware is used the material should be soaking for 6 h in diluted HCl (20% v/v) and then thoroughly rinsed with distilled water and dried before use.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.
4. Strongly method dependent.

BIBLIOGRAPHY

1. Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
2. Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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