Quantitative determination of total protein IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Proteins give an intensive violet-blue complex with copper salts in an alkaline medium. Iodide is included as an antioxidant.

 The intensity of the color formed is proportional to the total protein concentration in the sample^{1,4}.

CLINICAL SIGNIFICANCE

The proteins are macromolecular organic compounds, widely distributed in the organism. They act like structural and transport elements. The proteins of the serum are divided in two fractions, albumin and globulins

The determination of total proteins is useful in the detection of:

- High protein levels caused by hemoconcentration like in the dehydrations or increase in the concentration of specific proteins.
- Low protein level caused by hemodilution by an impaired synthesis or loss (as by hemorrhage) or excessive protein catabolism^{4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Sodium potassium tartrate	15 mmol/L
	Sodium iodide	100 mmol/L
	Potassium iodide	5 mmol/L
	Copper (II) sulphate	19 mmol/L
Biuret		

PRECAUTIONS

H314-Causes severe skin burns and eye damage. H412-Harmful to aquatic life with long lasting effects.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 546 nm \geq 0,22.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

 Serum or heparinized plasma¹:

Stability of the sample: 1 month at refrigerator (2-8°C).

REFERENCE VALUES¹

Adults: 6,6 – 8,3 g/dL

Newborn: 5,2 – 9,1 g/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	300
Test	TProt	Vol. R2	
Full Name	Total Protein	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	8	Decimal	0.01 Slope 1
Vol. Sample Increas.		Unit	g/dL Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	546	Reagent Blank	10-11
Sec. Wave.		React. Time	46-47

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	TP	No.	**
Full name	TOTAL PROTEIN	Print name	TP
Reac. Type	End Point	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	g/dL	Decimal	0,01
Reagent Blank	10 - 11	React. Time	46 - 47
Vol. Sample	8 ul	R1	300 ul
Increased		R2	
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **26 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS
Measuring range: From *detection limit* of 0,007 g/dL to *linearity limit* of 14 g/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (g/dL)	SD	Mean (g/dL)	SD
Mean (g/dL)	6,53	4,89	6,77	5,08
SD	0,01	0,01	0,07	0,05
CV (%)	0,21	0,24	1,05	0,94

Sensitivity: 1 g/dL = 0,0825 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)²: 0,97002

Regression equation y= 0,954x +0,511.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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