# LDL Cholesterol D

Enzymatic colorimetric. Liquid

### Quantitative determination of LDL cholesterol IVD




Store at 2-8°C

#### PRINCIPLE OF THE METHOD



Direct determination of serum LDLc (low-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation steps.

The assay takes place in two steps.

- 1° Elimination of lipoprotein no-LDL

- 2° Measurement of LDLc

The intensity of the color formed is proportional to the LDLc concentration in the sample.

#### CLINICAL SIGNIFICANCE

The LDLc particle is lipoproteins that transport cholesterol to the cells. Often called "bad cholesterol" because high levels are risk factor for coronary heart disease and are associated with obesity, diabetes and nephrosis<sup>1,2,9</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

#### REAGENTS

<b>R 1</b> Enzymes	PIPES Buffer pH 7,0	50 mmol/L
	Cholesterol esterase (CHE)	≥600 U/L
	Cholesterol oxidase (CHOD)	≥500 U/L
	Catalase	≥600 U/mL
	TOOS	2 mmol/L
<b>R 2</b> Enzymes	PIPES Buffer pH 7,0	50 mmol/L
	4 - Aminoantipyrine (4-AA)	4 mmol/L
	Peroxidase (POD)	≥4 KU/L

#### PREPARATION

**R 1 and R 2:** Are ready to use.

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

**R 1 and R 2:** Once opened is stable 4 weeks at 2-8°C.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.

#### ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

#### SAMPLES

Serum, heparinized plasma or EDTA plasma. If any sample show precipitates, centrifuge before using<sup>5</sup>.

Serum stable 6 days at 2-8°C. Do not freeze the samples.

#### QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

#### SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	210
Test	LDL	Vol. R2	70
Full Name	LDL Cholesterol	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	3	Decimal	0.1 Slope 1
Vol. Sample Increas.		Unit	mg/dL Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	570	Reagent Blank	0-0
Sec. Wave.		React. Time	43-76

#### SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	LDL	No.	**
Full name	LDL Cholesterol	Print name	LDL
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	605	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	0 - 0	React. Time	51 - 82
Vol. Sample	3 ul	R1	210 ul
Increased		R2	70 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

#### REFERENCE VALUES<sup>6,7,8</sup>

Optimal	< 100 mg/dL
Near or above optimal	100-129 mg/dL
Borderline high	130-160 mg/dL
High	> 160 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

#### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 10 mg/dL to *linearity limit* of 976 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

Media (mg/dL)	Intraserie (n= 20)		Interserie (n= 20)	
	31,4	67,8	32,1	68,1
SD	0,42	1,11	0,92	2,02
CV (%)	1,35	1,64	2,87	2,97

**Sensibility:** 1 mg/dL = 0,001784 (A).

**Accuracy<sup>10,11</sup>:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99123.

Regression equation: y = 0,914x + 1,58283.

The results of the performance characteristics depend on the analyzer used.

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