
Quantitative determination of Lipase IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The pancreatic lipase in presence of colipase, desoxycholate and calcium ions, hydrolyses the substrate 1-2-O-dilauryl-rac-glycero-3-glutaric acid-(6' -methylresorufin)-ester. The sequence of reactions involved in the enzymatic direct lipase determination is the following:

The rate of methylresorufin formation, measured photometrically, is proportional to the catalytic concentration of lipase present in the sample.

CLINICAL SIGNIFICANCE

Lipase (LPS) is a pancreatic enzyme necessary for the absorption and digestion of nutrients that catalyzes the hydrolysis of glycerol esters of fatty acids. Determination of LPS is used for diagnosis of diseases of pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct^{1,7,8}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 8.3	40 mmol/L
	Colipase	≥ 1 mg/L
	Desoxycholate	1,8 mmol/L
	Taurodesoxycholate	7,2 mmol/L
R 2 Substrate (micro-emulsion)	Tartrate pH 4,0	15 mmol/L
	Lipase Substrate	≥ 0,7 mmol/L
	Calcium chloride (CaCl ₂)	0,1 mmol/L

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 580 nm ≥ 1,00.
- R2 is a turbid orange-colored micro-emulsion, discard if turning to red.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment ^(Note 2).

SAMPLES

Serum or plasma with sodium citrate, EDTA or heparin¹.
Avoid repeated frozen and unfrozen. Stability: 2 days at 2-8°C.

REFERENCE VALUES¹

≤ 38 U/L (U/L methylresorufin at 37°C).
These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures:
SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).
If control values are found outside the defined range, check the instrument, reagents and technique for problems.
Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	200
Test	LIPASE	Vol. R2	40
Full Name	Lipase	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	2	Decimal	0.1 Slope 1
Vol. Sample Increas.		Unit	U/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	570	Reagent Blank	0-0
Sec. Wave.		React. Time	53-60

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	LIPASE	No.	**
Full name	LIPASE	Print name	LIPASE
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	570	Sec. Wave.	
Unit	U/L	Decimal	0.1
Reagent Blank	0 - 0	React. Time	60 - 67
Vol. Sample	2 ul	R1	200 ul
Increased		R2	40 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 5 U/L to linearity limit of 250 U/L.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	40,2	59,35	38,5
SD	0,410	0,875	1,10	1,25
CV (%)	1,02	1,47	2,86	2,13

Sensitivity: 1 U/L= 0,00059792 (A)

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 101 samples were the following:
Correlation coefficient (r)²:0,99732.

Regression equation: y= 0,50054x + 3,9443.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. In some storage conditions (i.e. storage at a temperature lower than the one indicated) a precipitate may appear in the vial that will not influence that the reagent performance; however, it is recommended to resuspend the product with a slight rotation.
2. In order to avoid contamination it is recommended to use disposable material.

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