
Quantitative determination of gamma-glutamyl transferase (γ-GT)

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Gamma-glutamyl transferase (γ-GT) catalyses the transfer of γ-glutamyl group from γ-glutamyl-p-nitroanilide to acceptor glycylglycine, according to the following reaction:

The rate of 2-nitro-5-aminobenzoic acid formation, measured photometrically, is proportional to the catalytic concentration of γ-GT present in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Gamma-glutamyl transferase (γ-GT) is a cellular enzyme with wide tissue distribution in the body, primarily in the kidney, pancreas, liver and prostate.

Measurements of gamma-glutamyl transferase (γ-GT) activity are used in the diagnosis and treatment of hepatobiliary diseases such as biliary obstruction, cirrhosis or liver tumours^{1,2,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	TRIS	100 mmol/L
Buffer	Glycylglycine	100 mmol/L
R 2	L-γ-glutamyl-3-carboxy-4-nitroanilide	3 mmol/L
Substrate		

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm $\geq 1,20$.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum¹. γ-GT is stable for at least 3 days at 2-8°C, 8 hours at 15-25°C and 1 month at -20°C.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
No	**	Vol. R1	240
Test	GGT	Vol. R2	60
Full Name	GGT	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	30	Decimal	1
Vol. Sample Increas.		Unit	U/L
Vol. Sample Dec		Slope	1
		Inter.	0
REACTION PARAMETERS			
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	412	Reagent Blank	0-0
Sec. Wave.		React. Time	50-70

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	GGT	No.	**
Full name	GGT	Print name	GGT
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	412	Sec. Wave.	
Unit	U/L	Decimal	0
Reagent Blank	0 - 0	React. Time	57 - 77
Vol. Sample	30 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **35 days**. After this period the Calibration must be performed again in order to obtain good results.

REFERENCE VALUES¹

	25°C	30°C	37°C
Women	4-18 U/L	5-25 U/L	7-32 U/L
Men	6-28 U/L	8-38 U/L	11-50 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 2 U/L to linearity limit of 300 U/L. If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	Mean	SD
Mean (U/L)	38,3	190	40,1	198
SD	0,39	0,53	0,82	2,30
CV (%)	1,03	0,28	2,05	1,16

Sensitivity: 1 U/L = 0,0008 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99990.

Regression equation: y= 1,334x - 1,493.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

- Gendler S. γ-GT. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1120-1123.
- Persijn J P et al. J Clin Chem Clin Biochem 1976; (14) 9: 421-427.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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