# Protein in urine and CSF

Pyrogallol red. Colorimetric

## Quantitative determination of total urinary and CSF protein IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Proteins react in acid solution with pyrogallol red and molybdate to form a colored complex.

The intensity of the color formed is proportional to the protein concentration in the sample<sup>1,2</sup>.

### CLINICAL SIGNIFICANCE

In healthy persons, the urine contains no protein or only a trace amount of protein; normally the glomeruli prevent passage of protein from the blood to the glomerular filtrate. Glomerular injury causes increased permeability to plasma proteins, resulting in proteinuria, which refers to the presence of protein in the urine.

A persistent finding of proteinuria is the single most important indication of renal disease.

Elevated concentration of protein in cerebro-spinal fluid (CSF) can be caused by infections and intracranial pressure<sup>1,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

R	Pyrogallol red	50 µmol/L
	Sodium molybdate	0,04 mmol/L

### PRECAUTIONS

R: H371. May cause damage to organs. Contains: Methanol (CH<sub>3</sub>OH).

Follow the precautionary statements given in MSDS and label of the product.

### PREPARATION

The reagent is ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 600 nm ≥ 0,30.

### ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

### SAMPLES

- Urine 24 h: Stability 8 days at 2-8°C.
- Cerebrospinal fluid (CSF): Stable 4 days at 2-8°C.

### REFERENCE VALUES<sup>5</sup>

Urine:	<100 mg/24h (<150 mg/24h in pregnancy)
CSF:	Children 300 -1000 mg/L
	Adults 150 - 450 mg/L

These values are for orientation purpose; each laboratory should establish its own reference range.

### QUALITY CONTROL

**A calibration using a urine protein calibrator is needed (Ref.: 1002450). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.**

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

## SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	300
Test	PTUR	Vol. R2	
Full Name	Urine Protein	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	6	Decimal	1 Slope 1
Vol. Sample Increas.		Unit	mg/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	605	Reagent Blank	10-11
Sec. Wave.		React. Time	46-47

## SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	uProt	No.	**
Full name	Urine Protein	Print name	uProt
Reac. Type	End Point	Direction	Increase
Pri. Wave.	605	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	10 - 11	React. Time	46 - 47
Vol. Sample Increased	6 ul	R1	300 ul
Decreased		R2	
Sample blank		R3	
		R4	

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit 9,44 mg/L to linearity limit of 4000 mg/L.

If the concentration is greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

	Intra-assay (n= 20)		
Mean (mg/L)	220	536	1014
SD	3,7	4,0	5,2
CV (%)	2,28	0,75	0,51

	Inter-assay (n= 20)		
	216	499	1018
	18,3	26,1	166,1
	7,35	5,22	16,43

**Sensitivity:** 1mg/L = 0,00026 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,9338

Regression equation: y = 0,4294x - 5,4159

The results of the performance characteristics depend on the analyzer used.

### NOTES

- In order to avoid possible contaminations, it is recommended cleaning the needle after this parameter determination. Due to its content the system can be easily contaminated.

### BIBLIOGRAPHY

- Orsonneau JL et al. An improved Pyrogallol Red-Molybdate Method for Determining Total Urinary Protein. Clin Chem 1989 (35):2233-2236.
- Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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### BIBLIOGRAFIA

1. Orsonneau JL et al. An improved Pyrogallol Red-Molybdate Method for Determining Total Urinary Protein. Clin Chem 1989; 35:2233-2236.
2. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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