HDL Cholesterol D

Direct. Enzymatic colorimetric




Quantitative determination of HDL cholesterol IVD

Store at 2-8°C



PRINCIPLE OF THE METHOD

Directly determination of serum HDLc (high-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation of the sample^{3,5}. The assay takes place in two steps.

- 1° Elimination of lipoprotein no-HDL

- 2° Measurement of HDLc

The intensity of the color formed is proportional to the HDLc concentration in the sample.

CLINICAL SIGNIFICANCE

HDL particles are high-density lipoproteins that transport cholesterol from the body tissues to the liver. Since HDL can remove cholesterol from the arteries and carry it back to the liver for their excretion, HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease and coronary artery disease.

A low HDL cholesterol levels, is considered a greater heart disease risk^{1,2,4}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 6,6	100 mM
	N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (HDAOS)	0,7 mM
	Cholesterol Esterase	≥ 800 U/L
	Cholesterol oxidase	≥ 500 U/L
	Catalase	≥ 300 U/L
	Ascorbic oxidase	≥ 3000 U/L
R 2	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 7,0	1,1 mmol/L
	4 - Aminoantipyrine (4-AA)	100 mM
	Peroxidase	≥ 3500 U/L

PREPARATION

- R 1 and R 2: Are ready to use.

STORAGE AND STABILITY^(Note 1)

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze the reagents.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum, heparinized plasma or EDTA plasma. If any samples show precipitates, centrifuge before using.

Stability of the sample: 6 days at 2-8°C and 1 year when stored at -70°C.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME			
Nº	**	Vol. R1	210		
Test	HDL	Vol. R2	70		
Full Name	HDL Cholesterol	Vol. R3			
Standard nº	1	Vol. R4			
SAMPLE VOLUME		RESULT SETUP			
Vol. Sample Stand.	2	Decimal	0.1	Slope	1
Vol. Sample Increases.		Unit	mg/dL	Inter.	0
Vol. Sample Dec					
REACTION PARAMETERS					
Reac. Type	Fixed Time	Direction	Increase		
Pri. Wave.	570	Reagent Blank	0-0		
Sec. Wave.		React. Time	43-76		

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	HDL	No.	**
Full name	HDL Cholesterol	Print name	HDL
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	570	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	0 - 0	React. Time	51 - 82
Vol. Sample	2 ul	R1	210 ul
Increased		R2	70 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

REFERENCE VALUES

	Men	Women
Low risk	> 50 mg/dL	> 60 mg/dL
Normal risk	35 - 50 mg/dL	45 - 60 mg/dL
High risk	< 35 mg/dL	< 45 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 5,0 mg/dL to *linearity limit* of 151 mg/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	28,0	76,1	27,5	75,3
SD	0,25	0,81	1,26	2,04
CV (%)	0,89	1,06	4,60	2,71

Sensibility: 1mg/L = 0,001399 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,938.

Regression equation: y = 0,9825x + 1,41606.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. The reagent 2 presents yellowish coloration due to the peroxidase, but it does not affect its functionality.

BIBLIOGRAPHY

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S.; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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BIBLIOGRAPHIE

1. National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
2. Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
3. Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
4. Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
5. Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S.; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref:MD1001096	Cont.	R1: 4 x 30 mL
		R2: 2 x 20 mL

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	Homens	Mulheres
Risco menor	> 50 mg/dL	> 60 mg/dL
Risco normal	35 - 50 mg/dL	45 - 60 mg/dL
Risco elevado	< 35 mg/dL	< 45 mg/dL

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BIBLIOGRAFIA

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih W.J, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref:MD1001096	Cont.	R1: 4 x 30 mL
		R2: 2 x 20 mL