Methadone Metabolite urine

DoA Liquid Test

Methadone Metabolite (EDDP) Assay for Urine IVD

Store at 2 - 8°C.

INTENDED USE

The Methadone Metabolite Enzyme Immunoassay (EIA) is a homogenous enzyme immunoassay system intended for use in the qualitative and semi-quantitative analysis of methadone metabolite in human urine.

The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method (1, 2). Clinical consideration and professional judgement should be exercised to any result, in order to establish the appropriate treatment or therapy.

PRINCIPLE OF THE METHOD

The methadone metabolite enzyme immunoassay is a homogeneous enzyme immunoassay (9) with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labelled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, methadone metabolite-labelled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to free drug; the unbound methadone metabolite-labelled G6PDH then exhibits its maximal enzyme activity.

Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

SIGNIFICANCE OF THE TEST

Methadone is a synthetic diphenylheptanonylamine opioid that has similar analgesic activity and potency as morphine when administered parenterally. However, unlike morphine, it reliably retains its effectiveness when given orally, and tolerance and physical dependency develop slowly (3,4).

Although methadone is prescribed to relieve chronic pain, its primary medical application is the detoxification and/or maintenance treatment of narcotic or heroin addiction (3-6). The abuse potential of methadone is comparable to that of morphine due to its resembling pharmacological activity (3, 5, 7).

Methadone is readily absorbed from gastrointestinal tract when ingested, and metabolized extensively in the liver. Initial N-demethylation results in normethadone metabolite, which rapidly undergoes cyclisation followed by dehydration to form the 2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine, commonly known as EDDP. Further N-demethylation yields a secondary metabolite, the 2-ethyl-5-methyl-3, 3-diphenyl-1-pyrroline (EMDP) (8). The metabolites are secreted in urine or bile along with unchanged drug.

Current methadone immunoassays can detect only the parent drug and are subject to "false positives" from adulteration of samples for drug of abuse testing, or "false negatives" for methadone compliance testing urine samples from fast metabolizers. As a result, confirmation of the presence of EDDP by LC or GC is often required. EDDP immunoassay detection will make compliance testing easy, and also rule out the possibility adulteration of unsupervised sample collections.

REAGENTS

Antibody/Substrate Reagent (R₁): Contains mouse monoclonal antibody to methadone metabolite, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide as preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains methadone metabolite-labelled glucose-6-phosphate dehydrogenase (G6PDH) in buffer with sodium azide as preservative.

Avoid prolonged exposures of the reagent at temperatures higher than 25°C.

CALIBRATORS AND CONTROLS

Negative Human Urine (Level 0): Contains negative human urine with sodium azide as preservative (Ref.: 933010 and 933015).

EDDP Calibrator (Ref.: 933040)

EDDP Calibrator Level 0: Contains negative human urine with sodium azide as preservative.

EDDP Calibrator Level 1: Contains 200 ng/mL EDDP in human urine with sodium azide as preservative.

EDDP Calibrator Level 2: Contains 300 ng/mL EDDP in human urine with sodium azide as preservative.

EDDP Calibrator Level 3: Contains 600 ng/mL EDDP in human urine with sodium azide as preservative.

EDDP Calibrator Level 4: Contains 1000 ng/mL EDDP in human urine with sodium azide as preservative.

EDDP Calibrator for Urine Cut-Off 300 ng/mL: Contains 300 ng/mL EDDP in human urine with sodium azide as preservative. (Ref.: 932955)

EDDP Control: Contains 225 and 375 ng/mL EDDP in human urine with sodium azide as preservative (Ref.: 935955).

PREPARATION

The reagents are ready to use. No reagent preparation is required. All assay components should be stored refrigerated when not in use.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
Nº	**	Vol. R1	160
Test	DoA	Vol. R2	60
Full Name	DoA (User defined)	Vol. R3	
Standard nº	User defined	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	16	Decimal	1
Vol. Sample Increas.		Unit	User defined
Vol. Sample Dec		Slope	1
		Inter.	0
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	340	Reagent Blank	0-0
Sec. Wave.		React. Time	56-70
JUDGEMENT CRITERIA			
Absorbance	-30000/30000	Lin. Range	User defined
Incre.Test		Lin. Limit	
Decre.Test		Subs. Limit	
Prozone (Rate-Antigen)		Q1	
PC		Q2	
ABS		Q3	
		Q4	
CALIBRATOR			
CALIBRATION TYPE: Spline or Linear two points (User defined)			

This parameter calibration is stable for 30 days.

WARNING AND PRECAUTIONS

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagents used in the assay contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the reagents beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix reagents from different manufacturers.
- Do not freeze reagents.

SPECIMEN COLLECTION

Urine samples may be collected in plastic or glass containers. Some plastics may adsorb drugs. Use a fresh urine specimen for the test. If the sample can not be analyzed immediately, it may be stored refrigerated for up to 3 days. For longer storage keep the sample frozen and then thaw before use. Samples should be brought to room temperature of 18-25°C for testing. Samples with high turbidity should be centrifuged before analysis. Urine samples within the normal pH range of 5-8 can be tested without any pre-treatment. Fresh and properly stored urine samples generally are within this range. Samples with pH out of the range should be adjusted to be within this range with 1M HCl or 1M NaOH before testing. Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and both samples should be forwarded to the laboratory for testing.

Handle all urine specimens as if they are potentially infectious.

INSTRUMENTATION REQUIRED

- Spin640 autoanalyzer
- Laboratory equipment.

CALIBRATION

For qualitative determinations, the reagent must be calibrated with the selected cut-off calibrator. For semi-quantitative determinations, the reagent can be calibrated with a 5 point calibration curve.

Cut-off (ng/mL)	Calibrator levels (ng/mL)		Control levels (ng/mL)
	QUALITATIVE	SEMI- QUANTITATIVE	
300	Cut-off (300)	Level 0 (0) Level 1 (200) Level 2 (300) Level 3 (600) Level 4 (1000)	Control -25% (225) Control +25% (375)

INTERPRETATION

For qualitative determinations, the cut-off calibrator (300ng/mL) of EDDP is used as a reference for distinguishing positive from negative samples. A sample with a change in absorbance ($\Delta A/min$) equal to, or greater than, that obtained with the cut-off calibrator is considered positive. A sample with a change in absorbance value lower than that obtained with the cut-off calibrator is considered negative.

For semi-quantitative determinations, a calibration curve with multiple calibrators is required. The concentration of the EDDP in the sample may then be estimated from the calibration curve.

QUALITY CONTROL

Good laboratory practices recommend the use of control specimens to ensure proper assay performance.

The calibration curve can be validated with the controls levels 225 and 375 ng/mL (Ref. 935955) or with commercial controls.

LIMITATIONS

1. A positive result from the assay indicates only the presence of EDDP
2. Refer to the intended use statement for details of recommended confirmation testing methods.
3. The test is designed for use with human urine only.

BIBLIOGRAPHY

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988).
3. Goodman, L.S., and A. Gilman. The Pharmacological Basis of Therapeutics, 4th edition, 380, The MacMillan Co., 1970
4. Roper, J.D., D. Garside, and B.A. Goldberger, "Opiates", in Contemporary Practice in Clinical Toxicology, 2nd edition, Leslie M. Shaw, editor-in-chief, AACCC, 2000.
5. Lacy, C., L.L. Armstrong, R.J. Lipsy, and L.L. Lance. Drug Information Handbook. Hudson, OH: Lexi-Comp, 1993.
6. Katzung, B.G., Basic and Clinical Pharmacology, 6th edition, Appleton & Lange, Norwalk, CT. 1995.
7. Physician's Desk Reference, 54th edition, 2711-2713, Medical Economics Company, Montvale, NJ. 2000.
8. Baselt, R.C. and R.H. Carvey, Disposition of Toxic Drugs and Chemicals in Man, 4th edition, Chemical Toxicology Institute, Foster City, CA. 1995.
9. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, vol. 53, No. 69, pp11970 (1988).
3. Goodman, L.S., and A. Gilman. The Pharmacological Basis of Therapeutics, 4th edition, 380, The MacMillan Co., 1970
4. Ropero, J.D., D. Garside, and B.A. Goldberger, "Opiates", in Contemporary Practice in Clinical Toxicology, 2nd edition, Leslie M. Shaw, editor-in-chief, AACC, 2000.
5. Lacy, C., L.L. Armstrong, R.J.Lipsy, and L.L. Lance. Drug Information Handbook. Hudson, OH: Lexi-Comp, 1993.
6. Katzung, B.G., Basic and Clinical Pharmacology, 6th edition, Appleton & Lange, Norwalk, CT. 1995.
7. Physician's Desk Reference, 54th edition, 2711-2713, Medical Economics Company, Montvale, NJ. 2000.
8. Baselt, R.C. and R.H. Carvey, Disposition of Toxic Drugs and Chemicals in Man, 4th edition, Chemical Toxicology Institute, Foster City, CA. 1995.
9. Rubenstein, K.E., R.S. Schneider, and E.F. Ullman, Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, 47, 846 (1972)

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Ref.: MD930078

Cont.

R1: 2 x 40 mL

R2: 2 x 15 mL