**Quantitative determination of α<sub>1</sub>-antitrypsin (α<sub>1</sub>-ATRYP) IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Anti- α<sub>1</sub>-antitrypsin antibodies when mixed with samples containing α<sub>1</sub>-antitrypsin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α<sub>1</sub>-antitrypsin concentration of the patient sample, that can be quantified by comparison from a calibrator of known α<sub>1</sub>-antitrypsin concentration.

**CLINICAL SIGNIFICANCE**

α<sub>1</sub>-antitrypsin is a glycoprotein synthesized in the hepatic parenchyma cells that circulates in the bloodstream. It is the second highest proteinase inhibitor in plasma after α<sub>2</sub>-macroglobulin. α<sub>1</sub>-antitrypsin is a strong reactor with elastase, skin collagenase, chymotrypsin, plasmin, and thrombin, and also shows inhibitory activity against fungal and leukocytic proteases.

α<sub>1</sub>-antitrypsin deficiency is an inherited disorder, and occurs when both parents pass on an abnormal gene (PIZ) to their child. This deficiency is associated with a very risk for development of pulmonary emphysema and diseases of the liver - neonatal cholestasis, hepatitis, cirrhosis and hepatocellular carcinoma-. α<sub>1</sub>-antitrypsin increases in inflammatory or necrosis process. Serum levels begin to rise after approximately 24 hours and peak at 3 or 4 days if the insult is acute and short-lived.

**REAGENTS**

<b>R 1</b> Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
<b>R 2</b> Antibody	Goat serum, anti-human α <sub>1</sub> -antitrypsin, pH 7,5. Preservative.
<b>Optional</b>	Ref: 1102003 PROT CAL.

**CALIBRATION**

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

**PREPARATION**

Reagents: Ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

**ADDITIONAL EQUIPMENT**

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

**SAMPLES**

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant.

Stable 7 days at 2-8°C or 3 months at -20°C.

Do not use highly hemolyzed or lipemic samples.

**REFERENCE VALUES**

Newborn: Between 124 - 348 mg/dL.

Adults: 90 – 200 mg/dL.

Each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is

available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

**SPIN640 APPLICATION**

<b>TEST INFORMATION</b>		<b>REAGENT VOLUME</b>		
Nº	**	Vol. R1	240	
Test	ATRYP	Vol. R2	60	
Full Name	Alpha1Antytripsin	Vol. R3		
Standard nº	6	Vol. R4		
<b>SAMPLE VOLUME</b>		<b>RESULT SETUP</b>		
Vol. Sample Stand.	2	Decimal	1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
<b>REACTION PARAMETERS</b>				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
<b>CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)</b>				
CALIBRATION TYPE: SPLINE				
<b>Nº CAL DIL</b>	<b>CONCENTRATION</b>	<b>SAMPLE DIL</b>	<b>VOL DIL</b>	<b>VOL</b>
0	CAL *0,1	10,0	90,0	2,0
1	CAL *0,25	30,0	90,0	2,0
2	CAL *0,49351	30,0	90,0	4,0
3	CAL *0,7354	30,0	90,0	6,0
4	CAL	0,0	0,0	2,0

**SPIN640Plus APPLICATION**

<b>EDIT PARAMETERS</b>			
Test	ATRYP	No.	**
Full name	ANTITRIPSIN	Print name	ATRYP
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	20-21	React. Time	63 - 64
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	
<b>CALIBRATION</b>	<b>Spline</b>	<b>Water</b>	<b>0</b>
	Calibrator	*0.1	3 10 150
	Calibrator	*0.25	3 25 175
	Calibrator	*0.489	3 30 90
	Calibrator	*0.735	6 30 90
	Calibrator	(*)Cal Value	3 -- --

The Calibration is stable until 31 days. After this period the Calibration must be performed again in order to obtain good results.

**PERFORMANCE CHARACTERISTICS**

**1. Measurement range:** Up to 500 mg/dL under the described assay conditions.

Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

**2. Detection Limit:** Values less than 16 mg/dL give non-reproducible results.

**3. Prozone effect:** No prozone effect was detected upon 1000 mg/dL.

**4. Sensitivity:** Δ 3,4 mA / mg/dL.

**5. Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

<b>EP5</b>	<b>CV (%)</b>		
	34,39 mg/dL	92,7 mg/dL	181,8 mg/dL
Total	4,2%	2,6%	2,8%
Within Run	0,8%	1,1%	1,6%
Between Run	3,8%	2,4%	2,3%
Between Day	1,6%	0%	0%

**6. Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Beckman Array 360 CE. 35 samples ranging from 70 to 250 mg/dL of α<sub>1</sub>antitrypsin were assayed. The correlation coefficient (r) was 0,92 and the regression equation y = 0,84x +26,5.

The results of the performance characteristics depend on the used analyzer.

**BIBLIOGRAPHY**

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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### BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS effects of drugs on clinical laboratory tests, 4th ed. AACCPres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACCPres, 1997.

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### APPLICATION AU SPIN640

TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	ATRYP	Vol. R2	60	
Full Name	Alpha1Antytrypsin	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	2	Decimal	1 Slope 1	
Vol. Sample Inceas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0,1	10,0	90,0	2,0
1	CAL *0,25	30,0	90,0	2,0
2	CAL *0,49351	30,0	90,0	4,0
3	CAL *0,7354	30,0	90,0	6,0
4	CAL	0,0	0,0	2,0

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### BIBLIOGRAPHIE

1. Clinical Guide to Laboratory Tests, Edited by NWTietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. PesceAJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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EP5	CV (%)		
		34,39 mg/dL	92,7 mg/dL
Total	4,2%	2,6%	2,8%
Within Run	0,8%	1,1%	1,6%
Between Run	3,8%	2,4%	2,3%
Between Day	1,6%	0%	0%

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### BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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Ref.: MD1102054

Cont.

R1. Solvente: 2 x 40 mL  
R2. Anticorpo: 1 x 20 mL