## Quantitative determination of complement C3 (C3) IVD

Store 2 - 8°C.

### PRINCIPLE OF THE METHOD

Anti-human C3 antibodies when mixed with samples containing C3, form insoluble complexes. These complexes cause an absorbance change, dependent upon the C3 concentration of the patient sample, that can be quantified by comparison from a calibrator of known C3 concentration.

### CLINICAL SIGNIFICANCE<sup>1</sup>

C3 is the functional link between classical and alternative pathways of activation and it is the most concentrate component of the complement system in human plasma. Hepatic cells synthesize C3, although bacterial endotoxins induce synthesis by monocytes and fibroblasts.

Concentration C3 increases as a consequence of an acute-phase response (trauma, surgery or inflammatory process), biliary obstruction and focal glomerulosclerosis. Decreasing C3 levels are consequence of a genetic deficiency that may increase the risk of infections particularly with encapsulated bacteria, or acquired deficiency that causes vascular disorders and severe infections.

### REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
Antibody (R2)	Goat serum, anti-human C3, pH 7,5. Preservative.
Optional	Cod: 1102003 PROT CAL.

### CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

### PREPARATION

**Reagents:** Ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not use.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

### ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

### REFERENCE VALUES<sup>5</sup>

Neonates: Between 70 - 196 mg/dL.

Adults: Between 90 - 180 mg/dL.

Each laboratory should establish its own reference range.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

### SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	C3	Vol. R2	60	
Full Name	COMPLEMENT C3	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	3.0
1	CAL *0.25	30.0	90.0	3.0
2	CAL *0.4951	30.0	90.0	6.0
3	CAL *0.7354	30.0	90.0	9.0
4	CAL	0.0	0.0	3.0

### SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	C3	No.	**
Full name	C3	Print name	C3
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	20-21	React. Time	63 - 64
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	
CALIBRATION			
Spline	Water	0	
Calibrator	*0.1	3	13 115
Calibrator	*0.25	3	35 105
Calibrator	*0.489	6	35 105
Calibrator	*0.735	9	35 105
Calibrator	(* )Cal Value	3	-- --

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

### PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 400 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** Values less than 1 mg/dL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 1500 mg/dL.
- Sensitivity:**  $\Delta$  8,86 mA. mg/dL (23,8 mg/dL),  $\Delta$  4,3 mA. mg/dL (190 mg/dL).
- Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	42,98 mg/dL	118,96 mg/dL	229,5 mg/dL
Total	6,6%	2,3%	3,1%
Within Run	0,9%	0,8%	0,8%
Between Run	3,7%	2,2%	1,8%
Between Day	5,4%	0%	2,4%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using an immunoturbidimetric method from Bayer. 48 samples ranging from 50 to 200 mg/dL of C3 were assayed. The correlation coefficient (r) was 0,96 and the regression equation  $y = 1,1x - 0,6$ .

The results of the performance characteristics depend on the used analyzer.

### BIBLIOGRAPHY

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Carroll MC. Annual Review of Immunology 1998; 16: 545-568.
- Lambiris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1: 16-45.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACCPres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCPres, 1997.

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**BIBLIOGRAFIA**

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Carroll MC. Annual Review of Immunology 1998; 16: 545-568.
- Lambiris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1: 16-45.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	C3	Vol. R2	60	
Full Name	COMPLEMENT C3	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	3.0
1	CAL *0.25	30.0	90.0	3.0
2	CAL *0.4951	30.0	90.0	6.0
3	CAL *0.7354	30.0	90.0	9.0
4	CAL	0.0	0.0	3.0

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**BIBLIOGRAPHIE**

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Carroll MC. Annual Review of Immunology 1998; 16: 545-568.
- Lambiris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1: 16-45.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres.

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### BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Yang Y et al. Curr Dir Autoimmun 2004; 7: 98-132.
- Borque L et al. Clin Biochem 1983; 16: 330-333.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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Ref.: MD1102094

Cont.

R1.Solvente: 2 x 40 mL

R2.Anticorpo: 1 x 20 mL